---

# All about RT

cDNA Synthesis

---

# Selection Chart of Bioneer's cDNA Synthesis Products

## by Application

### cDNA Synthesis Kits

Application	RT PreMix	<i>CycleScript</i> <sup>™</sup> RT PreMix	<i>RocketScript</i> <sup>™</sup> RT PreMix	<i>RocketScript</i> <sup>™</sup> Cycle RT PreMix	<i>RocketScript</i> <sup>™</sup> RT PreMix, RNase H Minus
Standard RT	√	√	√	√	√
High Efficiency RT		√		√	√
RT of Secondary Structured RNA			√	√	√
RT of Long kb RNA					√

### One-step RT-PCR Kits

Application	RT-PCR	<i>RocketScript</i> <sup>™</sup> RT-PCR PreMix	<i>RocketScript</i> <sup>™</sup> RT-PCR PreMix, RNase H Minus	<i>Dual-HotStart</i> <sup>™</sup> RT-PCR PreMix	<i>Dual-HotStart</i> <sup>™</sup> RT-PCR PreMix (with UDG)	RocketPlex RT-PCR PreMix
Standard RT-PCR	√	√	√	√	√	√
RT-PCR of Secondary Structured RNA		√	√	√	√	√
High Specificity & Sensitivity RT-PCR				√	√	
Multiplex RT-PCR				√	√	√
Prevent Air-contamination of PCR Products					√	
RT-PCR of Long kb RNA			√			

### One-step RT-qPCR Kits

Application	<i>Dual-HotStart</i> <sup>™</sup> RT-qPCR PreMix	<i>GreenStar</i> <sup>™</sup> RT-qPCR PreMix
RT-qPCR by Hydrolysis Probe Method	√	
RT-qPCR by dsDNA Binding Dye Method		√
High Specificity & Sensitivity RT-qPCR	√	√

### Enzymes

Application	<i>M-MLV</i> RTase	<i>CycleScript</i> <sup>™</sup> RTase	<i>RocketScript</i> <sup>™</sup> RTase	<i>RocketScript</i> <sup>™</sup> RTase, RNase H Minus
Standard RT	√	√	√	√
High Efficiency RT		√	√	√
RT of Secondary Structured RNA			√	√
RT of Long kb RNA				√

# Selection Chart of Bioneer's cDNA Synthesis Products

## by Specification

### cDNA Synthesis Kits

Product	Fragment Size	RNase H Activity	DNase Activity	RNase Activity	GC-rich
RT PreMix	≤ 9 kb	○	X	X	●●●
<i>CycleScript</i> ™ RT PreMix	≤ 9 kb	○	X	X	●●●
<i>RocketScript</i> ™ RT PreMix	≤ 10 kb	○	X	X	●●●
<i>RocketScript</i> ™ Cycle RT PreMix	≤ 10 kb	○	X	X	●●●
<i>RocketScript</i> ™ RT PreMix, RNase H Minus	≤ 12.5 kb	X	X	X	●●●●

### One-step RT-PCR Kits

Product	Fragment Size	5' → 3' Exonuclease	3' → 5' Exonuclease	Specificity	GC-rich	Leaves 3'-A
RT-PCR PreMix	≤ 5 kb	X	X	-	●●●	○
<i>RocketScript</i> ™ RT-PCR PreMix	≤ 6 kb	X	X	-	●●●	○
<i>RocketScript</i> ™ RT-PCR PreMix, RNase H Minus	≤ 12.5 kb	○	X	-	●●●	○
<i>Dual-HotStart</i> ™ RT-PCR PreMix	≤ 3 kb	○	X	●●●●●	●●●●●	○
<i>Dual-HotStart</i> ™ RT-PCR PreMix (with UDG)	≤ 3 kb	○	X	●●●●●	●●●●●	○
RocketPlex RT-PCR PreMix	≤ 1 kb	X	X	●●●●	●●●	○

### One-step RT-qPCR Kits

Product	Fragment Size	5' → 3' Exonuclease	3' → 5' Exonuclease	Specificity	GC-rich	Leaves 3'-A
<i>GreenStar</i> ™ RT-qPCR PreMix	≤ 200 bp	○	X	●●●●	●●●●	○
<i>Dual-HotStart</i> ™ RT-qPCR PreMix	≤ 200 bp	○	X	●●●●	●●●●	○

### Enzymes

Product	Fragment Size	RNase H Activity	DNase Activity	RNase Activity	GC-rich
<i>M-MLV</i> Reverse Transcriptase	≤ 9 kb	○	X	X	●●●
<i>CycleScript</i> ™ Reverse Transcriptase	≤ 9 kb	○	X	X	●●●
<i>RocketScript</i> ™ Reverse Transcriptase	≤ 10 kb	○	X	X	●●●
<i>RocketScript</i> ™ RTase, RNase H Minus	≤ 12.5 kb	X	X	X	●●●

# Technology


## 01

### Contents

- **cDNA synthesis Technology**
  - Cyclic Temperature Reverse Transcription (CTRT)
  - *RocketScript™* Reverse Transcriptase
- **One-Step RT-PCR & RT-qPCR Technology**
  - *Dual-HotStart™* Technology

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
## Conventional RT vs Bioneer's CTRT


## CTRTR technology-applied products

제품명	카탈로그 번호
AccuPower® CycleScript™ RT PreMix	K-2044 - K-2050-2
AccuPower® RocketScript™ Cycle RT PreMix	K-2201 - K-2210
CycleScript™ Reverse Transcriptase	E-3131 - E-3132

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## RocketScript™ technology-applied products

제품명	카탈로그 번호
AccuPower® RocketScript™ RT PreMix / Master Mix	K-2101 - K-2104 / K-2105
AccuPower® RocketScript™ Cycle RT PreMix / Master Mix	K-2201 - K-2210 / K-2216
AccuPower® RocketScript™ RT PreMix, RNase H Minus / Master Mix	K-2221 - K-2248 / K-2249
AccuPower® RocketScript™ RT-PCR PreMix / Master Mix	K-2501 - K-2504 / K-2505
AccuPower® RocketScript™ RT-PCR PreMix, RNase H Minus / Master Mix	K-2231 - K-2234 / K-2235
AccuPower® RocketPlex RT-PCR PreMix / Master Mix	K-2211 - K-2214 / K-2217
AccuPower® Dual-HotStart™ RT-PCR PreMix	K-6710 - K-6713
AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG)	K-6714 - K-6715
AccuPower® GreenStar™ RT-qPCR PreMix / Master Mix (2X)	K-6400 - K-6403
AccuPower® Dual-HotStart™ RT-qPCR PreMix / Master Mix (2X)	K-6704 - K-6707
RocketScript™ Reverse Transcriptase	E-3141 - E-3142
RocketScript™ Reverse Transcriptase, RNase H Minus	E-3161 - E-3162

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
- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

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- Enzyme: M-MLVRTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 9 kb

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- 20℃

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
Figure 1. Specific amplification of 5'- UTR region of HCV with AccuPower® RT PreMix.

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- Sequencing single and double-strand DNA or RNA
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

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- Enzyme: CycleScript™ RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 9 kb

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-20°C

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
Figure 1. Reaction conditions of CTRT compared with that of conventional RT.


The conventional fixed temperature RT reactions at 42°C and 55°C and cyclic temperature RT reactions at 37°C and 50°C using HeLa cell total RNA (800 ng each) were performed 800 bp, 1.5 kb, 2.0 kb, and 2.6 kb fragments of human transferrin receptor gene were amplified with each primers. AccuPower® CycleScript™ RT PreMix series (dT<sub>20</sub> & dN<sub>12</sub>) preferably shows good results at CTRT reaction condition.

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Note: Schematic representation of the 5' UTR of a gene, with complex secondary structure, at three different temperatures. Note that RocketScript™ RTase shows full activity at 70°C allowing it to synthesize the complete gene sequence where M-MLV and other reverse transcriptase's fail.

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

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- Enzyme: RocketScript™ RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 10 kb

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-20°C

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
Figure 1. Sensitivity comparison between AccuPower® RocketScript™ RT PreMix and M-MLV RTase.

Sensitivity results of AccuPower® RocketScript™ RT PreMix using GAPDH compared with conventional reverse transcriptases. Each 10 fg-100 ng of total RNA used for RT and the same amount of RT products used for electrophoresis.

- Lane 1: 10 fg human total RNA from HeLa cell
- Lane 2: 100 fg human total RNA from HeLa cell
- Lane 3: 1 pg human total RNA from HeLa cell
- Lane 4: 10 pg human total RNA from HeLa cell
- Lane 5: 100 pg human total RNA from HeLa cell
- Lane 6: 1 ng human total RNA from HeLa cell
- Lane 7: 10 ng human total RNA from HeLa cell
- Lane 8: 100 ng human total RNA from HeLa cell

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

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- Enzyme: RocketScript™ RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 10 kb

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-20°C

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
Figure 1. Complex RNA amplification results of AccuPower® RocketScript™ Cycle RT PreMix.

Each target gene (MYC, TAP) was amplified after performing reverse transcription with AccuPower® RocketScript™ Cycle RT PreMix. Rxn. conditions: Conventional 1 hr incubation at 42°C, 50°C, or 60°C, deactivation at 95°C for 5 min.

A: M-MLV Reverse Transcriptase

B: AccuPower® RocketScript™ Cycle RT PreMix with Oligo dT<sub>20</sub>

Lane 1: 100 ng Human total RNA from HeLa cell

Lane 2: 10 ng Human total RNA from HeLa cell

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

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- Enzyme: RocketScript™ RTase
- DNase activity: No
- RNase activity: No
- RNase H activity: No
- Fragment size: ~ 12.5 kb

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-20°C

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
Figure 1. Performance of the AccuPower® RocketScript™ RT PreMix, RNase H Minus in two-step RT-PCR.

Target mRNA ranged from 500 bp to 11 kb were reverse transcribed by using AccuPower® RocketScript™ RT PreMix, RNase H Minus. With each cDNA, PCR was performed by using AccuPower® ProFi Taq PCR PreMix (Cat. No. K-2631, Bioneer). Template sizes are indicated above the gel image.

M1: Lambda/Hind III marker (Cat. No. D-1050, Bioneer)

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- cDNA library construction
- Gene expression analysis

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- Enzyme: M-MLVRTase, Top DNA Polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 5 kb

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-20℃

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Figure 1. Sensitivity comparison between AccuPower® RT-PCR PreMix and Supplier I RT-PCR kit.

Each  $10^9$ - $10^3$  copy of PSTVd used for RT-PCR and the same amount of RT-PCR products used for electrophoresis.

Lane 1:  $10^9$  copy Lane 2:  $10^8$  copy Lane 3:  $10^7$  copy  
Lane 4:  $10^6$  copy Lane 5:  $10^5$  copy Lane 6:  $10^4$  copy  
Lane 7:  $10^3$  copy Lane 8: NTC

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 2. Comparison of reproducibility test for AccuPower® RT-PCR PreMix batch 1, 2, 3 and batch 4 products using serial diluted human total RNA.

Lane 1: 10 ng human total RNA from HeLa cell  
Lane 2: 1 ng human total RNA from HeLa cell  
Lane 3: 100 pg human total RNA from HeLa cell  
Lane 4: 10 pg human total RNA from HeLa cell

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- cDNA library construction
- Gene expression analysis

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- Enzyme: RocketScript™ RTase, Top DNA Polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 6 kb

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-20°C

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
Figure 1. Performance comparison between AccuPower® RocketScript™ RT-PCR PreMix and competitor RT-PCR kits.

Complex secondary structure RNA species amplification by RocketScript™ is outstanding compared to the leading competitor's reverse transcriptase. RT reactions were performed according to each manufacturer's recommendations.

- Lane 1: 10 ng of total RNA from HeLa cell
- Lane 2: 1 ng of total RNA from HeLa cell
- Lane 3: 100 pg of total RNA from HeLa cell
- Lane 4: 10 pg of total RNA from HeLa cell

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- Low copies detection
- RNA virus detection
- Gene expression analysis

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- Enzyme: RocketScript™ RTase (RNase H Minus) ProFi Taq DNA Polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: Yes
- 3' - A overhang: Yes
- RNase H activity: No
- Fragment size: ~ 12.5 kb

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- 20°C

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
Figure 1. Performance of the AccuPower® RocketScript™ RT-PCR PreMix, RNase H Minus.

AccuPower® RocketScript™ RT-PCR PreMix, RNase H minus was used to amplify cDNAs which size are between 2.5 kb to 12.5 kb. DNA sizes are indicated above the picture. Reactions were performed with 100 ng total RNA from HeLa cells.

M1: Lambda/Hind III marker (Cat. No. D-1050, Bioneer)

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA Polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 3 kb

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-20°C

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
Figure 1. Comparison of PCR amplification specificity between AccuPower® Dual-HotStart™ RT-PCR PreMix from Bioneer and other suppliers' Hotstart PCR master mix.

Lane 1: HCV 10<sup>7</sup>                      Lane 3: HCV 10<sup>5</sup>  
 Lane 2: HCV 10<sup>6</sup>                      Lane 4: HCV 10<sup>4</sup>


Figure 2. Comparison of PCR amplification sensitivity between AccuPower® Dual-HotStart™ RT-PCR PreMix from Bioneer and other suppliers' Hotstart PCR master mix.

Target: human PGK1  
 Lane 1: 10 ng of human RNA  
 Lane 2: 1 ng of human RNA  
 Lane 3: 100 pg of human RNA  
 Lane 4: 10 pg of human RNA

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA Polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 3 kb

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-20°C

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
Figure 1. Comparison of amplification quality using PCR products (not including uracil base or including uracil base) between AccuPower® Dual-HotStart™ RT-PCR PreMix and AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG).

Lanes 1: 10<sup>7</sup>      Lanes 2: 10<sup>6</sup>      Lanes 3: 10<sup>5</sup>  
 Lanes 4: 10<sup>4</sup>      Lanes 5: 10<sup>3</sup>      Lanes 6: 10<sup>2</sup>  
 Lanes 7: 10      Lanes 8: NTC

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- Multiplex RT-PCR
- Low copy detection
- Gene expression analysis

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- Enzyme: RocketScript™ RTase, HotStart Top DNA Polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 1 kb

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-20°C

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
Figure 1. Single RT-PCR and multiplex RT-PCR using AccuPower® RocketPlex RT-PCR PreMix.

M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 200 bp

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
-20°C

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
#### BIONEER

#### Company A

#### Company B

RNA (ng)	BIONEER	Company A	Company B
100	17.79	17.38	19.00
10	21.46	20.82	22.67
1	25.00	24.67	26.17
0.1	27.31	27.45	29.85
0.01	30.34	30.34	32.82
0.001	33.28	33.58	35.58
NTC	UD	37.32	36.28

Figure 1. Comparison of specificity between AccuPower® GreenStar™ RT-qPCR PreMix and Supplier's master mixture.

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- Gene expression profiling
- Target RNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA Polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes

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-20°C

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
Figure 1. High specificity of AccuPower® Dual-HotStart™ RT-qPCR PreMix.

Experiment with HCV target. 10 fold serial dilution of template RNA ( $10^6$  - 10 copies spiked in human total RNA).

Conventional hotstart qPCR always generate non-specific amplification at low template concentration, which deteriorate the sensitivity of qPCR. AccuPower® RocketScript™ RT-qPCR PreMix accurately amplifies target RNA without non-specific amplification, even at low concentration of template.

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