All about Bioneer's

Real-Time PCR

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Real-Time quantification of DNA and cDNA targets
Gene expression profiling
Gene functional analysis
Microbial & viral pathogen detection

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Enzyme: Top DNA Polymerase
5' → 3' exonuclease: No | 3' → 5' exonuclease: No | 3' - A overhang: Yes

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
Figure 1. Comparison of amplification efficiency between *AccuPower® GreenStar™* qPCR PreMix and other company's master mix products. Amplification curve of *AccuPower® GreenStar™* qPCR PreMix and other company's master mix kit. Lambda DNA primers were added into *AccuPower® GreenStar™* qPCR PreMix and other company's master mix kit. A series of Lambda DNA positive control diluents were tested. Reaction mixtures were prepared and qPCR was performed according to each company's protocol. All data were obtained using ABI 7500 Fast Real-time PCR machine (Applied Biosystems co.).


Figure 2. Real-Time PCR data of *AccuPower® GreenStar™* qPCR PreMix. *AccuPower® GreenStar™* qPCR PreMix provides at least 7 orders of magnitude in dynamic range (10 fg~10 ng /rxn).
A: Amplification curve of *AccuPower® GreenStar™* qPCR PreMix. Lambda DNA primers were added into *AccuPower® GreenStar™* qPCR PreMix. A series of lambda DNA positive control diluents were tested.
B: Standard curve of *AccuPower® GreenStar™* qPCR PreMix. $R^2 = 0.998$, PCR efficiency - 95%
C: Melting curve analysis of *AccuPower® GreenStar™* qPCR PreMix. The melting curve shows that only single amplified PCR product was obtained in all template range. All data were obtained using *Exicycler™* 96 Real-Time Quantitative Thermal Block (Cat. no. A-2060, Bioneer).

•Using a ABI 7500


Figure 3. Comparison of the specificity of dsDNA binding dye based Real-Time PCR. Amplification of a 90 bp target gene was detected using serially diluted LP (*Legionella Pneumoniae*) genomic DNA (10^5 dilution: $10^5 \sim 10^1$ copies) with *AccuPower® 2X GreenStar™* qPCR Master Mix and other commercial products. As shown in Figure, very small amount of primer dimers was appeared in *AccuPower® 2X GreenStar™* qPCR Master Mix than other kits.

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
Low copy viral/bacterial pathogen load determination in an earlier stage, Low copy mRNA amplification, Low copy target RNA quantification, RNA amplification for microarray and NGS

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5' → 3' exonuclease: Yes | 3' → 5' exonuclease: No | 3' - A overhang: Yes
Fragment size: 200 bp

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Gene expression profiling, Target DNA quantification, Microbial detection, Viral/bacterial pathogen load determination, Evaluation of primer pair performance for probe-based Real-Time PCR

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Enzyme: Taq DNA Polymerase

5' → 3' exonuclease: Yes | 3' → 5' exonuclease: No | 3' - A overhang: Yes

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
Figure 1. Comparison of detection sensitivity between AccuPower® DualStar™ qPCR PreMix and other company's master mix products.

West Nile Virus (WNV) primers and Hydrolysis-based probe were added into AccuPower® DualStar™ qPCR PreMix and other Company's master mix kits. A series of WNV positive control diluents were tested. Reaction mixtures were prepared and qPCR was performed according to each Company's protocol. All data were obtained using ABI 7500 Fast Real-Time PCR system (Applied Biosystems co.).

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Gene expression profiling, Target DNA quantification, Microbial detection, Viral/bacterial pathogen load determination, Evaluation of primer pair performance for probe-based Real-Time PCR

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Enzyme: HotStart Taq DNA Polymerase
5' → 3' exonuclease: Yes | 3' → 5' exonuclease: No | 3' - A overhang: Yes

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Figure 1. Comparison of amplification quality between AccuPower® Plus DualStar™ qPCR PreMix and other Company's Real-Time qPCR kit.

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Gene expression profiling, Target DNA quantification, Microbial detection, Viral/bacterial pathogen load determination, Evaluation of primer pair performance for probe-based Real-Time PCR

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Enzyme: HotStart Taq DNA Polymerase

5' → 3' exonuclease: Yes | 3' → 5' exonuclease: No | 3' - A overhang: Yes

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Figure 1. Efficiency of uracil DNA glycosylase using PCR products (including dN/dU base).

PCR product type	Linearity	Efficiency	10 ⁸	10 ⁷	10 ⁶	10 ⁵	10 ⁴	10 ³	10 ²	NTC
dNTPs	0.9988	98	12.81	16.41	19.72	22.84	26.37	29.76	33.48	UD
dNTPs including dUTPs	·	·	UD	UD						

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Gene expression profiling, Target DNA quantification, Microbial detection, Viral/bacterial pathogen load determination

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Enzyme: RocketScript™ RTase, HotStart Taq DNA Polymerase

5' → 3' exonuclease: Yes | 3' → 5' exonuclease: No | 3' - A overhang: Yes

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
Figure 1. High specificity of *AccuPower*® *Dual-HotStart*™ RT-qPCR PreMix.

Experiment with HCV target. 10 fold serial dilution of Template RNA (10^6 copies~10 copies) spiked in Human Total RNA. Conventional Hot start qPCR always generate non-specific amplification at low template concentration, which deteriorate the sensitivity of qPCR.

Dual-HotStart™ RT-qPCR accurately amplifies target RNA without non-specific amplification, even at low concentration of template.


Figure 2. High sensitivity of *AccuPower*® *Dual-HotStart*™ RT-qPCR PreMix.

Experiment with HIV target. 10-fold serial dilution of template RNA (10^{10} copies~10 copies) spiked in human total RNA).


Figure 3. PCR inhibitor study using *AccuPower*® *Dual-HotStart*™ RT-qPCR PreMix.

Human blood has various inhibitors which suppress PCR, such as hemoglobin and/or blood-EDTA. The blocking effect of these PCR inhibitors in RT-qPCR using *AccuPower*® *Dual-HotStart*™ RT-qPCR PreMix was examined, when added directly to the reaction at a final concentration of 10–10,000 ppm. Our results showed that *AccuPower*® *Dual-HotStart*™ RT-qPCR PreMix is either not affected (Hemoglobin) or affected less than (Blood-EDTA) Company A's product by PCR inhibitors in the RT-qPCR reaction.

PCR Inhibitor	BIONEER	Company A
	Totally inhibition (ppm)	
Blood-EDTA	1,000	100
Hemoglobin	*	*

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- Pacific Blue, Oregon Greens, AMCA, Bodipy 유도체
- FAM, JOE, TET, HEX, VIC, NED, TAMRA, ROX, Texas Red
- Cyanine 3, Cyanine 3.5, Cyanine 5, Cyanine 5.5

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High-Throughput Real-Time PCR System

- Superior 5-color Real-Time Quantitative PCR System

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
Figure 1. Excellent Uniformity.
Fluorescence data using 10^6 copies of IRF3 gene (FAM labeled) in each of 96 well positions. The average Ct of 96 well is 21.8 and the Ct variation range is 0.19.


Figure 2. Wide dynamic range.
Graph shows standard curve of tenfold serial dilutions of 10 copies to 10^{11} copies MMP9 gene (FAM labeled). The PCR efficiency generated by the standard curve is 103%.


Figure 3. Precise discrimination.
Fluorescence data from a series of 1.33-fold dilutions of TMV gene (10^6 copies) amplified using reporter dyes to check one target: FAM/TMV. The PCR efficiency generated by the standard curve is 101%.


Figure 4. Real 5-color Multiplexing.
5 target genes can be detected in a single tube (FAM:T. vaginalis, TET: M. Hominis, TAMRA: TMV, Texas Red: HSV type1, Cyanine5: HSV type2).

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
Figure 1. Excellent Uniformity.
qPCR result using 1×10^6 copies of Lambda DNA (FAM labeled) in each of 384 well positions. The average Ct of 384 well is 21.6 and the Ct variation range is 0.43.


Figure 2. Wide dynamic range.
Ct values of 10-fold diluted samples show a wide dynamic range of quantification. Fluorescence data from a series of 10-fold dilution of PGK1 DNA (10^{10} copies) amplified using reporter dyes to check on target: FAM/PGK1.


Figure 3. Precise discrimination.
Exicycler™ 384 provide sensitive detection and precise target discrimination down to 2-fold differences. Fluorescence data from a series of 2-fold dilution of CSF2 DNA (10^6 copies) amplified using reporter dyes to check one target: TET/CSF2.


Figure 4. Real 5-color Multiplexing.
5 target genes can be detected in a single tube with a minimum volume of 5 μ l (FAM:T. vaginalis, TET: M. Hominis, TAMRA: TMV, Texas Red: HSV type1, Cyanine5: HSV type2).

dsDNA Binding Dye Type Kit

카탈로그 번호	제품명
AccuPower® GreenStar™ qPCR PreMix & 2X Master Mix	
K-6200	AccuPower® GreenStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6201	AccuPower® GreenStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6202	AccuPower® GreenStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96)
K-6203	AccuPower® GreenStar™ qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6204	AccuPower® GreenStar™ qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, ABI7500, optical film included
K-6210	AccuPower® GreenStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6211	AccuPower® GreenStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6212	AccuPower® GreenStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96)
K-6213	AccuPower® GreenStar™ qPCR PreMix, 20 µl/rxn, 96-well plate, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6214	AccuPower® GreenStar™ qPCR PreMix, 20 µl/rxn, 96-well plate, 96 rxn, ABI7500, optical film included
K-6251	AccuPower® 2X GreenStar™ qPCR Master Mix, 50 µl/rxn, 100 rxn, 80X ROX Dye (0.1 ml X 1 ea)
K-6252	AccuPower® 2X GreenStar™ qPCR Master Mix, 50 µl/rxn, 200 rxn, 80X ROX Dye (0.1 ml X 1 ea)
K-6253	AccuPower® 2X GreenStar™ qPCR Master Mix, 50 µl/rxn, 100 rxn, without ROX Dye
K-6254	AccuPower® 2X GreenStar™ qPCR Master Mix, 50 µl/rxn, 200 rxn, without ROX Dye
AccuPower® GreenStar™ RT-qPCR PreMix & Master Mix	
K-6400	AccuPower® GreenStar™ RT-qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6403	AccuPower® GreenStar™ RT-qPCR Master Mix (2X), 2.5 ml, 100 rxn

Hydrolysis Probe Type Kit

카탈로그 번호	제품명
AccuPower® DualStar™ qPCR PreMix	
K-6100	AccuPower® DualStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6101	AccuPower® DualStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6102	AccuPower® DualStar™ qPCR PreMix, 20 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96)
K-6103	AccuPower® DualStar™ qPCR PreMix, 20 µl/rxn, 96-well plate, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6104	AccuPower® DualStar™ qPCR PreMix, 20 µl/rxn, 96-well plate, 96 rxn, ABI7500, optical film included
K-6110	AccuPower® DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6111	AccuPower® DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6112	AccuPower® DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96)
K-6113	AccuPower® DualStar™ qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6114	AccuPower® DualStar™ qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, ABI7500, optical film included
AccuPower® Plus DualStar™ qPCR PreMix & Master Mix	
K-6600	AccuPower® Plus DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6601	AccuPower® Plus DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6602	AccuPower® Plus DualStar™ qPCR PreMix, 50 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96), optical film included
K-6603	AccuPower® Plus DualStar™ qPCR Master Mix (2X), 2.5 ml, 100 rxn
AccuPower® Plus DualStar™ qPCR PreMix & Master Mix (with UDG)	
K-6605	AccuPower® Plus DualStar™ qPCR PreMix(with UDG), 50 µl/rxn, 8-tube strips, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6606	AccuPower® Plus DualStar™ qPCR PreMix (with UDG), 50 µl/rxn, 8-tube strips, 96 rxn, ABI7500, optical film included
K-6607	AccuPower® Plus DualStar™ qPCR PreMix (with UDG), 50 µl/rxn, 8-tube strips, 96 rxn, Opticon (CFX 96), optical film included
K-6608	AccuPower® Plus DualStar™ qPCR Master Mix (2X) (with UDG), 2.5 ml, 100 rxn
AccuPower® Dual-HotStart™ RT-qPCR PreMix & Master Mix	
K-6704	AccuPower® Dual-HotStart™ RT-qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, <i>Exicycler™</i> 96, optical film included
K-6705	AccuPower® Dual-HotStart™ RT-qPCR PreMix, 50 µl/rxn, 96-well plate, 96 rxn, ABI7500, optical film included
K-6706	AccuPower® Dual-HotStart™ RT-qPCR PreMix, 50 µl/rxn, 8-tube strips with cap, 96 rxn, Opticon (CFX 96), optical film included
K-6707	AccuPower® Dual-HotStart™ RT-qPCR Master Mix (2X), 2.5 ml, 100 rxn

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