

High-Throughput Automated Nucleic Acid Extraction System using *MagListo™* Kit (96 Samples)




MagListo™ X *ExiPrep™*96 Lite



- **Prevention of contamination**
 - UV-sterilization between extraction
- **Innovative magnetic rod system**
 - Superfast movement - Prep within 30 min
 - Multiple combination - 32, 64 or 96 preps
- **High yield and purity**
 - Multiple heating blocks
 - High efficiency of sample lysis/final elution
- Ease-of-use**
 - 7" - touch screen with built-in OS
 - Preprogrammed, sample-optimized protocol (User-modification possible)

*ExiPrep™*96 Lite can automatically extract nucleic acid using magnetic rods and *AccuNanoBead™* magnetic nanobeads. So, *MagListo™* series can be applied to *ExiPrep™* 96 Lite (A-5250) for high-throughput nucleic acid extraction.



Process of nucleic acid extraction using magnetic rods

Automation of Nucleic acid extraction using *ExiPrep™96 Lite*



- To apply *MagListo™* series to *ExiPrep™96 Lite* (A-5250), please refer to the guide below.
- For more information, please visit www.bioneer.com and refer to User's Guide of the product.

Step	Image	Description
Prepare buffer cartridge	<p>→ 96-Well dome plate</p>	Dispense adequate volume of buffers and samples to 96-well dome plate (Cat. No. 90060, 90061, 90062, 90063).
Set up	<p>→ Magnetic Rods → Buffer Cartridges Layer 3 → Buffer Cartridges Layer 2 → Buffer Cartridges Layer 1</p>	<ol style="list-style-type: none"> 1 Touch "Plate icon" to open "Buffer Cartridges Layer". 2 Set up the buffer cartridge on "Buffer Cartridges Layer" of <i>ExiPrep™96 Lite</i>.
	<p>Rod Cover Plate Magnetic Rod Cover</p>	<ol style="list-style-type: none"> 1 Touch "Plate icon" to open "Rod Cover Plate". 2 Set the magnetic rod cover to "Magnetic Rod Cover Plate" of <i>ExiPrep™96 Lite</i>.
Select protocol		<ol style="list-style-type: none"> 1 Touch "Standard Protocol". 2 Select adequate protocol for your sample. 3 Touch "select icon".
Run protocol		<ol style="list-style-type: none"> 1 Touch "Run icon". 2 Nucleic acid extraction will proceed automatically and it takes about 30 ~ 50 min depending on the protocol.
Purified nucleic acid	<p>→ Elution buffer cartridge</p>	<ol style="list-style-type: none"> 1 Purified nucleic acid will be in elution buffer cartridge (on the bottom of <i>ExiPrep™96 Lite</i>).