B_{RNAi}

- 01. siRNA Synthesis Service
- 02. miRNA Synthesis Service
- 03. *In vivo* siRNA Custom Service: SAMiRNA™

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
✉ E-mail: siRNA@bioneer.co.kr

01. siRNA Synthesis Service

<i>AccuTarget™</i> siRNAs	39
Genome-wide Predesigned siRNA & Primer	
<i>AccuTarget™</i> Genome-wide Predesigned siRNA Library	40
<i>AccuTarget™</i> Real-Time PCR Primers Library	41
Premade siRNA & Primer	
<i>AccuTarget™</i> Premade Human siRNA Sets	42
<i>AccuTarget™</i> Real-Time PCR Primers for Premade siRNA Sets	44
Custom siRNA	
<i>AccuTarget™</i> Custom Designed siRNA Synthesis	45
Control siRNA	
<i>AccuTarget™</i> Control siRNAs (Positive/Negative)	47
siRNA Oligonucleotide FAQs	51
Users Protocol	54

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
Turbo si-Designer - Bioneer's proprietary siRNA design algorithm

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
Figure 1. siRNA knockdown efficiency of siRNAs designed by Turbo si-Designer was analyzed by Northern blot and Real-Time PCR analysis.

A) siRNA knockdown efficiency of high NGIC score siRNAs.

B) siRNA knockdown efficiency of low NGIC score siRNAs.

AccuTarget™ Genome-wide Predesigned siRNA Library

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1. siRNA knockdown efficiency data: NC (AccuTarget™ Negative Control) and siRNA concentration at 100 nM
2. Transfection efficiency data: PC (AccuTarget™ GAPDH/GFP/Luciferase siRNA) and NC (AccuTarget™ Fluorescein-labeled Negative Control)

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- Functional genomics and proteomics research
- Gene expression studies
- Array analysis

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
Figure 2. siRNA knockdown efficiency of AccuTarget™ Genome-wide Predesigned siRNA.

AccuTarget™ Predesigned siRNAs are highly effective. To determine siRNA knockdown efficiency of predesigned siRNAs, HeLa cells were transfected with siRNAs at 100 nM concentration. By 24 hrs post-transfection, total RNA was isolated and the level of target mRNA was measured by qPCR. This data demonstrates the effectiveness of the Turbo si-Designer algorithm: 83.8% of tested siRNAs induced >70% siRNA knockdown and 38.1% of tested siRNAs elicited >90% knockdown.

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
A Fluorescence Analysis Data

B Melting Analysis Data


Figure 3. Real-Time PCR validation test of human oxidoreductase using AccuTarget™ Human Oxidoreductase Real-Time PCR Primer set.

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
AccuTarget™ Real-Time PCR Primer for Individual Gene

제품명	Reactions
AccuTarget™ Human Real-Time PCR Primer for Individual Gene	100
AccuTarget™ Human Real-Time PCR Primer for Individual Gene	200

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- Pathway analysis and target identification and validation
- Drug target HTS (High Throughput siRNA Screening)


Figure 4. Knockdown efficiency of AccuTarget™ siRNA Library.

To determine knockdown efficiency, HeLa cells were transfected with individual siRNAs at 100 nM concentration. By 24 hrs post-transfection, total RNA was isolated and the level of target mRNA was measured by qPCR. This data demonstrates the effectiveness of the Turbo si-Designer algorithm: 83.8% of tested siRNAs induced >70% knockdown and 38.1% of tested siRNAs elicited >90% knockdown.


Figure 5. siRNA knockdown efficiency of AccuTarget™ Human Validated Cell cycle siRNA Set. HeLa cells were transfected with siRNAs targeting 33 different genes at a concentration of 20 nM. Total RNA was isolated and the level of target mRNA was measured by qPCR. This data demonstrates the effectiveness of the Turbo si-Designer algorithm.

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AccuTarget™ Premade Human siRNA Sets			
카탈로그 번호	제품명	No. of Genes	No. of siRNAs
SHS-0010	AccuTarget™ Human Antioxidant siRNA set	38	114
SHS-0020	AccuTarget™ Human Apoptosis siRNA set	290	870
SHS-0250	AccuTarget™ Human Cancer siRNA set	1,157	3,471
SHS-0030	AccuTarget™ Human Caspase siRNA set	37	111
SHS-0040	AccuTarget™ Human Cell cycle siRNA set	112	336
SHS-0050	AccuTarget™ Human Cyclase siRNA set	21	63
SHS-0060	AccuTarget™ Human Cytochrome P450 siRNA set	52	156
SHS-0070	AccuTarget™ Human Deaminase siRNA set	22	66
SHS-0080	AccuTarget™ Human GPCR signaling pathway siRNA set	727	2,181
SHS-0090	AccuTarget™ Human Helicase siRNA set	114	342
SHS-0100	AccuTarget™ Human Isomerase siRNA set	104	312
SHS-0110	AccuTarget™ Human Kinase siRNA set	699	2,097
SHS-0120	AccuTarget™ Human Ligase siRNA set	272	816
SHS-0130	AccuTarget™ Human Lyase siRNA set	123	369
SHS-0140	AccuTarget™ Human Motor siRNA set	122	366
SHS-0150	AccuTarget™ Human NF-kB pathway siRNA set	37	111
SHS-0160	AccuTarget™ Human Nucleic acid binding siRNA set	2,573	7,719
SHS-0170	AccuTarget™ Human Oxidoreductase siRNA set	551	1,653
SHS-0180	AccuTarget™ Human Peptidase siRNA set	491	1,473
SHS-0190	AccuTarget™ Human Phosphatase siRNA set	188	564
SHS-0200	AccuTarget™ Human Receptor siRNA set	1,516	4,548
SHS-0210	AccuTarget™ Human Transferase siRNA set	1,428	4,284
SHS-0220	AccuTarget™ Human Transporter siRNA set	1,021	3,063
SHS-0230	AccuTarget™ Human Tubulin siRNA set	20	60
SHS-0240	AccuTarget™ Human Ubiquitin siRNA set	77	231

AccuTarget™ Real-Time PCR Primers for Premade siRNA Sets

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Figure 6. MALDI-TOF Mass Spectrometry Analysis of a custom siRNA.
All siRNAs are subjected to MALDI-TOF Mass Spectrometry to ensure its quality.


Figure 7. PAGE data of annealed double-stranded custom siRNA.
Complementary single-strand RNA strands were hybridized to form siRNA duplex and analyzed by 15% non-denaturing PAGE.
SS: single-strand RNA, DS: double-strand siRNA.

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AccuTarget™ Custom Designed siRNA Synthesis

○ Modification available for custom siRNA

5'- Modification	3'- Modification	Internal Modification	정제 방법	보장량 (nmole)
5' Fluorescein	3' Fluorescein	Phosphorothioate	HPLC	10
5' Phosphorylation	3' Phosphorylation	2'-OMe-rA, rC, rG, rU		
5' Biotin	3' Biotin	2'-F(A)		
5' Amine	3' Amine	2'-F(U)		20
5' TAMRA	3' TAMRA	2'-F(G)		
5' Thiol	3' Thiol	2'-F(C)		
5' Cyanine3	3' DABCYL	Inosine		50
5' Cyanine5	3' Cholesterol	Deoxy-abase		
5' PEG 2000	3' PEG 2000	Chimeric DNA		
5' Cyanine5.5	3' Cyanine5.5	Phosphorothioate		

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
A) Northern blotting

B) qPCR


Figure 8. Effects of Human GAPDH Positive Control siRNA.

HeLa cells were transfected separated with AccuTarget™ Human GAPDH Positive Control and Negative Control siRNA using lipofectamine 2000 (Invitrogen) at a final concentration of 100 nM. Total cellular RNA was isolated from transfected cells 24 hrs after transfection and subjected to northern blot and qPCR analysis. As can be seen from Figure 8-B, about 3% GAPDH mRNA remained.

2. GFP-siRNA


Figure 9. HeLa cells in a 24-well plate were co-transfected with 200 ng of CMV-GFP plasmid and 10 nM of GFP siRNA using lipofectamine 2000 transfection reagent.

Next day, the expression of GFP was observed by using a Nikon Eclipse TS100 epifluorescence microscope. In contrast to bright green fluorescence of GFP protein in NC-siRNA-transfected cells, no fluorescence was detected from GFP-siRNA-transfected cells, indicating efficient knockdown of GFP by using our positive control GFP-siRNA.

AccuTarget™ Control siRNAs (Positive / Negative)

3. Luciferase-siRNA


Figure 10. HeLa cells in a 6-well plate were co-transfected with 400 ng of CMV-luc plasmid and 10 nM of luciferase siRNA using lipofectamine 2000 transfection reagent.

Next day, cells were harvested and assayed for luciferase activity. As shown in Figure. 6, co-transfection with our positive control luciferase siRNA led to efficient knockdown of luciferase activity (85%-95% knockdown compared to luciferase activity of NC-siRNA-transfected cells).

4. Mouse Positive Control siRNA


Figure 11. Effects of Mouse Lamin A/C Positive Control siRNA.

Hepa1-6 cells were transfected separately with *AccuTarget™* Mouse Lamin A/C Positive Control and Negative Control siRNA using Lipofectamine RNAimax (Invitrogen) at a final concentration of 1 nM. Total cellular RNA was isolated from transfected cells 24 hrs after transfection and subjected to real-time PCR analysis.


Figure 12. Effects of Mouse PPib (Cyclophilin B) Positive Control siRNA.

Hepa 1-6 & NIH3T3 cells were transfected separately with *AccuTarget™* Mouse PPib (Cyclophilin B) Positive Control and Negative Control siRNA using Lipofectamine RNAimax (Invitrogen) at a final concentration of 5 & 1 nM. Total cellular RNA was isolated from transfected cells 24 hr after transfection and subjected to real-time PCR analysis.

AccuTarget™ Control siRNAs (Positive / Negative)

Negative control siRNA


Figure 13. HeLa cells transfected with FITC-labeled siRNA (Cat. No. SN-1021, Bioneer) was observed by confocal microscopy.

The fluorescent cells indicate that the target cells were successfully transfected with the siRNA.

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AccuTarget™ Positive Control siRNAs

카탈로그 번호	제품명	정제 방법	보장량 (nmole)
SP-1001	AccuTarget™ GAPDH Positive Control siRNA	Bio-RP	5
SP-1002			10
SP-1003			20
SP-1011		HPLC	5
SP-1012			10
SP-1013			20
SP-2001	AccuTarget™ GFP Positive Control siRNA	Bio-RP	5
SP-2002			10
SP-2003			20
SP-2011		HPLC	5
SP-2012			10
SP-2013			20
SP-3001	AccuTarget™ Luciferase Positive Control siRNA	Bio-RP	5
SP-3002			10
SP-3003			20
SP-3011		HPLC	5
SP-3012			10
SP-3013			20
SP-4001	AccuTarget™ Mouse Lamin A/C Positive Control	Bio-RP	5
SP-4002			10
SP-4003			20
SP-4011		HPLC	5
SP-4012			10
SP-4013			20
SP-5001	AccuTarget™ Mouse cyclophilin B Positive Control	Bio-RP	5
SP-5002			10
SP-5003			20
SP-5011		HPLC	5
SP-5012			10
SP-5013			20

AccuTarget™ Control siRNAs (Positive / Negative)

AccuTarget™ Negative Control siRNAs

카탈로그 번호	제품명	정제 방법	보장량 (nmole)
SN-1001	AccuTarget™ Negative Control siRNA	Bio-RP	5
SN-1002			10
SN-1003			20
SN-1011		HPLC	5
SN-1012			10
SN-1013			20
SN-1021	AccuTarget™ Fluorescein-labeled Negative Control siRNA	HPLC	5
SN-1022			10
SN-1023			20

AccuTarget™ Control siRNA Sets

카탈로그 번호	제품명	정제 방법	보장량 (nmole)
SS-1001	AccuTarget™ GAPDH Control siRNA Set	Bio-RP	5 nmole positive control + 2 nmole negative control
SS-1002	AccuTarget™ GFP Control siRNA Set		
SS-1003	AccuTarget™ Luciferase Control siRNA Set		
SS-1011	AccuTarget™ GAPDH Control siRNA Set	HPLC	
SS-1012	AccuTarget™ GFP Control siRNA Set		
SS-1013	AccuTarget™ Luciferase Control siRNA Set		

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Dilution Protocol

1. Briefly centrifuge tubes (or multi-well plates) containing siRNA to ensure that the siRNA pellet is located at the bottom of the tube.
2. Dissolve siRNAs to a convenient stock concentration using the recommended volume of DEPC-D.W. (or RNase-free water) shown in Table 1.
3. Pipette the solution up and down 3-5 times (or vortex briefly).
4. Briefly centrifuge tubes (or multi-well plates) containing siRNA to ensure that the solution is collected at the bottom of the tube.
5. Aliquot the siRNAs into small volumes and store at -20°C. siRNA is stable for 1 year under the specified storage condition. For best results, limit freeze-thaw events for each tube no more than five.

Table 1. Recommended siRNA resuspension volumes and concentrations.

siRNA Amount (nmol)	DEPC-D.W. volume(μl) for desired final concentration	
	100 μM stock	20 μM stock
10	100	500
20	200	1000
50	500	Exceeds tube volume
100	1000	

Transfection Protocol

- * We use the Lipofectamine™ RNAiMAX (Invitrogen; Cat. no. 13778) and HeLa cell for transfection procedure.
- * This protocol is fixed at 6-well plate *in vitro* culture condition (if you want to change this condition, you have to consider the relative surface area (Table 2) and Invitrogen protocol, when you are seeding the cells into the culture dish).

1. One day (24 hrs) before transfection, plate 3.0×10^5 HeLa cells in each well with 2.5 ml of growth medium without antibiotics such that they will be 50~60% confluent at the time of transfection.
2. Remove the growth medium from the 6-well plate before transfection. And add the 500 μl fresh growth medium without serum in each well.
3. For each well to be transfected, prepare siRNA duplex-Lipofectamine™ RNAiMAX complexes as follows.
 - 3-1. Dilute siRNA duplex (making final concentration as 5 nM~100 nM) in 250 μl growth medium (or Opti-MEM I Reduced Serum medium) without serum. Mix gently by vortex.
 - 3-2. Mix Lipofectamine™ RNAiMAX gently before use, then dilute 3.5 μl in 250 μl medium (or Opti-MEM I Reduced Serum medium) without serum. Incubate this solution 5 min at room temperature.
 - 3-3. Combine the diluted siRNA duplex with the diluted Lipofectamine™ RNAiMAX. Mix and incubate for 20 min at room temperature.
4. Add the mixture to each well containing HeLa cells, which result 1 ml as total volume. Mix gently by hand rocking the plate back and forth.
5. Incubate the cells for 5~6 hr at 37°C in CO₂ incubator.
6. Change the medium with fresh one containing serum and incubate the cells 24~48 hr until you are ready to assay for gene knockdown.

Table 2. The relative surface area of *in vitro* cell culture dish and culture media volume.

Culture Vessel	Relative surface area	Volume of plating medium
96-well	0.2	100 μl
48-well	0.4	200 μl
24-well	1	500 μl
6-well	5	2.5 ml
60 mm	10	5 ml
100 mm	30	10 ml

02. miRNA Synthesis Service

<i>AccuTarget™</i> miRNAs	56
miRNA Mimics & Inhibitors	
<i>AccuTarget™</i> miRNA mimics & inhibitors	57
Control miRNA	
<i>AccuTarget™</i> Control miRNAs (Positive / Negative)	59
miRNA mimic FAQs	61

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AccuTarget™ miRNA mimics & inhibitors

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
Figure 3. Performance of miRNA mimics positive and negative controls.

AccuTarget™ miRNA mimic Positive & Negative Controls were transfected at 20 nM using Lipofectamine™ RNAiMAX into HeLa cell lines and assessed for their ability to decrease target mRNA levels. Down-regulation of GAPDH was determined using the Real-Time quantitative RT-PCR at 48 hr post-transfection using Bioneer's Exicycler™ 96 qPCR instrument.

AccuTarget™ Control miRNAs (Positive / Negative)

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카탈로그 번호	제품명	정제 방법	보장량 (nmole)
SMC-1001	AccuTarget™ miRNA Housekeeping Positive Control (GAPDH)	Bio-RP	5
SMC-1002			10
SMC-1003			20
SMC-2001	AccuTarget™ miRNA mimic Negative Control #1		5
SMC-2002			10
SMC-2003			20
SMC-3001	AccuTarget™ miRNA mimic Negative Control #2		5
SMC-3002			10
SMC-3003			20
SMC-4001	AccuTarget™ Fluorescein-labeled miRNA mimic Negative Control siRNA #1	HPLC	5
SMC-4002			10
SMC-4003			20
SMC-5001	AccuTarget™ Fluorescein-labeled miRNA mimic Negative Control siRNA #2		5
SMC-5002			10
SMC-5003			20
SMC-2101	AccuTarget™ miRNA inhibitor Negative Control #1	Bio-RP	5
SMC-2102			10
SMC-2103			20
SMC-3101	AccuTarget™ miRNA inhibitor Negative Control #2		5
SMC-3102			10
SMC-3103			20
SMC-4101	AccuTarget™ Fluorescein-labeled miRNA inhibitor Negative Control siRNA #1	HPLC	5
SMC-4102			10
SMC-4103			20
SMC-5101	AccuTarget™ Fluorescein-labeled miRNA inhibitor Negative Control siRNA #2		5
SMC-5102			10
SMC-5103			20

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Reference

1. Erhard F, et al. (2014) Widespread context dependency of microRNA-mediated regulation. *Genome research* 24(6):906-919.
2. Carroll AP, Tooney PA, & Cairns MJ (2013) Context-specific microRNA function in developmental complexity. *Journal of Molecular Cell Biology* 5(2):73-84.
3. Gao F-B (2010) Context-dependent functions of specific microRNAs in neuronal development. *Neural Development* 5(1):25.
4. Huang Y, et al. (2011) Biological functions of microRNAs: a review. *Journal of physiology and biochemistry* 67(1):129-139.
5. Jin Y, Chen Z, Liu X, & Zhou X (2013) Evaluating the microRNA targeting sites by luciferase reporter gene assay. *Methods in molecular biology* (Clifton, N.J.) 936:117-127.

03. *In vivo* siRNA Custom Service: SAMiRNA™

<i>In vivo</i> siRNA Custom Service Using SAMiRNA™ Technology	63
SAMiRNA™ Technical Spec	65

In vivo siRNA Custom Service Using SAMiRNA™ Technology

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* NOTICE TO PURCHASER: LIMITED LICENSE

These products are for research purpose only. They are not for diagnostic or therapeutic purpose and are not to be administered to humans.

Use of SAMiRNA™ is covered by U.S. patents Nos. 8,779,114, 8,771,976, and their foreign counterparts or pending patents. The purchase of these products includes a limited, non-transferable immunity from suit under the forgoing patents for using only this amount of product solely for the purchaser's own internal research. No other patent rights to use this product for any other purpose or for commercial purpose, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, are conveyed expressly, by implication, or by estoppels. Further information on purchasing licenses may be obtained by contacting Bioneer at licensing@bioneer.com

SAMiRNA™ is a trademark of Bioneer Corporation.

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■ Genome-wide predesigned mouse, rat and human siRNA

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■ SAMiRNA™ nanoparticle synthesis for *in vivo* siRNA delivery and efficacy/biodistribution tests using animal models

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■ *In vitro* screen/validation of siRNAs of choice lead optimization service

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In vivo siRNA Custom Service Using SAMiRNA™ Technology

■ Quantitation of functional knockdown of target gene by Real-Time qPCR (AccuPower® qPCR Array Service)

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