
C DNA/RNA Amplification

- 01. DNA Amplification
- 02. RNA Amplification
- 03. Real-Time PCR
- 04. Customized PCR

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Selection Guide of AccuPower® PreMix Series

○ DNA Amplification

| Products | Applications | | | | | | | | |
|------------------------------------|--------------|--------------|---------------------------------|----------------------|--------------------|-------------------|----------------|-------------|---------------|
| | Standard PCR | HotStart PCR | Prevent Carryover Contamination | PCR for Gene Cloning | PCR for TA Cloning | High Fidelity PCR | Long Range PCR | GC Rich PCR | Multiplex PCR |
| AccuPower® PCR | √ | | | | √ | | | | |
| AccuPower® PCR (with UDG) | √ | | √ | | √ | | | | |
| AccuPower® Taq PCR | √ | | | | √ | | | | |
| AccuPower® HotStart PCR | √ | √ | | | √ | | | | |
| AccuPower® HotStart PCR (with UDG) | √ | √ | √ | | √ | | | | |
| AccuPower® GoldHotStart Taq PCR | √ | √ | | | √ | | | √ | |
| AccuPower® PyroHotStart Taq PCR | √ | √ | | | √ | | | | |
| AccuPower® HotStart Pfu PCR | √ | √ | | √ | | √ | | | |
| AccuPower® ProFi Taq PCR | √ | | | √ | √ | √ | √ | | |
| AccuPower® Pfu PCR | √ | | | √ | | √ | | | |
| AccuPower® Multiplex PCR | √ | √ | | | √ | | | | √ |
| AccuPower® Gold Multiplex PCR | √ | √ | | | √ | | | | √ |

Selection Guide

○ RNA Amplification

cDNA Synthesis Kits

| Products | Applications | | | |
|-----------------------------------------------------------------------------------|--------------|-----------------|--------------------------------|-------------------|
| | Standard RT | High Efficiency | RT of Secondary Structured RNA | RT of Long kb RNA |
| <i>AccuPower</i> [®] RT | √ | | | |
| <i>AccuPower</i> [®] <i>CycleScript</i> [™] RT | √ | √ | √ | |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT | √ | | √ | |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] Cycle RT | √ | √ | √ | |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT (RNase H Minus) | √ | √ | √ | √ |

One-step RT-PCR Kits

| Products | Applications | | | | | |
|----------------------------------------------------------------------------------------------|-----------------|------------------------------------|---------------------------------------|------------------|---------------------------------|-----------------------|
| | Standard RT-PCR | RT-PCR of Secondary Structured RNA | High Specificity & Sensitivity RT-PCR | Multiplex RT-PCR | Prevent Carryover Contamination | RT-PCR of Long kb RNA |
| <i>AccuPower</i> [®] RT-PCR | √ | | | | | |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT-PCR PreMix | √ | √ | | | | |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT-PCR PreMix (RNase H Minus) | √ | √ | | | | √ |
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-PCR | √ | √ | √ | | | |
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-PCR (with UDG) | √ | √ | √ | | √ | |
| <i>AccuPower</i> [®] <i>RocketPlex</i> RT-PCR | √ | √ | | √ | | |

○ Real-Time PCR

| Products | Applications | | | | |
|---------------------------------------------------------------------------------|--------------|---------------|--------------------------|-------------------------|---------------------------------|
| | HotStart RT | HotStart qPCR | qPCR (dsDNA Binding dye) | qPCR (Hydrolysis Probe) | Prevent Carryover Contamination |
| <i>AccuPower</i> [®] <i>GreenStar</i> [™] qPCR | | √ | √ | | |
| <i>AccuPower</i> [®] <i>DualStar</i> [™] qPCR | | √ | | √ | |
| <i>AccuPower</i> [®] Plus <i>DualStar</i> [™] qPCR | | √ | | √ | |
| <i>AccuPower</i> [®] Plus <i>DualStar</i> [™] qPCR (with UDG) | | √ | | √ | √ |
| <i>AccuPower</i> [®] <i>GreenStar</i> [™] RT-qPCR | √ | √ | √ | | |
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-qPCR | √ | √ | | √ | |

Selection Guide

Specifications of AccuPower® PreMix Series

○ PCR & One-step RT-PCR

| Product | Product Size | 5' → 3' Exonuclease | 3' → 5' Exonuclease | Specificity | Fidelity | GC-rich | Leaves 3'-A |
|--------------------------------------|--------------|---------------------|---------------------|-------------|----------|---------|-------------|
| PCR | ≤ 10 kb | No | No | ●●● | ●●● | ●●● | Yes |
| PCR (with UDG) | ≤ 10 kb | No | No | ●●● | ●●● | ●●● | Yes |
| Taq PCR | ≤ 10 kb | Yes | No | ●●● | ●●● | ●●● | Yes |
| Pfu PCR | ≤ 15 kb | No | Yes | ●●● | ●●●●● | ●●● | No |
| ProFi Taq PCR | ≤ 30 kb | Yes | Yes | ●●●● | ●●●● | ●●● | Yes |
| HotStart PCR | ≤ 12 kb | No | No | ●●●● | ●●● | ●●●● | Yes |
| HotStart PCR (with UDG) | ≤ 12 kb | No | No | ●●●● | ●●● | ●●●● | Yes |
| GoldHotStart Taq PCR | ≤ 5 kb | Yes | No | ●●●● | ●●● | ●●●●● | Yes |
| PyroHotStart Taq PCR | ≤ 5 kb | Yes | No | ●●●● | ●●● | ●●●● | Yes |
| HotStart Pfu PCR | ≤ 5 kb | No | Yes | ●●●● | ●●●●● | ●●● | No |
| Multiplex PCR | ≤ 1 kb | No | No | ●●●● | ●●● | ●●● | Yes |
| Gold Multiplex PCR | ≤ 1 kb | No | No | ●●●● | ●●● | ●●● | Yes |
| RT-PCR | ≤ 5 kb | No | No | - | - | ●●● | Yes |
| RocketScript™ RT-PCR | ≤ 6 kb | No | No | - | - | ●●● | Yes |
| RocketScript™ RT-PCR (RNase H Minus) | ≤ 12.5 kb | Yes | No | - | - | ●●● | Yes |
| Dual-HotStart™ RT-PCR | ≤ 3kb | Yes | No | ●●●●● | - | ●●●● | Yes |
| Dual-HotStart™ RT-PCR (with UDG) | ≤ 3kb | Yes | No | ●●●●● | - | ●●●● | Yes |
| RocketPlex RT-PCR | ≤ 1 kb | No | No | ●●●● | - | ●●● | Yes |

○ Reverse Transcription

| Product | Product Size | RNase H Activity | DNase Activity | RNase Activity | GC-rich |
|----------------------------------|--------------|------------------|----------------|----------------|---------|
| RT | ≤ 9 kb | Yes | No | No | ●●● |
| CycleScript™ RT | ≤ 9 kb | Yes | No | No | ●●● |
| RocketScript™ RT | ≤ 10 kb | Yes | No | No | ●●● |
| RocketScript™ Cycle RT | ≤ 10 kb | Yes | No | No | ●●● |
| RocketScript™ RT (RNase H Minus) | ≤ 12.5 kb | No | No | No | ●●● |

○ qPCR & One-step RT-qPCR

| Product | 5' → 3' Exonuclease | 3' → 5' Exonuclease | Specificity | GC-rich | Leaves 3'-A |
|--------------------------------|---------------------|---------------------|-------------|---------|-------------|
| GreenStar™ qPCR | No | No | ●●●● | ●●●● | Yes |
| DualStar™ qPCR | Yes | No | ●●●● | ●●●● | Yes |
| Plus DualStar™ qPCR | Yes | No | ●●●●● | ●●●● | Yes |
| Plus DualStar™ qPCR (with UDG) | Yes | No | ●●●●● | ●●●● | Yes |
| GreenStar™ RT-qPCR | Yes | No | ●●●● | ●●●● | Yes |
| Dual-HotStart™ RT-qPCR | Yes | No | ●●●●● | ●●●● | Yes |

Selection Guide

Selection Guide of Enzymes

| Products | | Applications | | | | | | | | |
|-------------------------|------------------------------------|--------------|-----------------|----------------------|--------------------|-------------------|----------------|-------------|---------------|------------|
| | | PCR | | | | | | | | |
| | | Standard PCR | HotStart PCR | PCR for Gene Cloning | PCR for TA Cloning | High fidelity PCR | Long Range PCR | GC Rich PCR | Multiplex PCR | Microbiome |
| DNA Polymerase | <i>Top</i> | √ | | | √ | | | | | |
| | <i>Taq</i> | √ | | | √ | | | | | |
| | <i>Pfu</i> | √ | | √ | | √ | | | | |
| | <i>ProFi Taq</i> | √ | | √ | √ | √ | √ | | | |
| | MicroBiome Assay <i>Taq</i> | √ | | | √ | | | | | √ |
| HotStart DNA Polymerase | <i>Top</i> | √ | √ | | √ | | | √ | √ | |
| | <i>Taq</i> | √ | √ | | √ | | | √ | √ | |
| Products | | RT | | | | | | | | |
| | | Standard RT | High Efficiency | Cyclic RT | Long Range RT | | | | | |
| Reverse Transcriptase | <i>M-MLV</i> | √ | | | | | | | | |
| | <i>CycleScript™</i> | √ | | √ | | √ | | | | |
| | <i>RocketScript™</i> | √ | | √ | | √ | | | | |
| | <i>RocketScript™</i> RNase H Minus | √ | | √ | | √ | | | √ | |

Specifications of Enzymes

| Product | Product Size | 5' → 3' Exonuclease | 3' → 5' Exonuclease | Specificity | Fidelity | GC-rich | Leaves 3'-A |
|--------------------------------------------|--------------|---------------------|---------------------|-------------|----------|---------|-------------|
| <i>Top</i> DNA Polymerase | ≤ 10 kb | No | No | ●●● | ●●● | ●●● | Yes |
| <i>Taq</i> DNA Polymerase | ≤ 10 kb | Yes | No | ●●● | ●●● | ●●● | Yes |
| HotStart DNA Polymerase | ≤ 12 kb | No | No | ●●●● | ●●● | ●●●● | Yes |
| HotStart <i>Taq</i> DNA Polymerase | ≤ 12 kb | Yes | No | ●●●● | ●●● | ●●●● | Yes |
| <i>ProFi Taq</i> DNA Polymerase | ≤ 30 kb | Yes | Yes | ●●●● | ●●●● | ●●●● | Yes |
| <i>Pfu</i> DNA Polymerase | ≤ 20 kb | No | Yes | ●●● | ●●●●● | ●● | No |
| MicroBiome Assay <i>Taq</i> DNA Polymerase | ≤ 8 kb | Yes | No | ●●● | ●●● | ●●● | Yes |
| <i>M-MLV</i> Reverse Transcriptase | ≤ 9 kb | - | - | - | - | ●●● | - |
| <i>CycleScript™</i> Reverse Transcriptase | ≤ 9 kb | - | - | - | - | ●●● | - |
| <i>RocketScript™</i> Reverse Transcriptase | ≤ 10 kb | - | - | - | - | ●●● | - |
| <i>RocketScript™</i> RTase RNase H Minus | ≤ 12.5 kb | - | - | - | - | ●●● | - |

01. DNA Amplification

Standard PCR Kit

| | |
|------------------------------------------------------------------|----|
| <i>AccuPower</i> [®] PCR PreMix & Master Mix | 74 |
| <i>AccuPower</i> [®] PCR PreMix (with UDG) | 76 |
| <i>AccuPower</i> [®] <i>Taq</i> PCR PreMix & Master Mix | 78 |

HotStart PCR Kit

| | |
|-------------------------------------------------------------------------------|----|
| <i>AccuPower</i> [®] HotStart PCR PreMix | 80 |
| <i>AccuPower</i> [®] HotStart PCR PreMix (with UDG) | 82 |
| <i>AccuPower</i> [®] <i>GoldHotStart Taq</i> PCR PreMix & Master Mix | 84 |
| <i>AccuPower</i> [®] <i>PyroHotStart Taq</i> PCR PreMix | 86 |
| <i>AccuPower</i> [®] HotStart <i>Pfu</i> PCR PreMix | 88 |

High Fidelity PCR Kit

| | |
|------------------------------------------------------------------|----|
| <i>AccuPower</i> [®] <i>Pfu</i> PCR PreMix & Master Mix | 90 |
|------------------------------------------------------------------|----|

Long & High Fidelity PCR Kit

| | |
|------------------------------------------------------------------------|----|
| <i>AccuPower</i> [®] <i>ProFi Taq</i> PCR PreMix & Master Mix | 92 |
|------------------------------------------------------------------------|----|

Multiplex PCR Kit

| | |
|-----------------------------------------------------------------|----|
| <i>AccuPower</i> [®] Multiplex PCR PreMix & Master Mix | 94 |
| <i>AccuPower</i> [®] Gold Multiplex PCR PreMix | 96 |

DNA Polymerase

| | |
|--------------------------------------------|-----|
| <i>Top</i> DNA Polymerase | 98 |
| <i>Taq</i> DNA Polymerase | 100 |
| <i>Pfu</i> DNA Polymerase | 102 |
| <i>ProFi Taq</i> DNA Polymerase | 104 |
| HotStart DNA Polymerase | 106 |
| HotStart <i>Taq</i> DNA Polymerase | 108 |
| MicroBiome Assay <i>Taq</i> DNA Polymerase | 110 |

Conventional PCR Instrument

AllInOneCycler[™] → Go to M. Instruments & Devices

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- Conventional PCR
- Primer extension
- TA cloning
- Gene sequencing

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- Enzyme: *Top DNA polymerase*
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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-20°C

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
Figure 1. Comparison of sensitivity test for PCR PreMix and other companies' products using serially diluted human gDNA.

Amplification of the human insulin receptor gene.

I: AccuPower® PCR PreMix

II: A company's *Taq* DNA polymerase

III: B company's *Taq* DNA polymerase

IV: C company's PCR PreMix

Lane 1: Human gDNA 10 ng

Lane 2: Human gDNA 1 ng

Lane 3: Human gDNA 100 pg

Lane 4: Human gDNA 10 pg

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

AccuPower® PCR PreMix & Master Mix

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
Figure 2. Comparison of processivity test for AccuPower® PCR PreMix

and other companies' products using lambda DNA.

Rxn. condition: 95°C 5 min, [94°C 30 sec, 57°C 30 sec] (30 cycles), 72°C 5 min.

A: AccuPower® PCR PreMix

B: S company's Taq premix type

C: T company's Taq premix type

Lane 1: 1 kb fragment of Lambda DNA

Lane 2: 2 kb fragment of Lambda DNA

Lane 3: 3 kb fragment of Lambda DNA

Lane 4: 4 kb fragment of Lambda DNA

M: 100 bp Plus DNA Ladder (Cat. No. D-1035, Bioneer)


Figure 3. Comparison of thermostability of AccuPower® PCR PreMix.

AccuPower® PCR PreMix is incubated at 95°C with various time.

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- Conventional PCR
- Primer extension
- TA cloning
- Gene sequencing
- Molecular diagnosis

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- Enzyme: *Taq* DNA polymerase
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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-20°C


Figure 1. Prevention of carryover contamination.

AccuPower® PCR PreMix (with UDG)

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
Figure 2. Comparison of sensitivity between *AccuPower*® PCR PreMix and *AccuPower*® PCR PreMix (with UDG).

Sensitivity test was operated using serial diluted human genomic DNA. Reaction mixture was incubated at 37°C for 2 min followed by 95°C for 5 min, 30 cycles of 20 sec at 95°C, 20 sec at 55°C, 30 sec at 72°C. *AccuPower*® PCR PreMix (with UDG) contains dUTP besides dATP, dGTP, dCTP and dTTP. The amount of DNA used to test is represented below.

Lane 1: 10 ng of human genomic DNA

Lane 2: 1 ng of human genomic DNA

Lane 3: 100 pg of human genomic DNA

Lane 4: 10 pg of human genomic DNA

M: *AccuLadder*™100 bp DNA Size Marker (Cat. No. D-1030-1, Bioneer)


Figure 3. Efficiency of uracil DNA glycosylase using PCR products (including uracil base).

Efficiency test of uracil DNA glycosylase was operated using serial diluted PCR products including uracil base. *AccuPower*® PCR PreMix was also tested for negative control. Reaction mixture was incubated at 37°C for 2 min, followed by 95°C for 5 min, 30 cycles of 20 sec at 95°C, 20 sec at 55°C, 30 sec at 72°C.

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- Conventional PCR
- Primer extension
- TA cloning
- Gene sequencing

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- Enzyme: Taq DNA polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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-20°C

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
Figure 1. Comparison of PCR amplification quality between AccuPower® Taq PCR PreMix from Bioneer and other companies' Taq PCR master mix.

Target Gene: Human insulin receptor gene.

The cycling conditions for AccuPower® Taq PCR PreMix were 95°C for 5 min, 30 cycles of 20 sec at 95°C, 20 sec at 55°C, and 30 sec at 72°C. PCR reactions using other companies' PCR master mix were performed according to each company's protocol.

Lane 1: 10 ng Human genomic DNA

Lane 2: 1 ng Human genomic DNA

Lane 3: 100 pg Human genomic DNA

Lane 4: 10 pg Human genomic DNA

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

AccuPower® Taq PCR PreMix & Master Mix

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
Figure 2. Comparison of amplification quality between AccuPower® Taq PCR PreMix and competitors' Taq PCR master mix.

Target Gene: IRGC (Immunity-related GTPase family, cinema).
 Bioneer reaction mixture was followed by 95°C for 5 min, 35 cycles of 20 sec at 95°C, 20 sec at 55°C, 30 sec at 72°C. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.
 Lane 1: 10 ng human genomic DNA
 Lane 2: 1 ng human genomic DNA
 Lane 3: 100 pg human genomic DNA
 Lane 4: 10 pg human genomic DNA
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 3. Comparison of long kb amplification between AccuPower® Taq PCR PreMix and competitors' Taq PCR master mix.

Bioneer reaction mixture was followed by 95°C for 5 min, 30 cycles of 20 sec at 95°C, 20 sec at 65°C, 8 min at 68°C. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.
 Lane 1: 3 kb fragment (Human tumor protein p53 gene)
 Lane 2: 4 kb fragment (Human beta globin region)
 Lane 3: 4.5 kb fragment (Human DNA cross-link repair 1A gene)
 Lane 4: 8 kb fragment (Human hemoglobin epsilon 1 gene)
 M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | |
|---------|-------------------------------|--------------------------------------------------|-------------------|
| K-2601 | AccuPower® Taq PCR PreMix | 0.2 ml thin-wall 8-tube strips with attached cap | 20 µl/rxn |
| K-2603 | | | 50 µl/rxn |
| K-2602 | | | 20 µl/rxn |
| K-2604 | | | 50 µl/rxn |
| K-2609 | AccuPower® Taq PCR Master Mix | 2.5 ml of 2X Mater mix solution | 1.25 ml x 2 tubes |
| K-2610 | | 25 ml of 2X Mater mix solution | 12.5 ml x 2 tubes |

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- Enzyme: *Top* DNA polymerase
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 12 kb

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-20°C


Figure 1. Enzyme-mediated HotStart PCR

AccuPower® HotStart PCR PreMix

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
Figure 2. Specificity comparison between AccuPower® PCR PreMix and AccuPower® HotStart PCR PreMix.

- Lane 1: P75/P73 primer set (139 bp)
- Lane 2: P55/P53 primer set (211 bp)
- Lane 3: P55/P63 primer set (447 bp)
- Lane 4: P75/P83 primer set (618 bp)
- Lane 5: P55/P73 primer set (1,082 bp)
- Lane 6: P65/P83 primer set (1,296 bp)
- Lane 7: P55/P83 primer set (1,561 bp)


Figure 3. Comparison of PCR amplification sensitivity between AccuPower® HotStart PCR PreMix from Bioneer and other suppliers' HotStart PCR kit.

- Lane 1: Human gDNA 10 ng
- Lane 2: Human gDNA 1 ng
- Lane 3: Human gDNA 100 pg
- Lane 4: Human gDNA 10 pg
- M: 100 bp Plus DNA ladder (Cat. No. D-1035, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|--------------------------------|--------------------------------------------------|-----------|-----------|
| K-5050 | AccuPower® HotStart PCR PreMix | 0.2 ml thin-wall 8-tube strips with attached cap | 96 tubes | 20 µl/rxn |
| K-5052 | | | | 50 µl/rxn |
| K-5051 | | | 480 tubes | 20 µl/rxn |
| K-5057 | | | | 50 µl/rxn |
| K-5053 | | 0.5 ml thin-wall tubes with attached cap | 100 tubes | 50 µl/rxn |

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- Enzyme: *Taq* DNA polymerase
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 12 kb

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-20°C


Figure 1. Prevention of carryover contamination.

AccuPower® HotStart PCR PreMix (with UDG)

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
Figure 2. Comparison of specificity between *AccuPower*® HotStart PCR PreMix and HotStart PCR PreMix (with UDG).

Specificity test was operated using 7 pairs of primers targeting P53 gene. Reaction mixture was incubated at 37°C for 2 min followed by 95°C for 5 min, 30 cycles of 20 sec at 95°C, 40 sec at 55°C, 1 min at 72°C.

The amount of DNA (human) used to test is 10 ng.

Lane 1: 139 bp

Lane 2: 211 bp

Lane 3: 447 bp

Lane 4: 618 bp

Lane 5: 1,082 bp

Lane 6: 1,296 bp

Lane 7: 1,561 bp

M: *AccuLadder*™100 bp DNA Size Marker (Cat. No. D-1030-1, Bioneer)


Figure 3. Efficiency of uracil DNA glycosylase using PCR product (including uracil base).

Efficiency test of uracil DNA glycosylase was operated using serial diluted PCR products including uracil base. *AccuPower*® HotStart PCR PreMix was also tested for negative control. Reaction mixture was incubated at 37°C for 2 min followed by 95°C for 5 min, 30 cycles of 20 sec at 95°C, 20 sec at 55°C, 30 sec at 72°C.

Lane 1: 10¹¹ copy

Lane 2: 10¹⁰ copy

Lane 3: 10⁹ copy

Lane 4: 10⁸ copy

Lane 5: 10⁷ copy

Lane 6: 10⁶ copy

Lane 7: 10⁵ copy

Lane 8: 10⁴ copy

Lane N: No template control

M: *AccuLadder*™100 bp DNA Size Marker (Cat. No. D-1030-1, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | |
|----------|---------------------------------------------------|-----------------------------------------------------|-----------|
| K-5050-1 | <i>AccuPower</i> ® HotStart PCR PreMix (with UDG) | 0.2 ml thin-wall tubes with attached cap, 20 µl/rxn | 96 tubes |
| K-5051-1 | | | 480 tubes |

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- Enzyme: GoldHotStart Taq DNA polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 5 kb (Human)

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-20°C

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
Figure 1. Comparison of PCR amplification efficiency between AccuPower® GoldHotStart Taq PCR PreMix from Bioneer and other suppliers' HotStart PCR master mix.

Target: Human insulin receptor gene.

The cycling conditions for AccuPower® GoldHotStart Taq PCR PreMix were 95°C for 5 min, 30 cycles of 95°C for 30 sec, 55°C for 30 sec and 72°C for 30 sec. PCR reactions using other suppliers' PCR master mix were performed according to each suppliers' protocol.

Lane 1: 10 ng of human gDNA

Lane 2: 1 ng of human gDNA

Lane 3: 100 pg of human gDNA

Lane 4: 10 pg of human gDNA

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

AccuPower® GoldHotStart Taq PCR PreMix & Master Mix

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
Figure 2. Comparison of PCR amplification specificity between *AccuPower® GoldHotStart Taq* PCR PreMix from Bioneer and other suppliers' HotStart PCR master mix.

Target: ApoE gene (The PCR product size is 268 bp). The cycling conditions for *AccuPower® GoldHotStart Taq* PCR PreMix were 95°C for 5 min, 35 cycles of 95°C for 30 sec, 57°C for 30 sec and 72°C for 30 sec. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.

Lane 1,2: 100 ng human genomic DNA

Lane 3,4: 10 ng human genomic DNA

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 3. Comparison of PCR amplification specificity between *AccuPower® GoldHotStart Taq* PCR PreMix from Bioneer and other suppliers' HotStart PCR master mix.

The cycling conditions for *AccuPower® GoldHotStart Taq* PCR PreMix were 95°C for 5 min, 35 cycles of 95°C for 30 sec, 65°C for 30 sec and 72°C for 30 sec. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.

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- Enzyme: Taq DNA polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 5 kb (Human)

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-20℃


Figure 1. Enzyme-mediated HotStart PCR

AccuPower® PyroHotStart Taq PCR PreMix

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
Figure 2. Comparison of PCR amplification specificity between *AccuPower® PyroHotStart Taq* PCR PreMix from Bioneer and other suppliers' hotstart PCR master mix.

PCR reactions were performed according to each supplier's protocol.

Target: Human PrP gene

The PrP gene was amplified from human genomic DNA with two different primer sets, separately. This data shows that *AccuPower® PyroHotStart Taq* PCR PreMix has higher amplification efficiency and specificity than other suppliers' HotStart PCR master mix.

Lane 1: 100 ng DNA PrP primer set (500 bp)

Lane 2: 10 ng DNA PrP primer set (500 bp)

Lane 3: 100 ng DNA PrP primer set (705 bp)

Lane 4: 10 ng DNA PrP primer set (705 bp)

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 3. Comparison of PCR amplification specificity between *AccuPower® PyroHotStart Taq* PCR PreMix from Bioneer and other suppliers' HotStart PCR master mix.

The ApoE gene was amplified from 100 ng of human genomic DNA (The PCR product size is 268 bp). This data shows that *AccuPower® PyroHotStart Taq* PCR PreMix has higher amplification efficiency and specificity than other suppliers' HotStart PCR master mix.

Lane 1: *AccuPower® PyroHotStart Taq* PCR PreMix

Lane 2: Supplier I HotStart *Taq* PCR premix

Lane 3: Supplier S HotStart *Taq* PCR master mix

Lane 4: Supplier T HotStart *Taq* PCR master mix

Lane 5: Supplier Q HotStart *Taq* PCR master mix

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | |
|---------|----------------------------------------|--------------------------------------------------|-----------|
| K-2611 | AccuPower® PyroHotStart Taq PCR PreMix | 96 tubes | 20 µl/rxn |
| K-2613 | | | 50 µl/rxn |
| K-2612 | | 480 tubes | 20 µl/rxn |
| K-2614 | | | 50 µl/rxn |
| | | 0.2 ml thin-wall 8-tube strips with attached cap | |


Figure 4. *AccuPower® PyroHotStart Taq* PCR PreMix has high amplification efficiency and specificity.

Specificity test was performed using 7 different sets of primers targeting the P53 gene. 10 ng of human genomic DNA was used for each PCR reaction. The cycling conditions were 95°C for 5 min, 30 cycles of 95°C for 20 sec, 55°C for 40 sec, and 72°C for 1 min, and 72°C for 5 min for final extension.

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- High fidelity amplification
- Gene cloning with blunt ends
- Site-directed mutagenesis
- High specificity PCR
- cDNA template PCR

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- Enzyme: Pfu DNA polymerase
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: Yes
- 3' - A overhang: No
- Fragment size: ~ 5 kb

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-20°C

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
Figure 1. AccuPower® HotStart Pfu PCR PreMix shows enhanced specificity compared to competitors.

Specificity test was performed using 7 different sets of primers targeting the p53 gene. 10 ng of human genomic DNA was used for each PCR reaction. The cycling conditions were 95°C for 5 min, 32 cycles of 95°C for 30 sec, 62°C for 40 sec, and 72°C for 1 min 30 sec, and 72°C for 5 min for final extension.

- Lane 1: P75/73 primer set (139 bp)
- Lane 2: P55/53 primer set (211 bp)
- Lane 3: P55/63 primer set (447 bp)
- Lane 4: P75/83 primer set (618 bp)
- Lane 5: P55/73 primer set (1,082 bp)
- Lane 6: P65/83 primer set (1,296 bp)
- Lane 7: P55/83 primer set (1,561 bp)
- M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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
Figure 2. AccuPower® HotStart Pfu PCR PreMix has high amplification efficiency.

Bioneer reaction mixture was followed by 95°C for 5 min, 35 cycles of 95°C for 20 sec, 65°C for 20 sec, and 68°C for 15 min, and 68°C for 5 min for final extension.

Lane 1: 2 kb fragment

Lane 2: 2.5 kb fragment

Lane 3: 3 kb fragment

Lane 4: 4.5 kb fragment

Lane 5: 5 kb fragment

M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)


Figure 3. Comparison of PCR amplification efficiency between AccuPower® HotStart Pfu PCR PreMix from Bioneer and other suppliers' PCR master mix

Target: human insulin receptor gene.

The cycling conditions for AccuPower® HotStart Pfu PCR PreMix were 95°C for 5 min, 30 cycles of 95°C for 30 sec, 55°C for 30 sec and 72°C for 2 min. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.

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- Gene synthesis
- Gene cloning
- Conventional PCR
- Primer extension
- Site directed mutagenesis
- High fidelity가 요구되는 실험

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- Enzyme: Pfu DNA polymerase
- 5' → 3' exonuclease activity: No
- 3' → 5' exonuclease activity: Yes
- 3' - A overhang: No
- Fragment size: ~ 15 kb

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-20°C

AccuPower® Pfu PCR PreMix & Master Mix

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
Figure 1. Template range & sensitivity of AccuPower® Pfu PCR PreMix for human DNA template.

Test of working range & sensitivity of AccuPower® Pfu PCR PreMix for human DNA template.

Lane 1: 100 ng Lane 2: 10 ng Lane 3: 1 ng
 Lane 4: 100 pg Lane 5: 10 pg Lane 6: Template negative
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 2. Amplification of lambda DNA of 1 kb to 10 kb with AccuPower® Pfu PCR PreMix.

Lane 1: 1 kb fragment Lane 2: 2 kb fragment
 Lane 3: 3 kb fragment Lane 4: 4 kb fragment
 Lane 5: 5 kb fragment Lane 6: 6 kb fragment
 Lane 7: 7 kb fragment Lane 8: 8 kb fragment
 Lane 9: 9 kb fragment Lane 10: 10 kb fragment
 M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|-------------------------------|--------------------------------------------------|-----------|-----------|
| K-2022 | AccuPower® Pfu PCR PreMix | 0.2 ml thin-wall 8-tube strips with attached cap | 96 tubes | 20 µl/rxn |
| K-2023 | | | | 50 µl/rxn |
| K-2024 | | 480 tubes | 20 µl/rxn | |
| K-2025 | | | 50 µl/rxn | |
| K-2027 | | 0.5 ml thin-wall 8-tube strips with attached cap | 100 tubes | 50 µl/rxn |
| K-2026 | AccuPower® Pfu PCR Master Mix | 1 ml of 2X master mix solution | | |

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- Primer extension
- Long range amplification from genomic DNA
- High amplification efficiency
- Excellent performance on difficult templates
- Amplification of low-copy targets
- High yield and high sensitivity PCR

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- Enzyme: ProFi Taq DNA polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: Yes
- 3' - A overhang: Yes
- Fragment size: ~ 30 kb

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-20°C

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
Figure 1. Comparison of PCR amplification efficiency between AccuPower® ProFi Taq PCR PreMix from Bioneer and other suppliers' PCR master mix.

cDNA synthesized from 10-fold serial-diluted human total RNA from 10 ng to 10 pg using AccuPower® RocketScript™ Cycle RT PreMix (Cat. No. K-2201, Bioneer) was used as a template for PCR amplification. The cycling conditions for AccuPower® ProFi Taq PCR PreMix were 95°C for 5 min, 33 cycles of 95°C for 20 sec, 55°C for 20 sec and 72°C for 30 sec. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.

Target: human GAPDH gene

Lane 1: 10 ng of human total cDNA

Lane 2: 1 ng of human total cDNA

Lane 3: 100 pg of human total cDNA

Lane 4: 10 pg of human total cDNA

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

AccuPower® ProFi Taq PCR PreMix & Master Mix

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
Figure 2. Comparison of PCR amplification sensitivity between *AccuPower® ProFi Taq* PCR PreMix from Bioneer and other suppliers' PCR master mix.

The cycling conditions for *AccuPower® ProFi Taq* PCR PreMix were 95°C for 5 min, 30 cycles of 95°C for 20 sec, 65°C for 20 sec and 68°C for 4 min. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol.

- Lane 1: 2 kb fragment (human tumor protein p53 gene)
- Lane 2: 3 kb fragment (human tumor protein p53 gene)
- Lane 3: 4.5kb fragment (human DNA cross-link repair 1A gene)
- Lane 4: 8 kb fragment (human hemoglobin epsilon 1 gene)
- M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)


Figure 3. Comparison of PCR amplification of long targets between *AccuPower® ProFi Taq* PCR PreMix from Bioneer and other suppliers' PCR master mix.

The cycling conditions for *AccuPower® ProFi Taq* PCR PreMix were 95°C for 5 min, 32 cycles of 95°C for 20 sec, 65°C for 40 sec, and 68°C for 15 min. PCR reactions using other suppliers' PCR master mix were performed according to each supplier's protocol. Human DNA was used as a template for PCR amplification.

- Lane 1: 11 kb fragment
- Lane 2: 13.5 kb fragment
- Lane 3: 17.6 kb fragment
- Lane 4: 21.4 kb fragment
- M1: Lambda DNA/*Hind* III marker (Cat. No. D-1050, Bioneer)
- M2: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|--------------------------------------------|---------------------------------------------------|-----------|-----------|
| K-2631 | <i>AccuPower® ProFi Taq</i> PCR PreMix | 0.2 ml thin-wall 8-tube strips with attached caps | 96 tubes | 20 µl/rxn |
| K-2633 | | | 50 µl/rxn | |
| K-2632 | | 480 tubes | 20 µl/rxn | |
| K-2634 | | | 50 µl/rxn | |
| K-2635 | <i>AccuPower® ProFi Taq</i> PCR Master Mix | 1 ml of 2X Master mix solution | | |

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- Enzyme: HotStart Top DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 1 kb

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-20°C

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
Figure 1. Single PCR and multiplex PCR using AccuPower® Multiplex PCR PreMix.

Each lane from left to right indicates the single and multiplex PCR product using AccuPower® Multiplex PCR PreMix.


a) 10-plex multiplex PCR

b) 20-plex multiplex PCR

M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)

AccuPower® Multiplex PCR PreMix & Master Mix

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10-plex

Figure 2. Comparison of amplification quality between *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit.

10-plex primers were added into *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit. A series of human genomic DNA diluents were tested.

Lane 1: Human genomic DNA 100 ng, Lane 2: Human genomic DNA 10 ng, Lane 3: Human genomic DNA 1 ng

All data were obtained using *MyGenie*™ 96 Gradient Thermal Block (Bioneer).


Supplier Q: Supplier Q's multiplex PCR kit

Supplier S: Supplier S's multiplex PCR kit

Supplier I: Supplier I's *Taq* DNA polymerase for multiplex PCR (0.5 U, added 2 mM MgCl₂)

Rxn. condition: 95°C for 10 min, followed by 35 cycles of 30 sec at 95°C, 30 sec at 65°C, 60 sec at 72°C

M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)

20-plex

Figure 3. Comparison of amplification quality between *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit.

20-plex primers were added into *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit. A series of human genomic DNA diluents were tested (1 ng~100 ng). All data were obtained using *MyGenie*™ 96 Gradient Thermal Block (Bioneer).

Supplier Q: Supplier Q's multiplex PCR kit

Supplier S: Supplier S's multiplex PCR kit

Supplier I: Supplier I's *Taq* DNA polymerase for multiplex PCR (0.5 U, added 2 mM MgCl₂)

Rxn. condition: 95°C for 10 min, followed by 35 cycles of 30 sec at 95°C, 30 sec at 57°C, 60 sec at 72°C

M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)


Figure 4. Comparison of amplification quality using Labchip™ between *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit.

a) Virtual Gel Image. Gel image illustrates data reproducibility of the LabChip™ 90 system.

b) Overlay of expression level using Bioneer's Multiplex PCR PreMix and other company's multiplex PCR kit. The electropherogram displays the data between 10 PCR products yield using 10-plex primer sets to illustrate the amplification efficiency.

c) The graph shows the total concentration of PCR products between *AccuPower*® Multiplex PCR PreMix and other suppliers' multiplex PCR kit.

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| 카탈로그 번호 | 제품명 | 규격 |
|---------|---------------------------------------------|--------------------------------|
| K-2111 | <i>AccuPower</i> ® Multiplex PCR PreMix | 96 tubes |
| K-2112 | | |
| K-2113 | | 480 tubes |
| K-2114 | | |
| K-2120 | <i>AccuPower</i> ® Multiplex PCR Master Mix | 1 ml of 2X Master mix solution |

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| Target | Application |
|------------------|-----------------------------------------------------------------|
| Human and Animal | STR analysis for determining genetic profiles in forensic cases |
| | Molecular diagnostic analysis |
| | Genotyping assay |
| | Qualitative and semi-qualitative gene expression assay |
| | Mutant screening |
| Plant | Transgenic organism analysis |
| | STR analysis |
| | Detection of pathogens/bacteria infection |
| | Qualitative and semi-qualitative gene expression assay |

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- Enzyme: Top DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 1 kb

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-20°C


Figure 1. Enzyme-mediated HotStart PCR

AccuPower® Gold Multiplex PCR PreMix

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
Figure 2. High specificity of *AccuPower*® Gold Multiplex PCR PreMix.

Each lane from left to right represents the progressive number of primer sets (1~20) included in *AccuPower*® Gold Multiplex PCR PreMix reactions.

Rxn. condition: 95°C for 10 min, followed by 30 cycles (a), 35 cycles (b) of 30 sec at 95°C, 30 sec at 57°C, 60 sec at 72°C


Figure 3. High specificity of *AccuPower*® Gold Multiplex PCR PreMix.

Each lane from left to right indicates the single and multiplex PCR product using *AccuPower*® Gold Multiplex PCR PreMix.

a) 10-plex multiplex PCR b) 20-plex multiplex PCR

Lane M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)

Rxn. condition: 95°C for 10 min, followed by 30 cycles (a), 35 cycles (b) of 30 sec at 95°C, 30 sec at 57°C, 60 sec at 72°C


Figure 4. Comparison of amplification quality between *AccuPower*® Gold Multiplex PCR PreMix and other suppliers' multiplex PCR kit.

6-plex a), 10-plex b), 20-plex c) primers were added into *AccuPower*® Gold Multiplex PCR PreMix and other suppliers' master mixture. A series of human genomic DNA diluents were tested (Lane 1: 100 ng, Lane 2: 10 ng, Lane 3: 1 ng). All data were obtained using *MyGene*™ 96 Gradient Thermal Block (Bioneer).

Supplier Q: Multiplex PCR master mix

Supplier S: Multiplex PCR master mix

Supplier I: *Taq* DNA polymerase for multiplex PCR (0.5U), added 2 mM MgCl₂

M: 25/100 bp Mixed DNA Ladder (Cat. No. D-1020, Bioneer)

Rxn. condition: 95°C for 10 min, followed by 35 cycles of 30 sec at 95°C, 30 sec at 60°C, 60 sec at 72°C

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| 카탈로그 번호 | 제품명 | 규격 |
|---------|--------------------------------------|-----------|
| K-2115 | AccuPower® Gold Multiplex PCR PreMix | 96 tubes |
| K-2117 | | |
| K-2116 | | 50 µl/rxn |
| K-2118 | | 480 tubes |
| | | 20 µl/rxn |
| | | 50 µl/rxn |

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- Real-Time quantification of DNA and cDNA targets using dsDNA binding dye
- Gene expression profiling
- Microbial & viral pathogen detection

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- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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- 10X Reaction buffer with (or without) MgCl₂: Tris (pH 9.0), 15 mM MgCl₂, etc
- 1X dilution buffer: 50% glycerol containing 50 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, stabilizers, pH 8.0
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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500 Units (5 U/μl)

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50% glycerol containing 50 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, stabilizers, pH 8.0

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-20°C

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One Unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 min at 72°C.

Top DNA Polymerase

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
Figure 1. Sensitivity test of *Top* DNA polymerase using Human genomic DNA. Each fragment was amplified from a template dilution series (10 ng to 10 fg DNA per reaction) using 1 U of each *Top* DNA polymerase.


Figure 2. Enzyme activity test of *Top* DNA Polymerase. *Top* DNA polymerase was serially diluted and used to amplify 10 ng of each human genomic DNA.

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| 카탈로그 번호 | 제품 규격 | | | | |
|----------|---------------------------|-----------------------------------|-----------------|-------------|-------------------------|
| | <i>Top</i> DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP | 20 mM MgCl ₂ |
| E-3100 | 500 U | 1 ml (with MgCl ₂) | 1 ml | 1 ml | - |
| E-3100-1 | | 1 ml (without MgCl ₂) | | 1 ml | 1 ml |
| E-3100-2 | | 1 ml (with MgCl ₂) | | - | - |
| E-3100-3 | | 1 ml (without MgCl ₂) | | - | 1 ml |
| E-3101 | 2,000 U | 1 ml (with MgCl ₂) | 1 ml x 4 ea | 1 ml x 4 ea | - |
| E-3101-1 | | 1 ml (without MgCl ₂) | | 1 ml x 4 ea | 1 ml x 4 ea |
| E-3101-2 | | 1 ml (with MgCl ₂) | | - | - |
| E-3101-3 | | 1 ml (without MgCl ₂) | | - | 1 ml x 4 ea |

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- Real-Time quantification of DNA and cDNA targets using Dual probe, dsDNA binding dye
- Gene expression profiling
- Microbial & viral pathogen detection

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- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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- 10X reaction buffer with (or without) MgCl₂: Tris-HCl, KCl, 15 mM MgCl₂, pH 9.0
- Dilution buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol, pH 8.0
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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500 units (5 U/μl)

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20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol pH 8.0

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-20°C

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One Unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 min at 72°C.

Taq DNA Polymerase

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
Figure 1. Using 1 Unit of *Taq* DNA Polymerase, the activity of the polymerase was tested on human genomic DNA (A), lambda genomic DNA (B) as template. Each template DNA was serially diluted by ten-folds, with different ranges.

Lane 1: 100 ng Template DNA Lane 2: 10 ng Template DNA
 Lane 3: 1 ng Template DNA Lane 4: 100 pg Template DNA
 Lane 5: 10 pg Template DNA Lane 6: 1 pg Template DNA
 Lane 7: 100 fg Template DNA
 M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)


Figure 2. Amplification of fragments ranging from 5 kb to 8 kb from template Lambda DNA 20 pg with 1 units of *Taq* DNA Polymerase.

Lane 1: 5 kb PCR product Lane 2: 6 kb PCR product
 Lane 3: 7 kb PCR product Lane 4: 8 kb PCR product
 M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품 규격 | | | | |
|----------|---------------------------|-----------------------------------|-----------------|-------------|-------------------------|
| | <i>Taq</i> DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP | 20 mM MgCl ₂ |
| E-2011 | 500 U | 1 ml (with MgCl ₂) | 1 ml | 1 ml | - |
| E-2011-1 | | 1 ml (without MgCl ₂) | | 1 ml | 1 ml |
| E-2011-2 | | 1 ml (with MgCl ₂) | | - | - |
| E-2011-3 | | 1 ml (without MgCl ₂) | | - | 1 ml |
| E-2013 | 2,000 U | 1 ml (with MgCl ₂) | 1 ml x 4 ea | 1 ml x 4 ea | - |
| E-2013-1 | | 1 ml (without MgCl ₂) | | 1 ml x 4 ea | 1 ml x 4 ea |
| E-2013-2 | | 1 ml (with MgCl ₂) | | - | - |
| E-2013-3 | | 1 ml (without MgCl ₂) | | - | 1 ml x 4 ea |

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- Gene synthesis
- PCR or Primer extension requested high fidelity
- Blunt-end PCR Cloning or mutagenesis requested high fidelity

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- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: Yes
- 3' - A overhang: No
- Fragment size: ~ 20 kb

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- 10X reaction buffer: MgSO₄, Tris-HCl, KCl, (NH₄)₂SO₄, Acetylated BSA pH 9.0
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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250 Units (2.5 U/μl)

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50 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, Stabilizers, 50 % Glycerol pH 8.2

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-20°C

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One Unit is defined at the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble material in 30 min at 72°C.

Pfu DNA Polymerase

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
Figure 1. Human DNA was amplified using 1 units of enzyme in 20 µl reaction volume.

Lane 1: 10 ng Lane 2: 1 ng
Lane 3: 100 pg Lane 4: 10 pg
M: 1 kb DNA ladder (Cat. No. D-1040, Bioneer)


Figure 2. Amplification of fragments ranging from 2 kb to 5 kb from template Lambda DNA 100 pg with 1 units of *Pfu* DNA Polymerase.

Lane 1: 2 kb fragment
Lane 2: 3 kb fragment
Lane 3: 4 kb fragment
Lane 4: 5 kb fragment
M: 1 kb DNA ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품 규격 | | | |
|----------|---------------------------|---------------------------------------|-----------------|-------------|
| | <i>Pfu</i> DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP |
| E-2015 | 250 U | 1 ml (with MgSO ₄) | 1 ml | - |
| E-2015-1 | | 1 ml | | 1 ml |
| E-2016 | 1,000 U | 1 ml (with MgSO ₄) x 4 ea | | - |
| E-2016-1 | | 1 ml x 4 ea | | 1 ml x 4 ea |

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- Primer extension
- long-range amplification from genomic DNA
- High amplification efficiency
- Excellent performance on difficult templates
- Amplification of low-copy targets
- High yield and high sensitivity PCR

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- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: Yes
- 3' - A overhang: Yes
- Fragment size: ~ 30 kb

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- 10X reaction buffer with (or without) MgCl₂: Tris-HCl, KCl, 15 mM MgCl₂, pH 9.0
- Dilution buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50 % Glycerol, pH 8.0
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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250 Units (5 U/μl)

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20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol pH 8.0

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-20°C

ProFi Taq DNA Polymerase

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
Figure 1. Comparison of PCR amplification efficiency between *ProFi Taq* DNA Polymerase from Bioneer and other suppliers' DNA polymerase.

The cycling conditions for *ProFi Taq* DNA Polymerase were 95°C for 5 min, 30 cycles of 95°C for 20 sec, 55°C for 20 sec and 72°C for 30 sec. PCR reaction using other suppliers' DNA polymerase were performed according to each supplier's protocol.

Target: human Insulin receptor gene.

Lane 1: 10 ng of human genomic DNA

Lane 2: 1 ng of human genomic DNA

Lane 3: 100 pg of human genomic DNA

Lane 4: 10 pg of human genomic DNA

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 2. Comparison of PCR amplification of long targets between *ProFi Taq* DNA Polymerase from Bioneer and other suppliers' DNA polymerase.

The cycling conditions for *ProFi Taq* DNA Polymerase were 95°C for 5 min, 32 cycles of 95°C for 20 sec and 68°C for 15 min. PCR reactions using other suppliers' DNA polymerase were performed according to each supplier's protocol. Human genomic DNA was used as a template for PCR amplification.

Lane 1: 11 kb fragment Lane 2: 13.5 kb fragment

Lane 3: 17.6 kb fragment Lane 4: 21.4 kb fragment

M1: Lambda DNA/*Hind*III marker (Cat. No. D-1050, Bioneer)

M2: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

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| 카탈로그 번호 | 제품 규격 | | | | |
|---------|---------------------------------|-----------------------------------|-----------------|-------------|-------------------------|
| | <i>ProFi Taq</i> DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP | 20 mM MgCl ₂ |
| E-2201 | 250 U | 1 ml (with MgCl ₂) | 1 ml | 1 ml | - |
| E-2202 | | 1 ml (without MgCl ₂) | | 1 ml | 1 ml |
| E-2203 | | 1 ml (with MgCl ₂) | | - | - |
| E-2204 | | 1 ml (without MgCl ₂) | | - | 1 ml |
| E-2205 | 1,000 U | 1 ml (with MgCl ₂) | 1 ml x 4 ea | 1 ml x 4 ea | - |
| E-2206 | | 1 ml (without MgCl ₂) | | 1 ml x 4 ea | 1 ml x 4 ea |
| E-2207 | | 1 ml (with MgCl ₂) | | - | - |
| E-2208 | | 1 ml (without MgCl ₂) | | - | 1 ml x 4 ea |

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
-20°C

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One Unit is defined as the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 min at 72°C.

HotStart DNA Polymerase

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
Figure 1. Multiplex PCR comparison of genomic DNA using 6 sets of primers and 2 different DNA polymerases

Lane 1: 750 bp fragment Lane 2: 590 bp fragment Lane 3: 450 bp fragment
 Lane 4: 360 bp fragment Lane 5: 260 bp fragment Lane 6: 150 bp fragment
 Lane 7: Multiplex PCR with primers used for Lane 1~6
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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| 카탈로그 번호 | 제품 규격 | | | | |
|----------|-------------------------|---------------------|-----------------|------------|-------------------------|
| | HotStart DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP | 20 mM MgCl ₂ |
| E-3150 | 250 U | 0.5 ml | 0.5 ml | 0.5 ml | 0.5 ml |
| E-3150-1 | | | | - | 0.5 ml |
| E-3151 | 1,000 U | 0.5 ml x 4 ea | | | - |

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- Real-Time quantification of DNA and cDNA targets using dsDNA binding dye
- HotStart PCR
- Multiplex PCR
- Automated PCR

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- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 10 kb

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- 10X reaction buffer: Tris-HCl, KCl, 15 mM MgCl₂, pH 9.0
- Dilution buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol, pH 8.0
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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250 Units (5 U/μl)

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20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol, pH 8.0

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-20°C

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One Unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 min at 72°C.

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
Figure 1. HotStart Taq DNA polymerase의 specificity test. Single and multiple PCR results in human genomic DNA p53 gene amplification.

Lane 1: 139 bp Lane 2: 211 bp
Lane 3: 447 bp Lane 4: 618 bp
Lane 5: 1,082 bp Lane 6: 1,296 bp
Lane 7: 1,561 bp
Lane 8: Multiplex PCR (139 bp, 447 bp, 618 bp)
M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

HotStart *Taq* DNA Polymerase


Figure 2. Specificity comparison between HotStart. *Taq* DNA polymerase and other competitor's products.
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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| 카탈로그 번호 | 제품 규격 | | | |
|----------|------------------------------------|--------------------------------------------|-----------------|---------------|
| | HotStart <i>Taq</i> DNA Polymerase | 10X Reaction Buffer with MgCl ₂ | Dilution Buffer | 10 mM dNTP |
| E-2017 | 250 U | 0.5 ml | | 0.5 ml |
| E-2017-3 | | | | - |
| E-2017-2 | 500 U | 0.5 ml x 2 ea | | |
| E-2017-1 | 1,000 U | 0.5 ml x 4 ea | | 0.5 ml x 4 ea |
| E-2017-4 | | | | - |

MicroBiome Assay *Taq* DNA Polymerase

Prevention of contamination by host cell *E. coli* genomic DNA

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- Routine PCR, multiplex PCR and qPCR
- Allele specific PCR
- 16S and 23S rRNA gene amplification
- Detection of bacteria in samples (e.g. blood)
- DNA labeling reaction & TA-cloning

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- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A overhang: Yes
- Fragment size: ~ 8 kb

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| 카탈로그 번호 | 제품 규격 | | | | |
|---------|-------------------------------|-----------------------------------|-----------------|------------|-------------------------|
| | MBA <i>Taq</i> DNA Polymerase | 10X Reaction Buffer | Dilution Buffer | 10 mM dNTP | 20 mM MgCl ₂ |
| E-3504 | 500 U | 1 ml (without MgCl ₂) | 1 ml | 1 ml | 1 ml |

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- 10X reaction buffer without MgCl₂: 200 mM Tris-HCl, 500 mM KCl Tween 20 0.1%
- Dilution buffer: 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, 50% Glycerol, pH 7.4
- 10 mM dNTPs mix: 2.5 mM of each dNTP
- MgCl₂ Solution: 20 mM

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500 Units (5 U/μl)

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20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol pH 8.0

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-20°C

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One unit is defined at the amount of enzyme that will incorporate 10 nmole of dNTP into acid-insoluble material in 30 min at 72°C.

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
Figure 1. High sensitivity and Low Contamination DNA of MicroBiome Assay *Taq* DNA polymerase.

- Lane 1: *E. coli* gDNA 100 pg
- Lane 2: *E. coli* gDNA 10 pg
- Lane 3: *E. coli* gDNA 1 pg
- Lane 4: *E. coli* gDNA 100 fg
- N: Non-Template

02. RNA Amplification

Standard RT & RT-PCR Kit

| | |
|----------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] RT PreMix & Master Mix | 112 |
| <i>AccuPower</i> [®] RT-PCR PreMix & Master Mix | 114 |

High Efficiency RT Kit

| | |
|--------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>CycleScript</i> [™] RT PreMix & Master Mix | 116 |
|--------------------------------------------------------------------------------------------|-----|

High Sensitivity RT & RT-PCR Kit

| | |
|---------------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT PreMix & Master Mix | 119 |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] Cycle RT PreMix & Master Mix | 121 |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT-PCR PreMix & Master Mix | 124 |

Long & High Efficiency RT Kit

| | |
|----------------------------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT PreMix, RNase H Minus & Master Mix | 126 |
| <i>AccuPower</i> [®] <i>RocketScript</i> [™] RT-PCR PreMix, RNase H Minus & Master Mix | 129 |

Multiplex RT-PCR Kit

| | |
|----------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>RocketPlex</i> RT-PCR PreMix & Master Mix | 131 |
|----------------------------------------------------------------------------------|-----|

HotStart RT & RT-PCR Kit

| | |
|------------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-PCR PreMix | 133 |
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-PCR PreMix (with UDG) | 135 |

Reverse Transcriptase

| | |
|-----------------------------------------------------------------------------|-----|
| <i>M-MLV</i> Reverse Transcriptase | 137 |
| <i>CycleScript</i> [™] Reverse Transcriptase | 138 |
| <i>RocketScript</i> [™] Reverse Transcriptase | 141 |
| <i>RocketScript</i> [™] Reverse Transcriptase, RNase H Minus | 143 |

Conventional PCR Instrument

AllInOneCycler[™] → Go to M. Instruments & Devices

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5' end mapping by primer extension analysis
- Real-Time PCR

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- Enzyme: M-MLV RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 9 kb

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-20°C

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Sensitivity & Reproducibility


Figure 1. Amplification of GAPDH target gene was detected using human total RNA (from 10 ng to 10 pg) with AccuPower® RT PreMix.

Lanes 1-4; 10 ng, 1 ng, and 100 pg, 10 pg of total RNA from HeLa cells, respectively.

Specific Amplification


Figure 2. Specific amplification of 5'-UTR region of HCV with AccuPower® RT PreMix.

Lane 1: negative control
 Lane 2,3: HCV positive serum
 Lane 4: HCV negative serum
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

AccuPower® RT PreMix & Master Mix

Reliability and Reproducibility Test


Figure 3. *AccuPower*® RT PreMix from each lot was tested to confirm reliability and reproducibility.

Human total RNA (panel A: actin, panel B: globin) and hog cholera virus RNA (panel C) were used as template. Following cDNA synthesis, *AccuPower*® PCR PreMix was used to amplify target genes.

Lane 1-4: Reliability test of each lot with *AccuPower*® RT PreMix
M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | | | | |
|----------|----------------------------------|--------------------------------------------------|------------------------------------------|----------------------|----------------------|----------------------|
| K-2041 | <i>AccuPower</i> ® RT PreMix | 0.2 ml thin-wall 8-strip tubes with attached cap | 96 tubes | 20 µl/rxn | | |
| K-2043 | | | | 50 µl/rxn (with dye) | | |
| K-2041-B | | | 480 tubes | | 20 µl/rxn | |
| K-2043-B | | | | | 20 µl/rxn (with dye) | |
| K-2040 | | | 0.5 ml thin-wall tubes with attached cap | 100 tubes | | 20 µl/rxn (with dye) |
| K-2042 | | | | | | 50 µl/rxn |
| K-2261-1 | | thin-wall 96-well flat plate | | | 10 µl/rxn | |
| K-2261-4 | | | | | 20 µl/rxn | |
| K-2261-2 | | | thin-wall 96-well full-skirted plate | | | 10 µl/rxn |
| K-2261-5 | | | | | | 20 µl/rxn |
| K-2261-3 | | | thin-wall 96-well semi-skirted plate | | | 10 µl/rxn |
| K-2261-6 | | | | | | 20 µl/rxn |
| K-2082-1 | | thin-wall 384-well full-skirted plate | | | 5 µl/rxn | |
| K-2082-2 | | | | | 10 µl/rxn | |
| K-2082-3 | | | | | 20 µl/rxn | |
| K-2263 | <i>AccuPower</i> ® RT Master Mix | 1 ml of 2X Master mix solution | | | | |

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- cDNA library construction
- Gene expression analysis

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- Enzyme: M-MLV RTase, Top DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 5 kb

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-20℃

AccuPower® RT-PCR PreMix & Master Mix

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Sensitivity


Figure 1. Sensitivity comparison between AccuPower® RT-PCR PreMix and Supplier I RT-PCR kit.

Each 10^9 copies ~ 10^3 copies of PSTVd used for RT/PCR and the same amount of RT/PCR products used for electrophoresis
 Lane 1: 10^9 copy Lane 2: 10^8 copy Lane 3: 10^7 copy
 Lane 4: 10^6 copy Lane 5: 10^5 copy Lane 6: 10^4 copy
 Lane 7: 10^3 copy Lane 8: NTC
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

Reproducibility


Figure 2. Comparison of reproducibility test for AccuPower® RT-PCR PreMix batch 1, 2, 3 and batch 4 products using serial diluted human total RNA.

Lane 1: 10 ng human total RNA from HeLa cell
 Lane 2: 1 ng human total RNA from HeLa cell
 Lane 3: 100 pg human total RNA from HeLa cell
 Lane 4: 10 pg human total RNA from HeLa cell
 M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)

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| 카탈로그 번호 | 제품명 | 규격 | | |
|----------|--------------------------|--------------------------------------------------|--------------------------------|-----------|
| K-2055 | AccuPower® RT-PCR PreMix | 0.2 ml thin-wall 8-strip tubes with attached cap | 96 tubes | 20 µl/rxn |
| K-2057 | | | 50 µl/rxn | |
| K-2055-B | | | 480 tubes | 20 µl/rxn |
| K-2057-B | | | | 50 µl/rxn |
| K-2056 | | 0.5 ml thin-wall tubes with attached cap | 100 tubes | 50 µl/rxn |
| K-2262-1 | | thin-wall 96-well flat plate | 10 µl/rxn | |
| K-2262-4 | | | 20 µl/rxn | |
| K-2262-2 | | thin-wall 96-well full-skirted plate | 10 µl/rxn | |
| K-2262-5 | | | 20 µl/rxn | |
| K-2262-3 | | thin-wall 96-well semi-skirted plate | 10 µl/rxn | |
| K-2262-6 | | | 20 µl/rxn | |
| K-2084-1 | | thin-wall 384-well full-skirted plate | 5 µl/rxn | |
| K-2084-2 | | | 10 µl/rxn | |
| K-2084-3 | | | 20 µl/rxn | |
| K-2264 | | AccuPower® RT-PCR Master Mix | 1 ml of 2X Master mix solution | |

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
AccuPower® CycleScript™ RT PreMix & Master Mix

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- Sequencing single and double-strand DNA or RNA
- RT-PCR
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- Enzyme: *CycleScript™* RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 9 kb

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-20℃

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
Figure 1. Reaction conditions of CTRT compared with that of conventional RT. The conventional fixed temperature RT reactions at 42°C and 55°C and cyclic temperature RT reactions at 37°C and 50°C using HeLa cell total RNA (800 ng each) were performed 800 bp, 1.5 kb, 2.0 kb, and 2.6 kb fragments of human transferrin receptor gene were amplified with each primers. *AccuPower® CycleScript™* RT PreMix series (dT₂₀ & dN₁₂) preferably shows good results at CTRT reaction condition.


Figure 2. Gene amplification test of transferrin receptor compared with other companies.

The reaction condition was performed according to each manufacturer's recommendation. All cDNAs were reacted with *AccuPower®* PCR PreMix (Cat. No. K-2012, Bioneer).

Lane 1~5: *AccuPower® CycleScript™* RT PreMix (dT₂₀) incubated at 55°C for 1 hr
 Lane 6~10: *AccuPower® CycleScript™* RT PreMix (dN₁₂) incubated at 55°C for 1 hr
 Lane 11~15: C company's RT product including dT primer incubated at 42°C for 1 hr
 Lane 16~20: C company's RT product including random primer incubated at 42°C for 1 hr

Lane 21~25: I company's RT product including dT primer incubated at 45°C for 1 hr
 Lane 26~30: I company's RT product including random primer incubated at 45°C for 1 hr

M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)


Figure 3. RT reaction at various temperature and short reaction time.

Rxn. condition: conventional 42°C 1 hr, 55°C 1hr & cyclic reaction 1: (37°C 2 min/50°C 3 min) X 12, cyclic reaction 2: (37°C 1 min/ 47°C 3 min/ 55°C 1 min) X 12. This product shows thermal stability.

Target: Human GAPDH, human β-actin

Lane 1, 5, 9, 13, 17, 21, 25 & 29: HeLa cell total RNA 100 ng

Lane 2, 6, 10, 14, 18, 22, 26 & 30: HeLa cell total RNA 10 ng

Lane 3, 7, 11, 15, 19, 23, 27 & 31: HeLa cell total RNA 1 ng

Lane 4, 8, 12, 16, 20, 24, 28 & 32: HeLa cell total RNA 100 pg

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer)


Figure 4. Amplification results at various reaction times.

CycleScript™ RT PreMix (dT₂₀) was performed cyclic temperature RT according to (37°C 2 min/50°C 3 min); X12, X4, X2 cycles. 10 min reaction is enough. HeLa cell total RNA templates were serially diluted such as 100 ng, 10 ng, 1 ng, 100 pg.

10 min: 2 times of 2 minutes at 37°C and 3 minutes at 50°C

20 min: 4 times of 2 minutes at 37°C and 3 minutes at 50°C

60 min: 12 times of 2 minutes at 37°C and 3 minutes at 50°C

Lane 1, 5, 9, 13, 17, 21: 100 ng of HeLa cell total RNA used for reaction

Lane 2, 6, 10, 14, 18, 22: 10 ng of HeLa cell total RNA used for reaction

Lane 3, 7, 11, 15, 19, 2: 1 ng of HeLa cell total RNA used for reaction

Lane 4, 8, 12, 16, 20, 24, 28, 32: 100 pg of HeLa cell total RNA used for reaction

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
Note: Schematic representation of the 5' UTR of a gene, with complex secondary structure, at three different temperatures. Note that RocketScript™ RTase shows full activity at 70°C allowing it to synthesize the complete gene sequence where M-MLV and other reverse transcriptase's fail.

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reactions
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- Enzyme: RocketScript™ RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 10 kb

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-20°C

AccuPower® RocketScript™ RT PreMix & Master Mix

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Thermostable Activity



Figure 1. Sensitivity comparison between AccuPower® RocketScript™ RT PreMix and M-MLV RTase.

Sensitivity results of AccuPower® RocketScript™ RT PreMix using GAPDH compared with conventional reverse transcriptases.

Each 10 fg~100 ng of total RNA used for RT and the same amount of RT products used for electrophoresis.

- Lane 1: 10 fg human total RNA from HeLa cell
- Lane 2: 100 fg human total RNA from HeLa cell
- Lane 3: 1 pg human total RNA from HeLa cell
- Lane 4: 10 pg human total RNA from HeLa cell
- Lane 5: 100 pg human total RNA from HeLa cell
- Lane 6: 1 ng human total RNA from HeLa cell
- Lane 7: 10 ng human total RNA from HeLa cell
- Lane 8: 100 ng human total RNA from HeLa cell

Sensitivity Test

| Concentration | RocketScript™ RT PreMix | Supplier Q | Supplier I |
|---------------|-------------------------|------------|------------|
| 10,000 | 23.91 | 25.63 | 24.43 |
| 1,000 | 27.33 | 28.92 | 28.03 |
| 100 | 30.62 | 32.42 | 30.88 |
| 10 | 33.63 | 35.43 | 33.95 |
| Efficiency | 104% | 103% | 108% |
| Linearity | 0.9999 | 0.9996 | 0.9995 |

Figure 2. Sensitivity comparison between AccuPower® RocketScript™ RT PreMix and Supplier RTases using Real-Time PCR.

Rxn. conditions: conventional 1 hr incubation at 60°C, deactivation at 95°C for 5 min

All cDNAs were amplified with AccuPower® DualStar™ qPCR PreMix (Cat. No. K-6110, Bioneer).

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| 카탈로그 번호 | 제품명 | 규격 | |
|---------|----------------------------------------|--------------------------------------------------|-----------|
| K-2101 | AccuPower® RocketScript™ RT PreMix | 0.2 ml thin-wall 8-tube strips with attached cap | |
| K-2103 | | | 96 tubes |
| K-2102 | | 480 tubes | 20 µl/rxn |
| K-2104 | | | 50 µl/rxn |
| K-2105 | AccuPower® RocketScript™ RT Master Mix | 1 ml of 2X Master mix solution | |

Enhanced Performance


Figure 3. Comparison of amplification efficiency between AccuPower® RocketScript™ RT PreMix (A) and competitors M-MLV RTase (B).

RocketScript™ is able to handle a wide range of sample concentrations and transcript lengths so your downstream applications are minimally affected by the reverse transcription step.

Lanes 1-3: 1,000 ng, 100 ng and 10 ng of total RNA from HeLa cells, respectively.

Note: Competitor products show inhibition at high input concentrations of total RNA

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
Note: Schematic representation of the 5' UTR of a gene, with complex secondary structure, at three different temperatures. Note that RocketScript™ shows full activity at 70°C allowing it to synthesize the complete gene sequence where M-MLV and other reverse transcriptase's fail.

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reactions
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- Enzyme: RocketScript™ RTase
- DNase activity: No
- RNase activity: No
- Fragment size: ~ 10 kb

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-20°C

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Thermostable Activity


Figure 1. Complex RNA amplification results of *AccuPower® RocketScript™* Cycle RT PreMix.

Each target gene (MYC, TFRC) was amplified after performing reverse transcription with *AccuPower® RocketScript™* Cycle RT PreMix.

Rxn. conditions: Conventional 1 hr incubation at 42°C, 50°C, or 60°C, deactivation at 95°C for 5 min


A: *M-MLV* Reverse Transcriptase

B: *AccuPower® RocketScript™* Cycle RT PreMix with Oligo (dT₂₀)

Lane 1: 100 ng Human total RNA from HeLa cell

Lane 2: 10 ng Human total RNA from HeLa cell

Cyclic Temperature Reverse Transcription

| Concentration (copies/rxn) | FTRT (Ct) | CTRT with 1 cycle (Ct) | CTRT with 10 cycles (Ct) |
|----------------------------|-----------|------------------------|--------------------------|
| 10,000 | 19.46 | 18.77 | 18.51 |
| 1,000 | 24.11 | 24.04 | 22.93 |
| 100 | 29.78 | 28.35 | 28.19 |
| 10 | 32.87 | 33.00 | 31.05 |

Figure 3. Low copy species enrichment by cycle.

Comparing FTRT (Fixed Temperature RT) to 1 and 10 cycle(s) of CTRT reveals progressive improvement of detection of cDNA in subsequent qPCR for CTRT as copies of template RNA decrease.

FTRT: 60 min incubation at 50°C followed by 5 min deactivation at 95°C

CTRT: Cycles of 37°C annealing 10 sec, 50°C cDNA synthesis 4 min, 55°C secondary structure melting and cDNA synthesis 30 sec

Target: Human GAPDH

Human total RNA from HeLa cells

qPCR with *AccuPower® GreenStar™* qPCR PreMix (Cat. No. K-6210, Bioneer)

Sensitivity & Full-length cDNA synthesis


Figure 2. Comparison of amplification efficiency between *AccuPower® RocketScript™* Cycle RT PreMix and competitors' RTases.

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M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

RT reaction condition is performed according to each manufacturer's recommendations.

(b) Full-Length cDNA synthesis test

RT reactions were performed according to each manufacturer's recommendation. All cDNAs were amplified with *AccuPower® HotStart* PCR PreMix (K-5050) from Bioneer

Note: Competitor products show inhibition at high input concentrations of total RNA.

Lane 1: 1 µg human total RNA from HeLa cell

Lane 2: 100 ng human total RNA from HeLa cell

Lane 3: 10 ng human total RNA from HeLa cell

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
Note: Schematic representation of the 5'UTR of a gene, with complex secondary structure, at three different temperatures. Note that RocketScript™ shows full activity at 70°C allowing it to synthesize the complete gene sequence where M-MLV and other reverse transcriptase's fail.

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- Enzyme: RocketScript™ RTase, ProFi Taq DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 6 kb

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-20°C

AccuPower® RocketScript™ RT-PCR PreMix & Master Mix

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Figure 1. Performance comparison between *AccuPower® RocketScript™* RT-PCR PreMix and competitor RT-PCR kits.

Complex secondary structure RNA species amplification by *RocketScript™* is outstanding compared to the leading competitor's reverse transcriptase. RT reactions were performed according to each manufacturer's recommendations.

Lane 1: 10 ng of total RNA from HeLa cells

Lane 2: 1 ng of total RNA from HeLa cells

Lane 3: 100 pg of total RNA from HeLa cells

Lane 4: 10 pg of total RNA from HeLa cells


Figure 2. *AccuPower® RocketScript™* RT-PCR PreMix shows enhanced performance compared to competitors.

One-step RT-PCR reactions were performed with 100 ng total RNA from HeLa cells using reagents and conditions specified in each manufacturer's protocol.

Lane 1: 3 kb Lane 2: 4.5 kb Lane 3: 5.2 kb

RT reaction condition is performed according to each manufacturer's recommendations.


Figure 3. Highly reproducible amplifications.

Amplification of an 500 bp target gene was detected using human total RNA (from 10 ng to 10 pg) with *AccuPower® RocketScript™* RT-PCR PreMix.

As shown in figure 3, highly reproducible amplifications were achieved within each Lot. set of triplicates

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- First-strand synthesis of cDNA from RNA molecules (RT)
- RT-PCR
- Random priming reactions
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- Enzyme: RocketScript™ RTase, RNase H Minus
- DNase activity: No
- RNase activity: No
- RNase H activity: No
- Fragment size: ~ 12.5 kb

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-20°C

AccuPower® RocketScript™ RT PreMix, RNase H Minus & Master Mix

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
Figure 1. Performance of the AccuPower® RocketScript™ RT PreMix, RNase H Minus in two-step RT-PCR.

Target mRNAs ranged from 500 bp to 11 kb were reverse transcribed by using AccuPower® RocketScript™ RT PreMix, RNase H Minus. With each cDNA, PCR was performed by using AccuPower® ProFi Taq PCR PreMix (Cat. No. K-2631, Bioneer). Template sizes are indicated above the gel image.

M1; 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

M2; Lambda DNA/Hind III Markers (Cat. No. D-1050, Bioneer).


Figure 2. Thermal stability using BIONEER AccuPower® RocketScript™ RT PreMix, RNase H Minus.

Synthesis of cDNA at 50°C-70°C for 30 mins using 100 ng, 10 ng, 1 ng and 100 pg of 600 bp RNA transcript as a template using RocketScript™ RNase H(-) RT PreMix. RocketScript™ RNase H(-) RT PreMix completed synthesis of 70°C.



Figure 3. High synthesis rate at 50°C using BIONEER AccuPower® RocketScript™ RT PreMix, RNase H Minus.

Synthesis of cDNA at 50°C or 42°C for 10, 20, 30 and 60 min using 1 µg of 9 kb RNA transcript as a template using Bioneer cDNA synthesis kit and other kits for first strand cDNA synthesis. Reaction products were resolved on a 1% alkaline agarose gel. Bioneer cDNA synthesis kit completed synthesis of 9 kb transcript in 10 min.

Lane 1: RT 10 min Lane 2: RT 20 min Lane 3: RT 30 min Lane 4: RT 60 min
M; 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

Supplier I- S III Reverse Transcriptase


Figure 4. Sensitivity comparison of the various target.


Rxn. conditions: Conventional 30 min incubation at 50°C deactivation at 95°C for 5 min

Lane 1: 10 ng of total RNA from HeLa cells

Lane 2: 1 ng of total RNA from HeLa cells

Lane 3: 100 pg of total RNA from HeLa cells

Lane 4: 10 pg of total RNA from HeLa cells

| Human RNA Concentration (pg) | Ct |
|------------------------------|-------|
| 10 ⁶ | 18.81 |
| 10 ⁵ | 21.96 |
| 10 ⁴ | 25.29 |
| 10 ³ | 29.16 |
| 10 ² | 32.43 |
| 10 | 35.70 |
| 1 | 39.01 |

Figure 5. Broad dynamic range of AccuPower® RocketScript™ RT PreMix, RNase H Minus.

First-strand cDNA was generated using the AccuPower® RocketScript™ RT PreMix, RNase H Minus. cDNA was amplified using the AccuPower® Plus DualStar™ qPCR PreMix on the Exicycler™ 96 Real-Time Quantitative Thermal Block from Bioneer. The standard curve illustrates high linearity ($R^2 = 0.999$) across a broad range of input RNA, suggesting that the relative representation of specific RNA transcripts is preserved in the cDNA pool regardless of the abundance of total RNA. Amplification was performed on 10-fold serial dilutions of HeLa total RNA (1 µg to 1 pg).

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- Low copies detection
- RNA virus detection
- Gene expression analysis

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- Enzyme: RocketScript™ RTase (RNase H Minus), ProFi Taq DNA polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: Yes
- 3' - A Overhang: Yes
- RNase H Activity: No
- Fragment size: ~ 12.5 kb

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-20°C

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
Figure 1. Performance of the AccuPower® RocketScript™ RT-PCR Pre-Mix, RNase H Minus.

AccuPower® RocketScript™ RT-PCR PreMix, RNase H minus was used to amplify cDNAs which size are between 2.5 kb to 12.5 kb. DNA sizes are indicated above the picture. Reactions were performed with 100 ng total RNA from HeLa cells.

M1: Lambda/HindIII marker (Cat. No. D-1050, Bioneer)

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
- Multiplex RT-PCR
- Low copy detection
- Gene expression analysis

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- Enzyme: RocketScript™ RTase, HotStart Top DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 1 kb

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-20°C

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
Figure 1. Single RT-PCR and multiplex RT-PCR using AccuPower® RocketPlex RT-PCR PreMix.

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AccuPower® RocketPlex RT-PCR Premix & Master Mix

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
Figure 2. Comparison of amplification quality between *AccuPower*® RocketPlex RT-PCR PreMix and other suppliers' RT-PCR kit.

10-plex primers were added into *AccuPower*® RocketPlex RT-PCR PreMix and other supplier's RT-PCR kit. A series of human total RNA diluents were tested. All data were obtained using *MyGenie*™ 96 Gradient Thermal Block (Bioneer). *AccuPower*® RocketPlex RT-PCR PreMix is able to perform reverse transcription reactions throughout a wide range of temperatures from 42°C to 70°C.

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- RNase H activity: Yes
- Fragment size: ~ 3 kb

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-20℃



Figure 1. Dual-HotStart™ Technology

AccuPower® Dual-HotStart™ RT-PCR PreMix

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
Figure 2. Comparison of PCR amplification sensitivity between AccuPower® Dual-HotStart™ RT-PCR PreMix from Bioneer and other suppliers' HotStart RT-PCR kit.

Target: human PGK1
 Lane 1: 10 ng of human RNA
 Lane 2: 1 ng of human RNA
 Lane 3: 100 pg of human RNA
 Lane 4: 10 pg of human RNA

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|-----------------------------------------|--------------------------------------------------|-----------|-----------|
| K-6710 | AccuPower® Dual-Hotstart™ RT-PCR PreMix | 0.2 ml thin-wall 8-tube strips with attached cap | 96 tubes | 20 µl/rxn |
| K-6711 | | | | 50 µl/rxn |
| K-6712 | | | 480 tubes | 20 µl/rxn |
| K-6713 | | | | 50 µl/rxn |

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: ~ 3 kb

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AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG)

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RNA template (not including uracil base)

RNA template (including uracil base)


Figure 2. Comparison of amplification quality using PCR products (not including uracil base or including uracil base) between AccuPower® Dual-HotStart™ RT-PCR PreMix and AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG).

Lane 1: 10^8 Lane 2: 10^7 Lane 3: 10^6 Lane 4: 10^5 Lane 5: 10^4 Lane 6: 10^3 Lane 7: 10^2 Lane N: NTC


Figure 3. Comparison of PCR amplification specificity between AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG) from Bioneer and other suppliers' Hotstart PCR master mix.

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- First strand cDNA synthesis from RNA
- RT-PCR
- and qRT-PCR

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- DNase activity : No
- RNase activity: No
- 3'-A overhang: No
- Strand displacement: Yes
- Fragment size: ~ 9 kb

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- 5X Reaction buffer: Tris-HCl, KCl, MgCl₂ (pH 8.1)
- 100 mM DTT
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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10,000 Units (200 U/μl)

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50% glycerol containing 20 mM Tris-HCl (pH 7.6), 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50 mM (NH₄)₂SO₄, Stabilizers.

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-20°C

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One Unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 min at 37°C using poly(A), oligo dT as template primer.

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| 카탈로그 번호 | 제품명 |
|---------|------------------------------------------------------------------------------------|
| E-3121 | M-MLV Reverse Transcriptase, 10,000 U, 5X Reaction Buffer, 100 mM DTT, 10 mM dNTPs |
| E-3122 | M-MLV Reverse Transcriptase, 50,000 U, 5X Reaction Buffer, 100 mM DTT, 10 mM dNTPs |

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
- First-strand synthesis of cDNA from RNA molecules
- RT-PCR
- Random priming reaction
- Library construction
- Probe labeling
- mRNA 5' end mapping by primer extension analysis

CycleScript™ Reverse Transcriptase

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- DNase activity : No
- RNase activity: No
- 3'-A overhang: No
- Strand displacement: Yes
- Fragment size: ~ 9 kb

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- 5X Reaction Buffer: Tris-HCl, KCl, MgCl₂ (pH 8.1)
- 100 mM DTT
- dNTPs mixture: 10 mM, each dNTP 2.5 mM

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10,000 Units (200 U/μl)

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50% glycerol containing 20 mM Tris-HCl (pH 7.6), 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50 mM (NH₄)₂SO₄, Stabilizers.

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
-20°C

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One Unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 min at 37°C using poly(A), oligo dT as template primer.

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CycleScript™ Reverse Transcriptase

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
Figure 1. Comparison of transferrin receptor gene amplification with different reverse transcriptases.

700 ng of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1 - 4: TFR (Transferrin receptor gene) amplified with *M-MLV*

Lane 5 - 8: TFR amplified with *CycleScript™*

Lane 9 - 12: TFR amplified with *CycleScript™*

Lane 13 - 16: TFR amplified with *M-MLV* from Company I

Lane 17 - 20: TFR amplified with *S-script* from Company S

Lane 21 - 24: TFR amplified with *S-script II* from Company I

Lane 25 - 28: TFR amplified with *S-script III* from Company I

Lane 29 - 32: TFR amplified with *O-script* from Company Q

Lane M.W.: 100 bp Plus DNA Ladder (Cat. No. D-1035, Bioneer)


Figure 2. Comparison of GAPDH gene amplification with different reverse transcriptases.

Each 10 ng, 1 ng, 100 pg, and 10 pg of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1 - 4: GAPDH amplified with *CycleScript™*

Lane 5 - 8: GAPDH amplified with *CycleScript™*

Lane 9 - 12: GAPDH amplified with *CycleScript™*

Lane 13 - 16: GAPDH amplified with *M-MLV* from Company I

Lane 17 - 20: GAPDH amplified with *S-script* from Company S

Lane 21 - 24: GAPDH amplified with *S-script* from Company S

Lane 25 - 28: GAPDH amplified with *S-script* from Company I

Lane 29 - 32: GAPDH amplified with *S-script* from Company I


Figure 3. Working temperature comparison of different reverse transcriptases.

Each 10 ng, 1 ng, 100 pg, and 10 pg of total RNA was used for reverse transcription and the same amount of amplified products were used for electrophoresis.

Lane 1 - 4: GAPDH amplified with *CycleScript™*

Lane 5 - 8: GAPDH amplified with *CycleScript™*

Lane 9 - 12: GAPDH amplified with *CycleScript™*

Lane 13 - 16: GAPDH amplified with *S-script* from company S

Lane 17 - 20: GAPDH amplified with *S-script* from company S

Lane 21 - 24: GAPDH amplified with *S-script* from company S

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| 카탈로그 번호 | 제품명 |
|---------|--------------------------------------------------------------------------------------------------|
| E-3131 | <i>CycleScript™</i> Reverse Transcriptase, 10,000 U, 10 mM dNTPs, 5X Reaction Buffer, 100 mM DTT |
| E-3132 | <i>CycleScript™</i> Reverse Transcriptase, 50,000 U, 10 mM dNTPs, 5X Reaction Buffer, 100 mM DTT |

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- Gene synthesis
- First-strand synthesis of cDNA from RNA molecules (Reverse Transcription)
- RT-PCR
- Random priming reactions
- Library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- DNase activity : No
- RNase activity: No
- 3'-A overhang : No
- Fragment size: ~10 kb

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- 5X Reaction Buffer: Tris-HCl, KCl, MgCl₂, Stabilizers (pH 8.5)
- 100 mM DTT
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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10,000 Units (200 U/μl)

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50% glycerol containing 20 mM Tris-HCl (pH 7.6), 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, Stabilizers.

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-20°C

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One Unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 min at 37°C using poly(A), oligo dT as template primer.

RocketScript™ Reverse Transcriptase

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

Figure 1. Complex RNA amplification.

Complex RNA amplification results of *RocketScript™* Reverse Transcriptase compared with that of conventional RTase.

Rxn. condition: conventional 42°C / 50°C / 60°C 1 hr, deactivation 95°C 5 min, this product shows thermal stability.


Lane A: *M-MLV* Reverse Transcriptase

Lane B: *RocketScript™* Reverse Transcriptase

Lane 1: 100 ng human total RNA from HeLa cell

Lane 2: 10 ng human total RNA from HeLa cell

Rocketscript™ RTase의 thermostable Activity

I Company RTase의 thermostable Activity


Figure 2. Thermostable stability check.

Amplification results of *RocketScript™* Reverse Transcriptase using myc compared with supplier I Reverse transcription.

Rxn. condition: Incubation at each temperature 45°C/50°C/55°C/60°C/65°C/70°C for 1 hr, deactivation at 95°C for 5 min

Target: human myc 495 bp

Lane 1: 100 ng human total RNA from HeLa cell

Lane 2: 10 ng human total RNA from HeLa cell

Lane 3: 1 ng human total RNA from HeLa cell

Lane 4: 100 pg human total RNA from HeLa cell

M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)

Figure 3. Comparison of amplification quality between *RocketScript™* RTase and competitor RTases.

Target gene expression Level.

Lane 1: 100 ng human total RNA from HeLa cell

Lane 2: 10 ng human total RNA from HeLa cell

Lane 3: 1 ng human total RNA from HeLa cell

Lane 4: 100 pg human total RNA from HeLa cell

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- Gene synthesis
- First-strand synthesis of cDNA from RNA molecules (Reverse Transcription)
- RT-PCR
- Random priming reactions
- Library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-Time PCR

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- DNase activity : No
- RNase activity: No
- RNase H activity: No
- Fragment size: ~12.5 kb

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- 5X Reaction Buffer: Tris-HCl, KCl, MgCl₂, Stabilizers (pH 8.3)
- 100 mM DTT
- 10 mM dNTPs mix: 2.5 mM of each dNTP

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10,000 Units (200 U/μl)

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50% glycerol containing 20 mM Tris-HCl (pH 7.6), 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, Stabilizers.

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-20°C

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One Unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into acid-precipitable material in 10 min at 37°C using poly(A), oligo dT as template primer.

RocketScript™ Reverse Transcriptase, RNase H Minus

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
Figure 1. Amplification results of *RocketScript*™ Reverse Transcriptase, RNase H Minus. (A) Sensitivity (B) Thermostability.

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer), Lane 1: 100 ng human total RNA from HeLa cell
 Lane 2: 10 ng human total RNA from HeLa cell, Lane 3: 1 ng human total RNA from HeLa cell
 Lane 4: 100 pg human total RNA from HeLa cell, Lane 5: 10 pg human total RNA from HeLa cell


Figure 2. Comparison of amplification quality among of *RocketScript*™ Reverse Transcriptase, RNase H Minus and other company.

(A) GAPDH 0.5 kb, (B) TUBB 1 kb.

M: 100 bp DNA Ladder (Cat. No. D-1030, Bioneer) Lane 1: 10 ng human total RNA from HeLa cell
 Lane 2: 1 ng human total RNA from HeLa cell, Lane 3: 100 pg human total RNA from HeLa cell
 Lane 4: 10 pg human total RNA from HeLa cell


Figure 3. Comparison of long kb amplification among of *RocketScript*™ Reverse Transcriptase, RNase H Minus and other company.

M1: Lambda/*Hind* III marker (Cat. No. D-1050, Bioneer), M2: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)
 Lane 1: 0.5 kb, Lane 2: 1 kb, Lane 3: 2.5 kb, Lane 4: 3.4 kb,
 Lane 5: 4.7 kb, Lane 6: 6 kb, Lane 7: 7.1 kb, Lane 8: 8 kb, Lane 9: 9 kb

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| 카탈로그 번호 | 제품명 |
|---------|---------------------------------------------------------------------------------|
| E-3161 | <i>RocketScript</i> ™ Reverse Transcriptase, RNase H Minus, 10,000 U (50 rxns) |
| E-3162 | <i>RocketScript</i> ™ Reverse Transcriptase, RNase H Minus, 50,000 U (250 rxns) |

03. Real-Time PCR

dsDNA Binding Dye Type Kit

| | |
|--------------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>GreenStar</i> [™] qPCR PreMix | 146 |
| <i>AccuPower</i> [®] 2X <i>GreenStar</i> [™] qPCR Master Mix | 148 |
| <i>AccuPower</i> [®] <i>GreenStar</i> [™] RT-qPCR PreMix & 2X Master Mix | 150 |


Hydrolysis Probe Type Kit

| | |
|------------------------------------------------------------------------------------------------------------|-----|
| <i>AccuPower</i> [®] <i>DualStar</i> [™] qPCR PreMix | 151 |
| <i>AccuPower</i> [®] Plus <i>DualStar</i> [™] qPCR PreMix & 2X Master Mix | 153 |
| <i>AccuPower</i> [®] Plus <i>DualStar</i> [™] qPCR PreMix & 2X Master Mix (with UDG) ... | 155 |
| <i>AccuPower</i> [®] <i>Dual-HotStart</i> [™] RT-qPCR PreMix & 2X Master Mix | 157 |

Real-Time PCR Instrument

Exicycler[™] 96/384 → Go to M. Instruments & Devices

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- Real-Time quantification of DNA and cDNA targets
- Gene expression profiling
- Gene functional analysis
- Microbial & viral pathogen detection

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- Enzyme: *Taq* DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes


Figure 1. Enzyme-mediated HotStart PCR (PyroHotStart 기술)

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High Specificity


Figure 2. Real-Time PCR data of AccuPower® GreenStar™ qPCR PreMix.

AccuPower® GreenStar™ qPCR PreMix provides at least 7 orders of magnitude in dynamic range (10 fg~10 ng /rxn).

A: Amplification curve of AccuPower® GreenStar™ qPCR PreMix.

Lambda DNA primers were added into AccuPower® GreenStar™ qPCR PreMix. A series of lambda DNA positive control diluents were tested.

B: Standard curve of AccuPower® GreenStar™ qPCR PreMix.

$R^2 = 0.998$, PCR efficiency = 95%

C: Melting curve analysis of AccuPower® GreenStar™ qPCR PreMix. The melting curve shows that only single amplified PCR product was obtained in all template range.

All data were obtained using Exicycler™ 96 Real-Time Quantitative Thermal Block (Cat. No. A-2060, Bioneer).


Figure 3. Comparison of amplification efficiency between AccuPower® GreenStar™ qPCR PreMix and other company's master mix.

Amplification curve of AccuPower® GreenStar™ qPCR PreMix and other company's master mix kit. Lambda DNA primers were added into AccuPower® GreenStar™ qPCR PreMix and other company's master mix kit. A series of Lambda DNA positive control diluents were tested. Reaction mixtures were prepared and qPCR was performed according to each company's protocol. All data were obtained using ABI7500 Fast Real-Time PCR machine (Applied Biosystems co.).

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| 카탈로그 번호 | 제품명 | 규격 | | | | |
|---------|-----------------------------------|---------------|---------------|-----------------------|-----------|-----------|
| K-6210 | AccuPower® GreenStar™ qPCR PreMix | Exicycler™ | 8-well strip | Optical film included | 96 rxn | 20 µl/rxn |
| K-6200 | | | | | | 50 µl/rxn |
| K-6213 | | 96-well plate | 20 µl/rxn | | | |
| K-6203 | | | 50 µl/rxn | | | |
| K-6211 | | ABI 7500 | 8-well strip | | 20 µl/rxn | |
| K-6201 | | | | 50 µl/rxn | | |
| K-6214 | | | 96-well plate | | 20 µl/rxn | |
| K-6204 | | | | 50 µl/rxn | | |
| K-6212 | | CFX96 | 8-well strip | | | 20 µl/rxn |
| K-6202 | | | | - | 50 µl/rxn | |

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- Real-Time quantification of DNA and cDNA targets
- Gene expression profiling
- Gene functional analysis
- Microbial & viral pathogen detection

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
- Enzyme: *Taq* DNA polymerase
- 5' → 3' exonuclease: No
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes

AccuPower® 2X GreenStar™ qPCR Master Mix

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
Equipment Compatibility

▶ Exicycler™ 96 of BIONEER

▶ Using a Bio-Rad IQ5

▶ Using a ABI7500


Figure 1. Comparison of the specificity of the dsDNA binding dye based Real-Time PCR.

Amplification of an 90 bp target gene was detected using serially diluted LP (Legionella Pneumoniae) genomic DNA (10n dilution: 10⁵~10 copies) with AccuPower® 2X GreenStar™ qPCR Master Mix. As shown in figure, very small amount of primer dimers was appeared in AccuPower® 2X GreenStar™ qPCR Master Mix than other kits.

All data were obtained using Exicycler™ 96 Real-Time Quantitative Thermal Block (Cat. No. A-2060, Bioneer).

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
Figure 2. Gene expression analysis.

AccuTarget™ Validated Real-Time PCR Primer Library for human is designed by Bioneer's bioinformatics tool and targeting for human genome. cDNA was synthesized using human PTGS2 target primer of those and human total RNA identically quantified from Hela cell and blood cell with AccuPower® CycleScript™ RT PreMix (Cat. No. K-2044, Bioneer). Gene analysis was carried out both Hela cell and blood cell by operating Real-Time PCR reaction (95°C 10 min, 1 cycle and 95°C 10 sec, 58°C 25 sec, 72°C 30 sec, 41 cycles) using the cDNA, AccuPower® 2X GreenStar™ PCR Master Mix and Exicycler™ 96 Real-Time Thermal Block (Cat. No. A-2060, Bioneer).

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
Figure 1. Comparison of specificity between AccuPower® GreenStar™ RT-qPCR PreMix and other supplier's master mixtures.

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| 카탈로그 번호 | 제품명 |
|---------|--------------------------------------------------------------------------------------------------------------|
| K-6400 | AccuPower® GreenStar™ RT-qPCR PreMix, 50 µl/rxn, 8-tube strips, 96 rxn, Exicycler™ 96, optical film included |
| K-6403 | AccuPower® GreenStar™ RT-qPCR Master Mix (2X), 2.5 ml, 100 rxn |

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- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

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- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes
- Fragment size: 200 bp

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- Gene expression profiling
- Target DNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination
- Evaluation of primer pair performance for probe-based Real-Time PCR

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- Enzyme: Taq DNA polymerase
- 5' → 3' exonuclease: Yes
- 3' → 5' exonuclease: No
- 3' - A overhang: Yes

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
Figure 2. Real-Time PCR data of AccuPower® DualStar™ qPCR PreMix.

AccuPower® DualStar™ qPCR PreMix provides at least 7 orders of magnitude in dynamic range ($10^7 \sim 10$ copies/rxn).

A: Amplification curve. West Nile Virus (WNV) primers and TaqMan-based probe were added into AccuPower® DualStar™ qPCR PreMix. A series of WNV positive control diluents were tested.

B: Standard curve. All data were obtained using Exicycler™ 96 Real-Time Quantitative Thermal Block (Cat. No. A-2060, Bioneer).

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
Figure 3. Data using various kinds of Real-Time PCR instruments.

AccuPower® DualStar™ qPCR PreMix is applicable to most of commercially available Real-Time quantitative PCR instruments. West Nile Virus (WNV) primers and TaqMan-based probe were added into AccuPower® DualStar™ qPCR PreMix. A series of WNV positive control diluents were tested.

A: Amplification curve and standard curve using ABI7500 Fast Real-Time PCR machine (Applied Biosystems Co.).

B: Amplification curve and standard curve using ABI7500 Real-Time PCR machine (Applied Biosystems Co.).

C: Amplification curve and standard curve using Opticon (CFX 96) Real-Time PCR machine (MJ Research, currently Bio-Rad Inc.).

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| 카탈로그 번호 | 제품명 | 규격 | | | | |
|---------|----------------------------------|---------------|-----------------------|-----------------------|-----------|-----------|
| K-6100 | AccuPower® DualStar™ qPCR PreMix | Exicycler™ | 8-well strip | Optical film included | 96 rxn | 20 µl/rxn |
| K-6110 | | | | | | 50 µl/rxn |
| K-6103 | | 96-well plate | Optical film included | 96 rxn | 20 µl/rxn | |
| K-6113 | | | | | 50 µl/rxn | |
| K-6101 | | ABI 7500 | 8-well strip | Optical film included | 96 rxn | 20 µl/rxn |
| K-6111 | | | | | | 50 µl/rxn |
| K-6104 | | 96-well plate | Optical film included | 96 rxn | 20 µl/rxn | |
| K-6114 | | | | | 50 µl/rxn | |
| K-6102 | | CFX96 | 8-well strip | - | 96 rxn | 20 µl/rxn |
| K-6112 | | | | | | 50 µl/rxn |


Figure 4. Comparison of detection sensitivity between AccuPower® DualStar™ qPCR PreMix and other suppliers' master mix.

West Nile Virus (WNV) primers and TaqMan-based probe were added into AccuPower® DualStar™ qPCR PreMix and other suppliers' master mixture. A series of WNV positive control diluents were tested. Reaction mixtures were prepared and qPCR was performed according to each supplier's protocol. All data were obtained using ABI7500 Fast Real-Time PCR machine (Applied Biosystems co.).

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- Gene expression profiling
- Target DNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination
- Evaluation of primer pair performance for probe-based Real-Time PCR

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- Enzyme: HotStart Taq DNA polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A Overhang: Yes

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-20°C

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Specificity


Figure 1. High specificity of AccuPower® Plus DualStar™ qPCR PreMix.


Sensitivity

| Template DNA copies | Ct Value | | |
|----------------------|----------|-----------|-------------|
| | BIONEER | A Company | B Company |
| 10 ³ | 32.6 | 34.72 | 31.28 |
| 10 ⁴ | 29.18 | 31.31 | 27.75 |
| 10 ⁵ | 25.66 | 27.61 | 24.44 |
| 10 ⁶ | 22.43 | 24.48 | 21.99 |
| 10 ⁷ | 19.22 | 20.98 | 19.99 |
| Efficiency/linearity | 99%/1 | 96%/1 | 100%/0.9999 |

Figure 2. Comparison of amplification quality between AccuPower® Plus DualStar™ qPCR PreMix and other suppliers' Real-Time qPCR kit.


AccuPower® Plus DualStar™ qPCR PreMix & 2X Master Mix

| Dye | Template DNA of copies | | | | | | | |
|-------------|------------------------|----------|----------|----------|----------|----------|----------|-----|
| | 1.00E+07 | 1.00E+06 | 1.00E+05 | 1.00E+04 | 1.00E+03 | 1.00E+02 | 1.00E+01 | NTC |
| FAM | 19.13 | 22.8 | 25.17 | 28.48 | 31.53 | 33.95 | 35.98 | UD |
| TET | 19.27 | 23.11 | 25.83 | 28.89 | 31.31 | 34.15 | 35.76 | UD |
| CY5 | 18.48 | 21.52 | 24.33 | 27.35 | 30.02 | 32.86 | 35.45 | UD |
| TEXAS_RED | 18.19 | 21.56 | 24.21 | 27.8 | 30.42 | 33.39 | 35.19 | UD |
| TAMRA (IPC) | 24.65 | 24.75 | 24.51 | 24.41 | 24.66 | 24.73 | 24.83 | UD |

Figure 3. Five-target multiplexing on the Exicycler™ 96 instrument using AccuPower® Plus DualStar™ qPCR PreMix.

Figure 3 shows amplification of a 5-target multiplex assay. The dyes used were FAM, TET, Cyanine 5, Texas Red and TAMRA, respectively. The data demonstrate that over a dilution series of input template, the AccuPower® Plus DualStar™ qPCR PreMix can successfully and reliably generate up to 5-target multiplex data on the Exicycler™.

| PCR inhibitor | | Bioneer | A company | B company |
|---------------|------------|--------------------------|-----------|-----------|
| | | Totally inhibition (PPM) | | |
| Humic acid | Soil | 10,000 | * | 10,000 |
| | Peat | 10,000 | 10,000 | 10,000 |
| | Leonardite | 10,000 | 10,000 | 10,000 |
| | Fulvic | 10,000 | 10,000 | 10,000 |
| Hemoglobin | | * | * | * |
| Blood EDTA | | 10,000 | 10,000 | 10,000 |

* : No inhibition

Figure 4. PCR inhibitor study using AccuPower® Plus DualStar™ qPCR PreMix.

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|-------------------------------------------|----------------------------------|------------------------------------------------|-----------|
| K-6100 | AccuPower® Plus DualStar™ qPCR PreMix | Exicycler™ | 8-well strip (96 rxn) optical film included | 50 µl/rxn |
| K-6110 | | ABI7500 | | |
| K-6103 | | CFX96 | | |
| K-6112 | AccuPower® Plus DualStar™ qPCR Master Mix | 2.5 ml of 2X Master mix solution | | |

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
AccuPower® Plus DualStar™ qPCR PreMix & 2X Master Mix (with UDG)

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| PCR product type | Linearity | Efficiency | 10 ⁸ | 10 ⁷ | 10 ⁶ | 10 ⁵ | 10 ⁴ | 10 ³ | 10 ² | NTC |
|-----------------------|-----------|------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----|
| dNTPs | 0.9988 | 98 | 12.81 | 16.41 | 19.72 | 22.84 | 26.37 | 29.76 | 33.48 | UD |
| dNTPs including dUTPs | - | - | UD | UD |

Figure 2. Efficiency of uracil DNA glycosylase using PCR products (including dN/dU base).

| | | | Ct Value | | | | | | | |
|-----------|-----------|------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----|--|
| Company | Linearity | Efficiency | 10 ⁷ | 10 ⁶ | 10 ⁵ | 10 ⁴ | 10 ³ | 10 ² | NTC | |
| Bioneer | 0.9999 | 98% | 15.73 | 19.14 | 22.99 | 25.91 | 28.96 | 33.1 | UD | |
| A company | 0.9994 | 97% | 16.86 | 19.98 | 23.55 | 27.26 | 30.05 | 34.61 | UD | |

Figure 3. Comparison of amplification quality between AccuPower® Plus DualStar™ qPCR PreMix (with UDG) and other supplier's Real-time qPCR kit.

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|------------------------------------------------------|----------------------------------|-----------------------|-----------|
| K-6605 | AccuPower® Plus DualStar™ qPCR PreMix (with UDG) | Exicycler™ 96-well plate | optical film included | 50 µl/rxn |
| K-6606 | | ABI7500 96-well plate | | |
| K-6607 | | CFX96 8-well strip | | |
| K-6608 | AccuPower® Plus DualStar™ qPCR Master Mix (with UDG) | 2.5 ml of 2X Master mix solution | | |

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- Gene expression profiling
- Target RNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination

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- Enzyme: RocketScript™ RTase, HotStart Taq DNA polymerase
- 5' → 3' exonuclease activity: Yes
- 3' → 5' exonuclease activity: No
- 3' - A Overhang: Yes
- RNase H activity: Yes

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-20°C



Figure 1. Dual-HotStart™ technology

AccuPower® Dual-HotStart™ RT-qPCR PreMix & 2X Master Mix

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
Figure 2. High sensitivity of AccuPower® Dual-HotStart™ RT-qPCR PreMix. Experiment with HIV target. 10 fold serial dilution of template RNA (10^{10} copies ~ 10 copies spiked in human total RNA).


Figure 3. High specificity of AccuPower® Dual-HotStart™ RT-qPCR PreMix. Experiment with HCV target. 10 fold serial dilution of template RNA (10^6 copies ~ 10 copies spiked in human total RNA). Conventional hotstart qPCR always generate non-specific amplification at low template concentration, which deteriorate the sensitivity of qPCR. AccuPower® Dual-HotStart™ RT-qPCR Pre-Mix accurately amplifies target RNA without non-specific amplification, even at low concentration of template

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| 카탈로그 번호 | 제품명 | 규격 | | |
|---------|-------------------------------------------|----------------------------------|------------------------------------------------|-----------|
| K-6704 | AccuPower® Dual-HotStart™ qPCR PreMix | Exicycler™ | 8-well strip (96 rxn) optical film included | 50 µl/rxn |
| K-6705 | | ABI7500 | | |
| K-6706 | | CFX96 | | |
| K-6707 | AccuPower® Dual-HotStart™ qPCR Master Mix | 2.5 ml of 2X Master mix solution | | |

04. Customized PCR

| | |
|--------------------------------------------------------------------------|-----|
| Customized <i>AccuPower</i> [®] PCR/RT/qPCR PreMix Service..... | 160 |
|--------------------------------------------------------------------------|-----|

Customized AccuPower® PCR/RT/qPCR PreMix Service

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