
M Instruments & Devices

- 01. Conventional PCR: *AllInOneCycler™*
- 02. Real-Time PCR: *Exicycler™ 96 / Exicycler™ 96 Fast / Exicycler™ 384 / Exicycler™ HD / Injector™ / ExiStation™ HT*
- 03. Protein Synthesis and Purification: *ExiProgen™*
- 04. DNA/RNA Preparation: *ExiPrep™16 Plus / ExiPrep™48 / ExiPrep™96 Lite / ExiCracker™ / ExiBeater™*
- 05. Electrophoresis: *Agaro-Power™ / DUALED Blue/White Transilluminator*
- 06. Microbial Culture: *HT-MegaGrow®*
- 07. Vortexing & Spin-down: *ExiSpin™ / ExiSpin™ 96*

o Instruments & Devices

☎ Phone: 1588-9788 (ext.4 → 3)

✉ Email: Instrument-support@bioneer.co.kr

01. Conventional PCR

<i>AllInOneCycler</i> [™] PCR System	418
---	-----

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Figure 8. AllInOneCycler™의 앞, 뒤, 측면 디자인.

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Physical Specifications		
Dimension (cm)	24.6(W) x 36.9(D) x 23.2(H)	
Weight (kg)	8.4	
Display	7-inch Touch LCD screen	
Control	Embedded Board	
Power consumption	100-240 V~, Max 8.5 A, 50/60 Hz, 630 VA	
Operating temperature	15~35℃	
Operating humidity	20~80%, no condensation	
Port	2 USB ports for data storage	
Operating specifications		
Method of heating / cooling	Peltier element	
Temperature range	4.0~99.9℃	
Gradient range	20~95℃	
Block type	Fast Block	Normal Block
Max. heating ramp rate	9.5℃/sec	6.5℃/sec
Max. cooling ramp rate	7.7℃/sec	4.5℃/sec
Temperature accuracy	± 0.3℃	
Temperature uniformity	± 0.3℃	
Lid temperature	90~110℃	
Ramp rate control	1~100%	

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Exicycler™ 384 Real-Time Quantitative Thermal Block 427

Exicycler™ HD Ultimate Multiplex Real-Time PCR Platform 430

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
Figure 2. Fluorescence data using 10^6 copies of IRF3 gene (FAM labeled) in each of 96 well positions. The average Ct of 96 well is 21.8 and the Ct variation range is 0.19.

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

Figure 3. Comparison of operating time and ramping time between the previous model and the newest model of Exicycler™ 96.

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FAM TET TAMRA Texas Red Cyanine5

Figure 4. 5-target genes can be detected in a single tube. (FAM: T. vaginalis, TET: M. Hominis, TAMRA: TMV, Texas Red: HSV type1, Cyanine5: HSV type2)

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
Figure 5. Exicycler™ 96 has superior well-to-well uniformity.

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
Figure 6. Graph shows standard curve of tenfold serial dilutions of 10 copies to 10^{11} copies MMP9 gene (FAM labeled). The PCR efficiency generated by the standard curve is 103%.


Figure 7. Fluorescence data from a series of 1.33-fold dilutions of TMV gene (10^6 copies) amplified using reporter dyes to check one target: FAM/TMV. The PCR efficiency generated by the standard curve is 101%.

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
Figure 8. Exicycler™ 96 Analysis Software.

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Physical Specifications		
Dimension (cm)	35.5(W) x 54(D) x 47(H)	
Weight (kg)	39	
Sample Capacity/ Size	Fast block	Opaque White 96-well Low Profile PCR Plate 0.1 ml Opaque White 8-strip Low profile PCR Tube
	Normal block	Opaque White 96-well PCR Plate/0.2 ml Opaque White 8-strip PCR Tube
Sample Volume (μl)	Fast block	10~50 (recommended 20 μl)
	Normal block	20~100 (recommended 50 μl)
Power Consumption	100-240 VAC, 50/60 Hz, Max 800 VA (Fuse: 250 V, F10AL, 2 ea)	
Optical Specifications		
Light Source	Short arc lamp	
Wattage	120 W	
Sensor	16-bit 2D CCD	
Excitation Filter / Emission Filter	5 Sets	
Operating Specifications		
Method of Heating / Cooling	Peltier element	
Temperature Range	4.0~99.9℃	
Temperature Accuracy	± 0.3℃	
Temperature Uniformity	± 0.3℃	
Ramping Rate	Fast block	Max 5.0℃/sec
	Normal block	Max 4.5℃/sec
Lid Temperature	90~120℃	
Gradient Range	20~95℃	
Temperature Increment Range	0.1~2.0℃	
Time Increment Range	1~60 sec	
Ramp Rate Control Range	1~100%	
Operating Temperature	15~35℃	
Operating Humidity	20~80%, no condensation	
Computer Specifications		
Operating System	Windows 7&10 (32/64-bit OS)	
Processor Speed	Intel Dual Core E2160 (1.8 GHz) or higher	
Memory	1 GB or higher	
Communication Port	USB 2.0 high speed	
Screen Resolution	1280X1024 or higher	

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
Figure 2. qPCR result using 1×10^6 copies of Lambda DNA (FAM labeled) in each of 384 well positions. The average Ct of 384 well is 21.36, the Ct variation range is 0.39, the CV value is 0.3%.


Figure 3. qPCR results using 1×10^6 copies of Lambda DNA (FAM labeled) in each of 384 well positions. (a) Bioneer's Exicycler™ 384; Ct average is 21.65, Ct variation is 0.43 and CV value is 0.28%. (b) competitor A; Ct average is 22.16, Ct variation is 1.79 and CV value is 1.13%. (c) competitor B; Ct average is 22.74, Ct variation is 1.69 and CV value is 1.09%.

Exicycler™ 384 Real-Time Quantitative Thermal Block

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
Figure 4. 5-target genes can be detected in a single tube with a minimum volume of 5 µl (FAM: T. vaginalis, TET: M. Hominis, TAMRA: TMV, Texas Red: HSV type1, Cyanine5: HSV type2)

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
Figure 5. Exicycler™ 384 has superior well-to-well uniformity.

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
Figure 6. Ct values of 10-fold diluted samples show a wide dynamic range of quantification. Fluorescence data from a series of 10-fold dilution of PGK1 DNA(10^{10} copies) amplified using reporter dyes to check on target: FAM/PGK1.


Figure 7. Exicycler™ 384 provide sensitive detection and precise target discrimination down to 2-fold differences Fluorescence data from a series of 2-fold dilution of CSF2 DNA(10^8 copies) amplified using reporter dyes to check one target: TET/CSF2.

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Exicycler™ 384 Real-Time Quantitative Thermal Block

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Physical Specifications	
Dimension (cm)	35.5(W) x 54(D) x 47(H)
Weight (kg)	41
Sample Capacity/ Size	384-well plate
Sample Volume (μl)	5~20 (recommended 10 μl)
Power Consumption	100~240 VAC, 50/60 Hz, 800VA Max
Optical Specifications	
Light Source	Short arc lamp (120W)
Sensor	16-bit 2D CCD
Excitation Filter / Emission Filter	5 Sets
Operating Specifications	
Method of Heating / Cooling	Peltier element
Temperature Range	4.0~99.9°C
Temperature Accuracy	± 0.3°C
Temperature Uniformity	± 0.3°C
Ramping Rate	Max. 4.5°C /sec
Lid Temperature	90~120°C
Gradient Range	20~95°C
Temperature Increment Range	0.1~2.0°C
Time Increment Range	1~60 sec
Ramp Rate Control Range	1~100%
Operating Temperature	15~35°C
Operating Humidity	20~80%, no condensation
Computer Specifications	
Operating System	Windows 7&10 (32/64-bit OS)
Processor Speed	Intel Dual Core E2160 (1.8 GHz) or higher
Memory	1 GB or higher
Communication Port	USB 2.0 high speed
Screen Resolution	1280X1024 or higher

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- Cancer mutation & expression study
- Pathogen study, genotyping, drug-resistance study for Infectious diseases
- Forensic science
- Personalized medicine
- Precision medicine
- Food safety, and bio-defense

Exicycler™ HD Ultimate Multiplex Real-Time PCR Platform

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Physical Specifications	
Dimension (cm)	35.5(W) x 54(D) x 47(H)
Weight (kg)	41
Sample Capacity/ Size	HD-qPCR panel (max 1,536 well)
Sample Volume (μl)	2
Power Consumption	100-240 VAC, 50/60 Hz, 800VA Max
Optical Specifications	
Light Source	Short arc lamp (120W)
Sensor	16-bit 2D CCD
Excitation Filter / Emission Filter	5 Sets
Operating Specifications	
Method of Heating / Cooling	Peltier element
Temperature Range	4.0~99.9°C
Temperature Accuracy	± 0.3°C
Temperature Uniformity	± 0.3°C
Ramping Rate	Max. 4.5°C /sec
Lid Temperature	90~120°C
Temperature Increment Range	0.1 - 9.9°C
Time Increment Range	1~60 sec
Ramp Rate Control Range	1~100%
Operating Temperature	15 ~ 35°C
Operating Humidity	20 ~ 80%, no condensation
Computer Specifications	
Operating System	Windows 7 이상
Processor Speed	Intel Dual Core E2160 (1.8 GHz) or higher
Memory	1 GB or higher
Communication Port	USB 2.0 high speed
Screen Resolution	1280X1024 or higher

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Operating Specifications

Servo Motor	0 ~ 3,000 RPM
Pressure Sensor	-101 kPa ~ 0 kPa
Oil-less Rocking Piston Pump	Max Vacuum: 7,500 mmHg, Flow Rate: 70 lpm
Operating Temperature	15~35°C
Operating Humidity	20~80 %, no condensation
Operating System	Built-in
User Interface Display	4.3 inch 480 x 272 TFT LCD With Touch Screen

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
Physical Specifications		
Dimension (cm)	218.2(W) x 97(D) x 195.5(H) (Fold state of the Signal Tower Lamp)	
	218.2(W) x 97(D) x 220(H) (Unfold state of the Signal Tower Lamp)	
Weight (kg)	650	
Operating Temperature	15 ~ 35℃	
Operating Humidity	20 ~ 80%, no condensation	
Input Voltage	110-120 VAC or 220-240 VAC (Voltage Selection Switch on the back)	
Frequency	50/60 Hz	
Power	Max. 9A/12A	
PC Interface	TCP/IP(x2), USB-B(x2), USB-A(x2), IEEE1394(x2)	
Computer specifications (x2)		
Computer	For Control and Analysis	For Analysis
Communication Port	USB 2.0 high speed port (Minimum 6 ports) IEEE 1394 port (Minimum 2 ports)	USB 2.0 high speed port (Minimum 4 ports)
Operating System	Windows 7 이상	
Processor Speed	Intel I5 (2.8GHz) or higher	
Memory	4 GB RAM or higher	
Screen Resolution	1280 x 800 screen resolution	
network	TCP/IP	
Hard Disk Drive	500 GB minimum	
External Air compressor specifications (Voltage Options)		
Dimension(cm)	25(W) x 50(D) x 55(H)	25(W) x 50(D) x 55(H)
Voltage/ Frequency	220-240 VAC, 50/60 Hz	110-120 VAC, 50/60 Hz
Current Drawn(50/60Hz)	2.6/2.8A	5.1/5.7A
Air Tank Capacity	10 L	10 L
Fuses	250V, F10AL, 2ea	250V, F10AL, 2ea

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03. Protein Synthesis and Purification

<i>ExiProgen™</i>	436
-------------------------	-----

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CAT protein (His-tagged protein)'s expression/purification with ExiProgen™

M: AccuLadder™ Protein Size Marker (Low) (Cat. No. D-2020, Bioneer)

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His-tagged proteins expression/purification with ExiProgen™

M: AccuLadder™ Protein Size Marker (Low) (Cat. No. D-2020, Bioneer)

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Lanes 1, 3, 5, 7, 9, 11, 13, 15 were extracted with 1×10^6 cells of cultured HeLa cell and lanes 2, 4, 6, 8, 10, 12, 14, 16 were extracted with ddH₂O as a negative control in DNA extraction. Note all samples have similar yields. Purity was also tested and was consistently between 1.8 and 2.0 (not shown).

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- Enzyme engineering
- Protein evolution
- Synthetic biology
- Bio-energy R&D
- Protein drug R&D

○ Available Experiments with Bioneer's Kits

- Protein Synthesis/Purification Kit
- Genomic DNA Extraction Kit
- Total RNA Extraction Kit
- Viral DNA/RNA Extraction Kit
- Fragment DNA Purification Kit

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카탈로그 번호	제품명
A-5041	ExiProgen™
카탈로그 번호	Accessories
A-5041-A	Accessories Set for ExiProgen™ (A2,3,4,5,6,7,10)
A-5041-A1	Multi Puncher (option)
A-5041-A2	Setup Tray
A-5041-A3	Disposable Tip Rack
A-5041-A4	Elution Tube Rack
A-5041-A5	Reaction Block (For Protein synthesis)
A-5041-A6	Waste Tray
A-5041-A7	Hole Puncher (6-hole)
A-5041-A9	AC Adapter for ExiProgen™, ExiPrep™16 Dx
A-5041-A10	Contamination Shield
A-5041-A12	Protection Cover Separation Tool
카탈로그 번호	Protein Synthesis and Purification Kits for ExiProgen™
K-7300	ExiProgen™ EC Protein Synthesis Kit, 16 rxns
K-7301	ExiProgen™ EC Protein Synthesis Kit, 32 rxns
K-7302	ExiProgen™ EC Protein Synthesis Kit, 96 rxns
K-7310	ExiProgen™ EC-Maxi Protein Synthesis Kit, 8 rxns
K-7320	ExiProgen™ EC-Tagfree Protein Synthesis Kit, 8 rxns
K-7330	ExiProgen™ EC-Disulfide Protein Synthesis Kit, 8 rxns
K-7220	ExiProgen™ His-tagged Protein Purification Kit, 16 rxns
K-7221	ExiProgen™ His-tagged Protein Purification Kit, 32 rxns
K-7240	ExiProgen™ Dialysis Kit, 16 rxns
KA-3001	ExiProgen™ Consumable SET, 16 rxns
카탈로그 번호	Nucleic Acid Extraction Kits for ExiProgen™
K-3225	ExiPrep™ Tissue Genomic DNA Kit
K-3200-CB	ExiPrep™ Beef Genomic DNA Kit
K-3200-CR	ExiPrep™ Rice Genomic DNA Kit
K-4211	ExiPrep™ Plus Blood Genomic DNA Kit
K-4214	ExiPrep™ Plus Bacteria Genomic DNA Kit
K-4215	ExiPrep™ Plus Plant Genomic DNA Kit
K-4217	ExiPrep™ Plus Seed Genomic DNA Kit
K-4241	ExiPrep™ Plus Total RNA Kit
K-4244	ExiPrep™ Plus Plant total RNA Kit
K-4271	ExiPrep™ Plus Viral DNA/RNA Kit

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
○ Genomic DNA extraction from cultured mammalian cell (HeLa)


Figure 1. Lanes 1, 3, 5, 7, 9, 11, 13, 15 were extracted with 1×10^6 cells of cultured HeLa cell and lanes 2, 4, 6, 8, 10, 12, 14, 16 were extracted with ddH₂O as a negative control in DNA extraction. Note that all the samples have similar yields.

[Sample volume: 1×10^6 cells, Average yield: 4-8 µg]

○ RNA extraction


Figure 2. Comparison of Total RNA extracted from HeLa cell ($\sim 4 \times 10^6$) using ExiPrep™ Plus Total RNA Kit (Bioneer, automatic) and a competitor's kit (Competitor Q, manual/single column). DNase was not treated.

Lane 1-3: Bioneer ExiPrep™ Plus Total RNA Kit, mean yield; 38.6 µg, mean purity; 2.06

Lane 4-5: Competitor Q Total RNA Extraction Kit, mean yield; 39.7 µg, mean purity; 2.05

Lane 6-8: Bioneer ExiPrep™ Plus Total RNA Kit, mean yield; 46.1 µg, mean purity; 2.03

Lane 9-10: Competitor Q Total RNA Extraction Kit, mean yield; 49.7 µg, mean purity; 2.06

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- Gene Expression Study
- Genetic Engineering
- GMP & Food poisoning test
- Biological terror detection

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Dimension (cm)	32(W) x 53.5(D) x 48.7(H)
Weight (kg)	22
Temperature Range	15~35℃
Humidity Range	20~80%, no condensation
Operating System	Stand alone (Built-in)
User Interface	320 x 240 touch screen graphic LCD
Electrical (Voltage/Frequency)	100-240 V, 50/60 Hz
UV Sterilization	15-minute cycle
Communications	TCP/IP
Heating Temperature	40~90℃

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

Operating Specifications

Magnet & Heating Block	30~90°C
Temperature Controlled Block	4~90°C
Operating Temperature	15~35°C
Operating Humidity	20~80%, no condensation
Operating System	Built-in
Communication	TCP/IP
User Interface Display	13.3 inch Full HD (1920 x 1080), Touch
Data Storage	USB 3.0 (Front x2), USB 2.0 (Rear x2)
Output Display	HDMI Port (x1)

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


Mini-scale nucleic acid extraction


Midi-scale nucleic acid extraction

Three sets of rods and covers for 96 samples

Two sets of rods and covers for 64 samples

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Lane 2: *ExiPrep*™96 Lite

Lane 3: Company Q


Figure 1. Animal tissue genomic DNA extraction using *MagListo*™ Genomic DNA Extraction Kit. DNA extraction from 30 mg of Beef.


Figure 2. Blood genomic DNA extraction using *MagListo*™ Genomic DNA Extraction Kit. DNA extraction from 200 μl of human whole blood.


Figure 3. Cultured cell genomic DNA extraction using *MagListo*™ Genomic DNA Extraction Kit. DNA extraction from 1×10⁶ of HeLa cell.

Plasmid DNA Extraction


Figure 4. Plasmid DNA extraction using *ExiPrep*™96 Lite

Lane 1-12: *MagListo*™ 5M Plasmid DNA Extraction Kit with *ExiPrep*™96 Lite

Lane 13-14: *AccuPrep*® Plasmid Extraction Kit

Average yield: 8 μg of pBlueScript plasmid in DH5α (OD₆₀₀=2.7) using *ExiPrep*™96 Lite

ExiPrep™96 Lite

Application on cfDNA Extraction (*MagListo*™ cfDNA Extraction Kit)


Figure 5. Cell-free DNA extracted using *MagListo*™ cfDNA Extraction kit and *ExiPrep*™96 Lite. Cell-free DNA was isolated from various volume of plasma and serum using *ExiPrep*™ 96 Lite DNA was visualized using Agilent 5200 Fragment Analyzer System.

Application on FFPE DNA Extraction



Figure 6. FFPE DNA extracted by using Manual method and *ExiPrep*™96 Lite.

- (A) Gel electrophoresis and NanoDrop measurement of FFPE DNA. FFPE DNA was isolated with *MagListo*™ Genomic DNA Extraction kit by manual method (lane 1-3) and *ExiPrep*™96 Lite (lane 4-11).
- (B) qPCR quantification of FFPE DNA extracted by manual method and *ExiPrep*™96 Lite. Mouse housekeeping gene *Cox6* primers and probe were used for qPCR.

Genomic DNA Extraction

Lane 1: *MagListo*™ (manual type, Bioneer), Lane 2: *ExiPrep*™96 Lite, Lane 3: Company Q, M: 1 kb DNA Ladder (Cat. No. D-1040, Bioneer)


Figure 7. Animal tissue RNA extraction using *MagListo*™ Universal RNA Extraction Kit. RNA extraction from 20mg of *Mus musculus* liver. DNase was not treated.


Figure 8. Cultured cell RNA extraction using *MagListo*™ Universal RNA Extraction Kit. RNA extraction from 1x10⁶ of HeLa cell. DNase was not treated.



Figure 9. Plant RNA extraction using *MagListo*™ Universal RNA Extraction Kit. RNA extraction from 100mg of *Brassica oleracea* var. *italica*. DNase was not treated.

Protein Purification

A Purification of AcGFP using Ni-NTA magnetic beads

M: Protein Size Marker (Cat. No. D-2010, Bioneer)


Lane 1-17: Purification samples with DNA

Lane NC: Purification sample without DNA

Lane EP: Expression sample with DNA

Lane EN: Expression sample without DNA

B Purification of scFv using Protein L magnetic beads

M: Protein Size Marker (Cat. No. D-2010, Bioneer)


Lane 1-46: Purification samples with DNA

Lane NC: Purification sample without DNA

Lane EP: Expression sample with DNA

Lane EN: Expression sample without DNA


C Purification of IgG using Protein G magnetic beads

M: Protein Size Marker (Cat. No. D-2010, Bioneer)

Lane 1-16: Purification samples

Lane L: Loading sample

M: Protein Size Marker (Cat. No. D-2010, Bioneer)

Lane 1-14: Purification samples

Lane L: Loading sample

Figure 6. SDS-PAGE analysis of purified protein

After purification, various proteins were loaded onto a 12% SDS-PAGE and Stained with coomassie Brilliant Blue. Average concentration of purified proteins;

(A) AcGFP: 18 µg/rxn, (B) scFv: 15 µg/rxn, (C-1) IgG: 62 µg/rxn, (C-2) IgG: 30 µg/rxn,

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
Physical Specifications	
Dimension (cm)	40(W) x 57.5(D) x 54.6(H)
Weight (kg)	46
Voltage / Frequency	100-240V~, 50/60Hz
Power	500 VA Max (Fuse: 250V, F6:3AL, 2ea)
Operating Specifications	
Heating Block	30~90°C
Temperature Controlled Block	4~90°C
Operating Temperature	15~35°C
Operating Humidity	20~80%, no condensation
Operating System	Built-in
Communication	TCP/IP
User Interface Display	7-inch Touch Screen
Data Storage	USB 2.0 (x2)

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

Dimensions (cm)	30(W) x 23(D) x 28(H)
Weight (kg)	6.5
No. of wells	16 wells
Diameter of wells	8 mm

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Dimension (cm)	30(W)X42(D)X18(H)
Power	400W
Power Supply	220V / 50 Hz
Travel	35 mm
Vibration Frequency Setting	10~2,100 minute
Sample Capacity	1~2 ml
Sample Size	≤ 8 mm
Output Particle Size	~ 3 μm
Typical Crushing Time	2 min
Weight(kg)	~ 20
Number of Grinding Platforms	2
Grinding Tank Type	Milling tank with screw-on screw cap 1.5-50 ml, optional
Grinding Kit Material	Hard steel, Teflon nylon 0.5 / 1.5mL / 2mL (centrifuge tube / PCR tube)
Grinding Kit Size	25 ml / 50 ml (grinding tank)
Grinding Ball Material	Agate, stainless steel, zirconia, tungsten carbide, ceramic
Grinding Adapter	Adapter 2 × 48 holes; 5 ml adapter 2 × 24 holes.
Crushing Time Setting	Digital display 0 seconds-9,999 minutes

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- 핵산 확인: Cloning, PCR, gel cutting, purification, sequencing
- 단백질 확인: SDS-PAGE, western blot X-ray film

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
Dimension (cm)	18.5(W) x 22(D) x 3(H)
Viewing Surface Dimension (cm)	12(W) x 18(D)
LED Wavelength (nm)	470 nm
Lid / Filter	Amber (580 nm)
LED Life (hours)	> 30,000
Certifications	CE/ETL
Input Voltage	100-240 VAC
Input Current	2.0 A
Automatic Power-off	5 min
Unit Weight (kg)	2.27
The number of LED lamps	36 blue-light LED lamps/ 24 white-light LED lamps
Package Set	Gel-cutting knife, replacement blade, mini darkroom

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DUALED Blue/White Transilluminator

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번호	구성
①	DUALED Blue/White Transilluminator base
②	Amber filter cover
③	Blue uniform plate
④	White uniform plate
⑤	Power cord
⑥	Gel-cutting knife and replacement blade
⑦	Mini darkroom

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Blue Light Mode		White Light Mode
Agarose Gel Cutting	Agarose Gel Imaging	SDS-PAGE Gel Imaging

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DUALLED Blue/White Transilluminator FAQs

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
Figure 1. Cell Growth Curve.

Table.1 Culture Conditions.

HT-MegaGrow®	37°C, 1,200 rpm, no aeration, oxygenation (No delay, 1-minute interval, 1 second on)
Plate	96 Deep well Plate, round-bottom, square-wells, plate lids (Cat. No.90061, Bioneer)
Media	500 µl/well of 1x Terrific Broth _ Ampicillin (100 µg/ml)
Strain	<i>E. coli</i> , DH5a carrying pBluescript (no insert)
Innoculation	5 µl pre-culture added to each wells with multichannel pipettor
Growth Time	12 hours

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■ Fully Automated Spin-Mix-Spin Technology

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
1st spin: Spin down sample
Vortexing of samples
2nd spin: Spin down mixed sample

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Dimensions (cm)	19(W) x 12.5(D) x 23.5(H)
Weight (kg)	2.7
Spin Regulation	1,000~3,500 rpm (increment 100 rpm)
Spin Timer	1 sec~30 min
Vortexing Strength	Soft, Medium, Hard
SMS- Cycle Regulation	1~999 cycles
Power Supply	AC 24 V, 1,250 mA

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