

[Cat. No.] Please refer to the **Ordering Information**

Introduction

RNA-Guided Endonucleases (RGENs) are ready-to-use Cas9 recombinant proteins for your genome editing experiments. *AccuTool™* Cas9 recombinant proteins are recombinant *Streptococcus pyogenes* Cas9 (WT) protein, purified from *E. coli*. The Cas9 recombinant proteins and custom-designed guide RNA (aRGENs) form a stable ribonucleoprotein (RNP) complex, which can be used on cells that are difficult to be transfected like primary cells. The Cas9 Nickase (D10A) generates a DNA nick rather than a double-strand breaks (DSBs), and the Cas9 Dead (D10A/H840A) has no endonuclease activity but remains its ability to bind onto the target sequences. Sniper Cas9 improves specificity and on-target activity by reducing off-targeting. The reagents can be used for *in vitro* and *in vivo* without additional processing.




Figure 1. Cas9 recombinant proteins

Applications

- Genome editing
- Drug discovery: CRISPR library screening, target validation
- Bioprocessing: Cell line engineering
- Agriculture: Plant breeding

Components

Components	Amount
Recombinant SpCas9 (WT or Nickase or Dead or Sniper) protein	50 µg
10X Reaction buffer	1 ml
1X Dilution buffer	1 ml

* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

Specifications

- *AccuTool™* Cas9 recombinant proteins are recombinant *Streptococcus pyogenes* Cas9 (WT) protein, purified from *E. coli*.

- *AccuTool™* Cas9 recombinant proteins concentration: 1 mg/ml.

Buffer Composition

10X Reaction buffer	1 M NaCl, 500 mM Tris-HCl, 100 mM MgCl ₂ , 1 mg/ml BSA, pH 7.9
---------------------	---

Storage Buffer

Cas9 recombinant protein is supplied in 10 mM Tris-HCl, 300 mM NaCl, 0.1 M EDTA, 1 mM DTT, 50% glycerol, and stabilizer, pH 7.4.

Quality Assurance

- Protein purity: >95% by SDS-PAGE with Coomassie Blue Staining.
- RNase Activity: No
- DNase Activity: No
- Protease Activity: No

Storage

- *AccuTool™* Cas9 recombinant proteins 50 µg (1 mg/ml)/vial are delivered at -20°C.
- Store the appropriate aliquot of proteins at or below -70°C. Do not store in a frost-free freezer.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols.

Ordering Information

Description	Cat. No.	
<i>AccuTool™</i> Recombinant SpCas9 WT protein	50 µg 50 µg x 2 ea 50 µg x 5 ea	ATS-0010 ATS-0011 ATS-0012
<i>AccuTool™</i> Recombinant SpCas9 Nickase (D10A) protein	50 µg 50 µg x 2 ea 50 µg x 5 ea	ATS-0016 ATS-0017 ATS-0018
<i>AccuTool™</i> Recombinant SpCas9 dead (D10A/H840A) protein	50 µg 50 µg x 2 ea 50 µg x 5 ea	ATS-0019 ATS-0020 ATS-0021
<i>AccuTool™</i> Recombinant Sniper Cas9 protein	50 µg 50 µg x 2 ea 50 µg x 5 ea	ATS-0025 ATS-0026 ATS-0027

Copyright 2021 BIONEER Corporation. All Rights Reserved.

Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements, and other changes to its products, services, specifications, or product descriptions at any time without notice.

Explanation of Symbols



Experimental Procedures

Steps		Procedure Details																					
Digestion of target sequence <i>in vitro</i> using CRISPR sgRNAs																							
1	 Preparation of reaction mixture	1. Set up the reaction mixture as below. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;">Components</th> <th style="width: 30%;">Amount</th> <th style="width: 40%;"></th> </tr> </thead> <tbody> <tr> <td>Cas9 Protein</td> <td>500 ng</td> <td>(100-1,000 ng)</td> </tr> <tr> <td>CRISPR sgRNA</td> <td>250 ng</td> <td>(100-700 ng)</td> </tr> <tr> <td>Targeting substrate</td> <td>100-150 ng</td> <td>PCR product</td> </tr> <tr> <td></td> <td>80 ng</td> <td>(Plasmid)</td> </tr> <tr> <td>10X Reaction buffer</td> <td>1 µl</td> <td></td> </tr> <tr> <td>D.W.</td> <td>To 10 µl</td> <td></td> </tr> </tbody> </table>	Components	Amount		Cas9 Protein	500 ng	(100-1,000 ng)	CRISPR sgRNA	250 ng	(100-700 ng)	Targeting substrate	100-150 ng	PCR product		80 ng	(Plasmid)	10X Reaction buffer	1 µl		D.W.	To 10 µl	
Components	Amount																						
Cas9 Protein	500 ng	(100-1,000 ng)																					
CRISPR sgRNA	250 ng	(100-700 ng)																					
Targeting substrate	100-150 ng	PCR product																					
	80 ng	(Plasmid)																					
10X Reaction buffer	1 µl																						
D.W.	To 10 µl																						
2	 Incubation of reaction mixture	2. Incubate the reaction mixture at 37°C for 1 hr. 3. Add 4 µg of RNase and incubate for 15 min at 37°C. 4. Add 1 µl of STOP solution to the reaction mixture and incubate for 15 min at 37°C. [STOP solution: 30% glycerol, 1.2% SDS, 250 mM EDTA (pH 8.0)] 5. Analyze on 2% agarose gel.																					
General guidelines for the application of CRISPR RNP to cultured cells by lipid-based transfection (Lipofection)																							
<ul style="list-style-type: none"> The amounts of the reagents given in the protocol below are for one well of a 24-well plate. For other reaction formats, scale the amounts of reagents up or down accordingly. 																							
1	 Preparation of transfection reagent	1. Add Cas9 RNP complex (0.5 µg of Cas9 Nuclease and 250 ng of CRISPR sgRNA) to 50 µl Opti-MEM I Reduce Serum Medium. 2. In a separate tube, dilute the transfection reagent by adding 4 µl of the Lipofectamine 2000 transfection reagent to 50 µl of Opti-MEM I Reduce Serum Medium and mix gently. 3. Incubate for 5 min at room temperature. 4. Add the diluted transfection reagent to the tube containing Cas9 protein/gRNA RNP complexes and mix gently. 5. Incubate at room temperature for 20 min to allow the formation of Cas9/gRNA-lipid complexes.																					
2	 Cell transfection	6. Add the Cas9/gRNA-lipid complexes to the 1x10 ⁵ NIH3T3 cells to be transfected. Swirl the plates gently to allow the mixing of the transfection mixture with the medium. 7. Incubate the plate at 37°C in a humidified CO ₂ incubator in a cell culture incubator for 2-3 days, and proceed with sample assay to determine the genome editing efficiency by T7E1 assay or targeted deep sequencing.																					
	 Option	Refer to the electroporation method for highly efficient transfection only.																					