AccuRapid™ TA Cloning Kit

Cat. No. K-7170

AccuRapid[™] TA Cloning Kit

User Guide

K-7170

Version No.: 1 (2022-02-23)

Please read all the information in booklet before using the unit

BIONEER Corporation

**8-11, Munpyeongseo-ro, Daedeok-gu, Daejeon
34302, Republic of Korea**

Tel: 1588-9788

Email: sales@bioneer.co.kr

www.bioneer.com

Intended Use

AccuRapid[™] TA Cloning Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Safety Warning and Precaution

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information, please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to BIONEER within 30 days.

Quality Management System ISO 9001 Certified

Every aspect of our quality management system from product development, production to quality assurance and supplier qualification meets the world-class standards.

Trademark

AccuRapid[™] is a trademark of BIONEER Corporation.

Copyright

Copyright 2022 BIONEER Corporation. All Rights Reserved.

Notice

BIONEER Corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Contents

| | |
|---|-----------|
| Product Information | 1 |
| Components | 1 |
| Storage | 1 |
| Specifications | 1 |
| Introduction | 2 |
| Product Description | 2 |
| pBHA-T Vector | 3 |
| Experimental Procedures | 5 |
| Materials to be Prepared by Users | 5 |
| Preparation of Insert | 6 |
| Ligation Reaction | 7 |
| Transformation | 8 |
| Results Analysis | 9 |
| Troubleshooting | 11 |
| Ordering Information | 13 |
| Related Products | 13 |
| Explanation of Symbols | 14 |

Product Information

Components

| Components | Amount |
|--|-----------------|
| | K-7170 (20 rxn) |
| pBHA-T vector (25 ng/μl) | 40 μl |
| Control Insert (70 ng/μl) | 20 μl |
| <i>AccuRapid</i> [™] 2X Reaction buffer | 100 μl |
| T4 DNA Ligase (200 U/μl) | 20 μl |

* **Unit definition:** One weiss unit (200 U) is defined as the amount of enzyme required to ligate 90% of *Hind* III fragments of 1 μg of lambda DNA in total volume of 20 μl at 16°C for 30 min.

Storage

All components of *AccuRapid*[™] TA Cloning Kit should be stored at -20°C. To minimize the decrease in cloning efficiency, avoid repeated freeze and thaw cycles of pBHA-T vector and *AccuRapid*[™] 2X Reaction buffer. Therefore, it is recommended to aliquot the solution and store each tube for single reaction.

Specifications

| | <i>AccuRapid</i> [™] TA Cloning Kit |
|-------------------------------|--|
| Reaction time | 15 min |
| Reaction temperature | 25°C |
| Need for PCR/Gel purification | Recommendation |
| Selection marker | Ampicillin |

Introduction

Product Description

TA Cloning is a method that directly clones PCR product by ligating deoxyriboadenosine (dA)-tailed DNA at the 3' end, tailed by *Taq* DNA Polymerase, to a vector with deoxyribothymidine (dT) added to its 3' end. *Taq* DNA Polymerase and other non-proofreading DNA polymerases have terminal transferase activity, which adds dA to the 3' end of the PCR product, resulting in an A-overhang structure. These PCR products can be directly inserted into a pre-linearized T-vector (pBHA-T vector) through ligation. It is possible to perform cloning with various sizes of inserts within 15 min of ligation reaction, and ampicillin or blue/white screening are possible through the ampicillin selection marker and *LacZ* gene in the vector. Therefore, TA cloning can be applied without using primers containing restriction enzyme sites or restriction enzyme treatment, so it is particularly useful when the restriction enzyme sites are limited. In addition, there are multiple restriction enzyme sites around the cloning site in the vector, making it easy to analyze the results or clone onto other vectors. BIONEER's *AccuRapid*™ TA Cloning Kit is a product that can easily and quickly clone the PCR product amplified with a DNA polymerase such as *Taq* DNA Polymerase that can add dA to the 3' end into a linearized pBHA-T vector.

Step 1: Preparation of insert DNA


Figure 1. Schematic protocol for *AccuRapid*™ TA Cloning Kit.

pBHA-T Vector

Features

1. *LacZ* α gene

This gene encodes amino acids 1-146 of β -galactosidase. By complementing the ω -fragment of the host *lacZ* gene, it presents β -galactosidase activity, enabling blue/white screening.

2. Ampicillin resistance gene

It contains an ampicillin resistance gene and is used for transformant selection.

3. M13 Forward and Reverse primer site

It is used to analyze the sequences to confirm if the gene of interest is inserted in the plasmid.


4. T7 promoter and SP6 promoter

It is used to analyze sequences to confirm if the gene of interest is inserted in the plasmid or to transcribe RNA from inserted sequence *in vitro*.

5. *Bam*H I restriction site

It is used to digest the inserted DNA from the recombinant plasmid.

Vector Map

| pBHA-T vector | Position (Bases) |
|-----------------------------------|------------------|
| M13 Reverse primer site | 1-17 |
| <i>LacZα</i> gene | 13-441 |
| SP6 promoter | 35-53 |
| T7 promoter | 236-254 |
| M13 Forward primer site | 261-277 |
| Ampicillin resistance gene | 1,003-1,863 |
| pUC <i>ori</i> (high copy number) | 2,034-2,622 |

Experimental Procedures

Materials to be Prepared by Users

- *AccuPower*[®] *Taq* PCR PreMix (Cat. No. K-2601) or *AccuPower*[®] *ProFi Taq* PCR PreMix (Cat. No. K-2631)
- *AccuPower*[®] PCR PreMix (Cat. No. K-2012)
- *AccuPrep*[®] PCR/Gel Purification Kit (Cat. No. K-3038)
- Thermal cycler (*AllInOneCycler*[™] PCR system, A-2041)
- Chemically or electrocompetent cell
- IPTG (Cat. No. C-8001)
- X-Gal (Cat. No. C-8002)
- LB agar plate with ampicillin
- LB broth or SOC media
- M13 Forward primer, M13 Reverse primer

Preparation of Insert

1. Amplify insert from template DNA using *AccuPower® ProFi Taq* PCR PreMix (Cat. No. K-2631) and *AllInOneCycler™* PCR system (Cat. No. A-2041).

*** Note:** For the amplification of insert DNA, DNA polymerase without 3'→5' exonuclease activity must be used. In addition, *AccuPower®* PCR PreMix (Cat. No. K-2012) and *AccuPower® Taq* PCR PreMix (Cat. No. K-2601) can be used, but if you need proofreading function, we recommend using *AccuPower® ProFi Taq* PCR PreMix.

2. Add template DNA, primers, and nuclease-free water into *AccuPower® ProFi Taq* PCR PreMix tubes to make a total volume of 20 µl or 50 µl. Do not calculate the dried pellet.

• Preparation of reaction mixture

| Components | 20 µl reaction | 50 µl reaction |
|-----------------------------|--------------------|--------------------|
| Template DNA (1-500 ng) | Variable (1-10 µl) | Variable (1-25 µl) |
| Forward primer (10 pmol/µl) | 0.5-2 µl | 1-5 µl |
| Reverse primer (10 pmol/µl) | 0.5-2 µl | 1-5 µl |
| Nuclease-free water | Variable | Variable |
| Total volume | 20 µl | 50 µl |

3. Dissolve the vacuum-dried pellet by vortexing, and briefly spin down. Perform the reaction under the following conditions.

| Step | Temperature | Time | Cycles |
|------------------|-------------|-----------|-----------|
| Pre-denaturation | 95°C | 5 min | 1 cycle |
| Denaturation | 95°C | 15-20 sec | 30 cycles |
| Annealing | 45-65°C | 15-30 sec | |
| Extension | 72°C | 1 min/kb | |
| Final extension | 72°C | 5-10 min* | 1 cycle |

*** Note:** For efficient A-tailing of PCR products, it is recommended to proceed with final extension at 72°C for more than 5 min.

4. After electrophoresis of the amplified insert, gel elution is performed with *AccuPrep®* PCR/Gel Purification Kit (Cat. No. K-3038).

*** Note:** Using an unpurified PCR product can significantly decrease the cloning efficiency, so it is recommended to purify the target DNA band through gel purification.

* **Note:** Commercially available elution buffer contains EDTA, which can reduce ligation efficiency. It is recommended to use the elution buffer (EA buffer) without EDTA of BIONEER's *AccuPrep*® PCR/Gel Purification Kit or distilled water.

Ligation Reaction

1. Thaw pBHA-T vector and *AccuRapid*™ 2X Reaction buffer on ice before ligation reaction.
2. Prepare reaction mixture for ligation reaction as follows.

| Components | Positive | Sample |
|---------------------------------------|----------|-----------|
| pBHA-T vector (25 ng/μl) | 2 μl | 2 μl |
| PCR product | - | 0.5-2 μl* |
| Control Insert (70 ng/μl) | 1 μl | - |
| <i>AccuRapid</i> ™ 2X Reaction buffer | 5 μl | 5 μl |
| T4 DNA Ligase (200 U/μl) | 1 μl | 1 μl |
| Nuclease-free water | 1 μl | Variable |
| Total volume | 10 μl | 10 μl |

* **Note:** Calculate the required amount of PCR product using the following formula. In general, 0.5-2 μl of PCR product is used, and it is recommended to add in a molar ratio of 5 : 1 (insert: vector) or higher for a more efficient ligation reaction.

$$\text{PCR product (ng)} = \frac{\text{pBHA-T vector (ng)} \times \text{size of PCR product (bp)}}{\text{Size of pBHA-T vector (bp)}} \times \text{molar ratio of } \frac{\text{PCR product}}{\text{pBHA-T vector}}$$

* **Example:** The required PCR product concentration, when performing with 1,000 bp of PCR product and with molar ratio (insert : vector) = 5 : 1.

$$\text{PCR product (ng)} = \frac{50 \text{ ng} \times 1,000 \text{ bp}}{\text{around } 3,000 \text{ bp}} \times \frac{5}{1} = \sim 83 \text{ ng}$$

3. Mix the reaction mixture by tapping and briefly spin down.
4. Incubate the reaction mixture at 25°C for at least 15 min.

* **Note:** Sufficient colonies can be obtained in only 15 min of reaction, but if you want to increase the cloning efficiency, it is recommended to increase the reaction time up to 1 hr.
5. After the reaction, reaction mixture should be stored either on ice or at -20°C until transformation.

Transformation

1. Take competent cells from deep freezer and slowly thaw on ice.
2. Add 2-5 μ l of ligation reaction mixture into 50 μ l of competent cells and incubate on ice for 20 min.
* **Note:** Do not pipetting or vortexing.
3. Heat shock the cells at 42°C for 60 sec. Then, immediately place tubes on ice and incubate 2 min.
4. Add 1 ml of pre-heated LB broth (not included antibiotics) or SOC media in 37°C and incubate at 37°C shaking incubator (200 rpm) for 1 hr.
5. Spread 50 μ l of 0.1 M IPTG (Cat. No. C-8001) and X-Gal (50 mg/ml) (Cat. No. C-8002) onto each LB agar plates containing ampicillin (50-100 μ g/ml) and incubate at 37°C.
6. After the incubation, centrifuge at 3,000 rpm for 5 min, remove the supernatant, add 200 μ l of LB broth, and gently resuspend the cell pellet by tapping.
7. Transfer to preheated LB agar plate (included antibiotics) in step 5, and spread completely.
8. Incubate overnight at 37°C in an incubator.

Results Analysis

1. After transformation, check colonies in the cultured plate.

* **Note:** Colony numbers can be different depending on the size, conditions of competent cells, transformation conditions, etc.

2. Select at least 10 white colonies and add colony, primers, and nuclease-free water into *AccuPower*[®] PCR PreMix (Cat. No. K-2012) tubes to make a total volume of 20 µl or 50 µl. Do not calculate the dried pellet.

* **Note:** Please refer to p.10 of table for primer sequences.

• Preparation of reaction mixture

| Components | 20 µl reaction | 50 µl reaction |
|---------------------------------|----------------|----------------|
| Template DNA | Colony | Colony |
| M13 Forward primer (10 pmol/µl) | 0.5-2 µl | 1-5 µl |
| M13 Reverse primer (10 pmol/µl) | 0.5-2 µl | 1-5 µl |
| Nuclease-free water | Variable | Variable |
| Total volume | 20 µl | 50 µl |

3. Dissolve the vacuum-dried pellet by vortexing, and briefly spin down. Perform the reaction under the following conditions.

| Step | Temperature | Time | Cycles |
|------------------|-------------|----------|--------------|
| Pre-denaturation | 95°C | 5 min | 1 cycle |
| Denaturation | 95°C | 20 sec | 25-30 cycles |
| Annealing | 45-65°C | 20 sec | |
| Extension | 72°C | 1 min/kb | |
| Final extension | 72°C | 3-5 min | 1 cycle |

4. Perform electrophoresis to check whether insert DNA is inserted.

* Expected size of DNA in experimental group: 249 bp[†] + size of insert DNA

† Length between M13 Forward primer and M13 Reverse primer

* Expected size of DNA in control group: 649 bp

5. When performing sequencing, T7 promoter, SP6, M13F(-40), M13R(-40), M13F(-20), and M13R(-20) primers can be used.

| Primer Name | Sequences (5'→3') | Length (bases) |
|-------------|-------------------------------|----------------|
| T7 promoter | TAA TAC GAC TCA CTA TAG GG | 20 |
| SP6 | ATT TAG GTG ACA CTA TAG | 18 |
| M13F(-40) | GTT TTC CCA GTC ACG AC | 17 |
| M13R(-40) | CAG GAA ACA GCT ATG AC | 17 |
| M13F(-20) | GTAAA CGA CGG CCA GT | 17 |
| M13R(-20) | GCG GAT AAC AAT TTC ACA CAG G | 22 |

Troubleshooting

1. No colonies produced

| Cause | Solution |
|--|---|
| Decreased transformation efficiency of competent cells | Check the expiration date and efficiency of competent cells. |
| Incorrect concentration of antibiotics in LB agar plates | The pBHA-T vector of this product uses the ampicillin selection marker. Check the ampicillin concentration in LB agar plates, and it is recommended to use antibiotics produced within 1 month. |

2. Most colonies are blue or light blue.

| Cause | Solution |
|--|---|
| Low expression of the <i>lacZ</i> gene or partial disruption of the <i>lacZ</i> gene by insert DNA | If the length of the insert DNA is less than 300 bp, it is recommended to analyze the light blue colonies as they may also contain insert DNA. |
| Non-formation of A-tailing at the 3' end of the PCR product | DNA polymerase such as <i>Pfu</i> DNA Polymerase with 3'→5' exonuclease activity does not have 3' A-tailing function. Therefore, it is recommended to use <i>Taq</i> DNA Polymerase without 3'→5' exonuclease activity. |
| Long-term storage of PCR products | It is recommended to use immediately the PCR product after amplification and purification before ligation. |
| Low molar ratio of insert DNA and vector | It is recommended to use insert DNA in a molar ratio of 5 times or more for ligation reaction. |
| Inhibited ligation reaction | High salt concentrations may inhibit the ligation reaction. It is recommended to minimize the amount (volume) of PCR product. |

Inappropriate host cell for blue/white screening

Check genotype of host cell whether host bacterial strains have *lacZ* mutation.

3. White colonies do not contain insert DNA.

| Cause | Solution |
|--|---|
| Use of plasmid DNA as template DNA for PCR amplification | If the plasmid DNA used as a template uses the ampicillin antibiotics marker, the template plasmid DNA may be inserted during transformation. After the PCR, remove the plasmid DNA through gel purification or <i>Dpn I</i> treatment. |
| Inclusion of non-specific products in PCR product | After amplifying the insert DNA by PCR, electrophoresis to gel purification of the target band. |
| Formation of satellite colonies | When ampicillin is used as an antibiotic marker, satellite colonies can be formed around large white colonies. Select as large white colonies as possible. |
| Degraded 3' T-overhangs on the vector | Repeated freezing and thawing of the pBHA-T vector can degrade the 3' T overhangs. If the 3' T-overhangs are damaged, it undergoes self-ligation which results in damage to the reading frame of the <i>lacZ</i> gene and thus no activity can be seen. Avoid repeated freeze and thaw cycles and it is recommended to store in small aliquots. |

4. Most colonies are blue or light blue.

| Cause | Solution |
|-------------------------------------|--|
| No X-Gal or IPTG in LB agar plates. | Check whether both IPTG and X-Gal are properly added to in LB agar plates. |

Ordering Information

| Description | | Cat. No |
|----------------------------------|---------|---------|
| <i>AccuRapid™</i> TA Cloning Kit | 20 rxns | K-7170 |

Related Products

| Description | Cat. No |
|--|----------|
| Oligo Synthesis (Primer) Service | S-1001 |
| Standard Sequencing Service | S-3010-1 |
| <i>AccuPower®</i> PCR PreMix & Master Mix | K-2012 |
| <i>AccuPower® Taq</i> PCR PreMix & Master Mix | K-2601 |
| <i>AccuPower® ProFi Taq</i> PCR PreMix & Master Mix | K-2631 |
| <i>AccuPrep®</i> Nano-Plus Plasmid Mini Extraction Kit | K-3111 |
| <i>AccuPrep®</i> PCR/Gel Purification Kit | K-3038 |
| T4 DNA Ligase | E-3061 |
| IPTG | C-8001 |
| X-Gal | C-8002 |
| Agarose, 100 g | C-9100 |
| 50X TAE ,500 ml | C-9004 |
| <i>GreenStar™</i> Nucleic Acid Staining Solution I | C-9036 |
| <i>AllInOneCycler™</i> PCR System | A-2041 |
| <i>Agaro-Power™</i> System | A-7020 |
| DUALED Blue/White Transilluminator | A-6020 |

Explanation of Symbols

| | | | |
|--|---|--|---|
|  <p>Batch Code</p> |  <p>Catalog Number</p> |  <p>Caution</p> |  <p>Consult Instructions For Use</p> |
|  <p>Contains Sufficient for <n> tests</p> |  <p>Manufacturer</p> |  <p>Research Use Only</p> |  <p>Temperature Limitation</p> |
|  <p>Use-by Date</p> | | | |

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BIONEER Corporation - HQ
BIONEER Inc. - USA Branch
BIONEER Corp. - European Branch