

[Cat. No.] Please refer to the **Ordering Information**

### Introduction

AccuPower® PCR PreMix is the powerful technology for convenient and easy performance of DNA amplification. This product contains vacuum-dried components including *Top* DNA Polymerase, dNTPs, reaction buffer, stabilizer, and tracking dye. It simplifies preparation of reaction mixture by adding template DNA and primers without any extra process. After the reaction, samples can be applied directly on agarose gel for analysis.

### Applications

- Conventional PCR
- Primer extension
- TA cloning
- Gene sequencing

### Features & Benefits

- **Stability:** Included stabilizer enables delivery at room temperature and provides increased stability compared to solution-type products.
- **User-friendly:** Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform PCR by adding template DNA and primers.
- **Sensitivity:** Excellent sensitivity and amplification efficiency even with small amounts of DNA.
- **Reproducibility:** Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

### Composition

Composition	Concentration
<i>Top</i> DNA Polymerase	1 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM
Reaction buffer with 1.5 mM MgCl <sub>2</sub>	1X
Stabilizer and tracking dye	O

### Specifications

<i>Top</i> DNA Polymerase	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 10 kb

### Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

### Precautions

- This enzyme is specifically optimized for increasing base incorporation rate by inactivating 5' to 3' exonuclease activity. Therefore, this product is not recommended to use for real-time PCR using hydrolysis probe method.

### Online Resources



Korean



English

Visit our **product page** for additional information and protocols


### Ordering Information

Description		Cat. No.	
0.2 ml thin-wall 8-tube strips with attached cap	96 tubes	20 µl/rxn	K-2012
		50 µl/rxn	K-2013
	480 tubes	20 µl/rxn	K-2016
		50 µl/rxn	K-2017
0.5 ml thin-wall microtubes	100 tubes	50 µl/rxn	K-2011
	thin-wall 96-well	flat plate	10 µl/rxn
		20 µl/rxn	K-2260-4
full-skirted plate		10 µl/rxn	K-2260-2
		20 µl/rxn	K-2260-5
semi-skirted plate		10 µl/rxn	K-2260-3
		20 µl/rxn	K-2260-6
thin-wall 384-well	full-skirted plate	5 µl/rxn	K-2080-1
		10 µl/rxn	K-2080-2
		20 µl/rxn	K-2080-3





### Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

### Explanation of Symbols



**Experimental Procedures**

Steps		Procedure Details																																
1	 <b>Preparation of reaction mixture</b>	<p>1. Add template DNA, primers, and nuclease-free water into <i>AccuPower®</i> PCR PreMix tubes make to a total volume of 20 µl or 50 µl. Do not include the dried pellet.</p> <ul style="list-style-type: none"> <li>Amount of template</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Template DNA</th> <th colspan="2">Amount of template</th> </tr> <tr> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Bacteriophage λ, Plasmid DNA</td> <td>100 fg-200 ng</td> <td>100 fg-500 ng</td> </tr> <tr> <td>Total genomic DNA</td> <td>1-500 ng</td> <td>1 ng-1 µg</td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA</td> <td>Variable (1-10 µl)</td> <td>Variable (1-25 µl)</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td>0.5-2 µl</td> <td>1-5 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td>0.5-2 µl</td> <td>1-5 µl</td> </tr> <tr> <td>Nuclease-free water</td> <td>Variable</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>20 µl</td> <td>50 µl</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried blue pellet by pipetting or vortexing, and briefly spin down.</p>				Template DNA	Amount of template		20 µl reaction	50 µl reaction	Bacteriophage λ, Plasmid DNA	100 fg-200 ng	100 fg-500 ng	Total genomic DNA	1-500 ng	1 ng-1 µg	Components	20 µl reaction	50 µl reaction	Template DNA	Variable (1-10 µl)	Variable (1-25 µl)	Forward primer (10 pmol/µl)	0.5-2 µl	1-5 µl	Reverse primer (10 pmol/µl)	0.5-2 µl	1-5 µl	Nuclease-free water	Variable	Variable	Total volume	20 µl	50 µl
Template DNA	Amount of template																																	
	20 µl reaction	50 µl reaction																																
Bacteriophage λ, Plasmid DNA	100 fg-200 ng	100 fg-500 ng																																
Total genomic DNA	1-500 ng	1 ng-1 µg																																
Components	20 µl reaction	50 µl reaction																																
Template DNA	Variable (1-10 µl)	Variable (1-25 µl)																																
Forward primer (10 pmol/µl)	0.5-2 µl	1-5 µl																																
Reverse primer (10 pmol/µl)	0.5-2 µl	1-5 µl																																
Nuclease-free water	Variable	Variable																																
Total volume	20 µl	50 µl																																
2	 <b>Incubate reactions in a thermal cycler</b>	<p>3. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>20 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>45-65°C</td> <td>20 sec</td> <td>25-35 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>0.5-1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> </tbody> </table>				Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing	45-65°C	20 sec	25-35 cycles	Extension	72°C	0.5-1 min/kb		Final extension	72°C	3-5 min	1 cycle					
Step	Temperature	Time	Cycles																															
Pre-denaturation	95°C	5 min	1 cycle																															
Denaturation	95°C	20 sec																																
Annealing	45-65°C	20 sec	25-35 cycles																															
Extension	72°C	0.5-1 min/kb																																
Final extension	72°C	3-5 min	1 cycle																															
3	 <b>Analyze with gel electrophoresis</b>	<p>4. After the reaction, maintain the reaction mixture at 4-8°C. The samples can be stored at -20°C until use.</p> <p>5. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																																
	 <b>Option</b>	<ul style="list-style-type: none"> <li>If primer's T<sub>m</sub> value is more than 65°C or PCR product size is more than 5 kb, follow the conditions as below.</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>20 sec</td> <td></td> </tr> <tr> <td>Annealing/Extension</td> <td>68°C</td> <td>1 min/kb</td> <td>30-35 cycles</td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> </tbody> </table>				Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing/Extension	68°C	1 min/kb	30-35 cycles	Final extension	72°C	3-5 min	1 cycle									
Step	Temperature	Time	Cycles																															
Pre-denaturation	95°C	5 min	1 cycle																															
Denaturation	95°C	20 sec																																
Annealing/Extension	68°C	1 min/kb	30-35 cycles																															
Final extension	72°C	3-5 min	1 cycle																															