

[Cat. No.]      **K-2263**

## Introduction

AccuPower® RT Master Mix is a ready-to-use mixture for cDNA synthesis including *M-MLV* Reverse Transcriptase, reaction buffer, and RNase inhibitor. It simplifies preparation of reverse transcription reaction mixture by adding template RNA and primers without any extra process. Furthermore, by using the RNase H<sup>+</sup> of *M-MLV* Reverse Transcriptase, template RNA is removed after the cDNA synthesis to minimize the effect of residual template RNA in PCR.

## Applications

- First-strand synthesis of cDNA from RNA molecules (RT)
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

## Features & Benefits

- User-friendly: Reactants are included in a tube, it allows any user simply perform cDNA synthesis by adding template RNA and primers.
- Stability: Included stabilizer and RNase inhibitor provides high resistance to degradation.
- Reproducibility: Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

## Composition

Composition	Concentration
<i>M-MLV</i> Reverse Transcriptase	200 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 μM
5X Reaction buffer	1X
DTT	0.25 mM
Stabilizer	1X

## Specifications

<i>M-MLV</i> Reverse Transcriptase	
DNase activity	No
RNase activity	No
Fragment size	Up to 9 kb

## Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

## Online Resources



Korean



English

Visit our **product page** for additional information and protocols

## Ordering Information

Description	Cat. No.
1 ml of 2X Master Mix solution	1 ml x 1 ea
	K-2263





## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols

LOT Batch Code	Biological Risks	REF Catalog Number	Caution
Consult Instructions For Use	Contains Sufficient for <n> tests	Do not Re-use	Manufacturer
RUO Research Use Only	Temperature Limitation	Use-by Date	

## Experimental Procedures

Steps		Procedure Details																								
1	 <b>Primer annealing</b>	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> <li>Amount of template RNA and primers</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>0.5-1 µg</td> <td>1-2 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td>0.05-0.1 µg</td> <td>0.1-0.2 µg</td> </tr> <tr> <td rowspan="3">Primers</td> <td>Random primer</td> <td>100 pmol</td> <td>200 pmol</td> </tr> <tr> <td>Gene specific primer</td> <td>10-30 pmol</td> <td>20-50 pmol</td> </tr> <tr> <td>Oligo dT</td> <td>100 pmol</td> <td>200 pmol</td> </tr> </tbody> </table> <p>2. Incubate the mixture at 70°C for 5 min and place it on ice directly.</p>	Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	0.5-1 µg	1-2 µg	Poly(A) RNA	0.05-0.1 µg	0.1-0.2 µg	Primers	Random primer	100 pmol	200 pmol	Gene specific primer	10-30 pmol	20-50 pmol	Oligo dT	100 pmol	200 pmol			
Components		20 µl reaction	50 µl reaction																							
Template RNA	Total RNA	0.5-1 µg	1-2 µg																							
	Poly(A) RNA	0.05-0.1 µg	0.1-0.2 µg																							
Primers	Random primer	100 pmol	200 pmol																							
	Gene specific primer	10-30 pmol	20-50 pmol																							
	Oligo dT	100 pmol	200 pmol																							
2	 <b>Thaw reagents</b>	<p>3. Thaw <i>AccuPower</i>® RT Master Mix on ice and mix thoroughly before use. Then, briefly spin down components.</p> <p>4. Dispense appropriate volumes of <i>AccuPower</i>® RT Master Mix into PCR tubes (not provided). Use 10 µl and 25 µl of 2X Master Mix for 20 µl reaction and 50 µl reaction, respectively.</p>																								
3	 <b>cDNA synthesis</b>	<p>5. Transfer the incubated primer annealing mixture to PCR tubes including master mix, and then fill up the nuclease-free water to make a total volume of 20 µl or 50 µl.</p> <p>6. Mix the reaction mixture by pipetting or vortexing, and briefly spin down.</p> <p>7. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>42°C</td> <td>60 min</td> </tr> <tr> <td>Heat inactivation</td> <td>95°C</td> <td>5 min</td> </tr> </tbody> </table> <p>8. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>	Step	Temperature	Time	cDNA synthesis	42°C	60 min	Heat inactivation	95°C	5 min															
Step	Temperature	Time																								
cDNA synthesis	42°C	60 min																								
Heat inactivation	95°C	5 min																								
	 <b>Option</b>	<ul style="list-style-type: none"> <li>If PCR is needed, transfer 2-5 µl of reaction mixture (synthesized cDNA) into <i>AccuPower</i>® PCR PreMix tubes (K-2012, not provided), and perform the reaction under the following conditions.</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>20 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>45-65°C</td> <td>20 sec</td> <td>25-35 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>0.5-1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* <b>Note:</b> For maximum yield and specificity, temperatures and cycling times should be optimized for each new template DNA or primers.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing	45-65°C	20 sec	25-35 cycles	Extension	72°C	0.5-1 min/kb		Final extension	72°C	3-5 min	1 cycle
Step	Temperature	Time	Cycles																							
Pre-denaturation	95°C	5 min	1 cycle																							
Denaturation	95°C	20 sec																								
Annealing	45-65°C	20 sec	25-35 cycles																							
Extension	72°C	0.5-1 min/kb																								
Final extension	72°C	3-5 min	1 cycle																							