

[Cat. No.] Please refer to the **Ordering Information**

Introduction

AccuPower[®] RT PreMix contains vacuum-dried components essential for cDNA synthesis including *M-MLV* Reverse Transcriptase, reaction buffer, and RNase inhibitor. It simplifies preparation of reverse transcription reaction mixture by adding template RNA and primers without any extra process. Furthermore, by using the RNase H⁺ of *M-MLV* Reverse Transcriptase, template RNA is removed after the cDNA synthesis to minimize the effect of residual template RNA in PCR.

Applications

- First-strand synthesis of cDNA from RNA molecules (RT)
- Random priming reaction
- cDNA library construction
- Probe labeling
- mRNA 5'-end mapping by primer extension analysis
- Real-time PCR

Features & Benefits

- User-friendly: Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform cDNA synthesis by adding template RNA and primers.
- Stability: Included stabilizer and RNase inhibitor provides high resistance to degradation.
- Reproducibility: Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

Composition

Composition	Concentration
<i>M-MLV</i> Reverse Transcriptase	200 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM
5X Reaction buffer	1X
DTT	0.25 mM
Stabilizer	1X

Specifications

<i>M-MLV</i> Reverse Transcriptase	
DNase activity	No
RNase activity	No
Fragment size	Up to 9 kb

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols

Ordering Information


Description	Cat. No.
20 µl/rxn	K-2041
96 tubes	
0.2 ml thin-wall 8-tube strips with attached cap	50 µl/rxn (with dye) K-2043
20 µl/rxn	K-2041-B
480 tubes	
0.5 ml thin-wall tubes with attached cap	50 µl/rxn (with dye) K-2043-B
20 µl/rxn (with dye)	K-2040
100 tubes	
50 µl/rxn (with dye)	K-2042
10 µl/rxn	K-2261-1
flat plate	20 µl/rxn K-2261-4
thin-wall 96-well	10 µl/rxn K-2261-2
full-skirted plate	20 µl/rxn K-2261-5
10 µl/rxn	K-2261-3
20 µl/rxn	K-2261-6
5 µl/rxn	K-2082-1
thin-wall 384-well	10 µl/rxn K-2082-2
full-skirted plate	20 µl/rxn K-2082-3

* **Note:** The RT series does not contain tracking dyes for electrophoresis, with the exception of some products.

Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols



Experimental Procedures

Steps		Procedure Details																											
1	 Primer annealing	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> Amount of template RNA and primers <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2" style="text-align: center;">Template RNA</td> <td>Total RNA</td> <td style="text-align: center;">0.5-1 µg</td> <td style="text-align: center;">1-2 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td style="text-align: center;">0.05-0.1 µg</td> <td style="text-align: center;">0.1-0.2 µg</td> </tr> <tr> <td rowspan="3" style="text-align: center;">Primers</td> <td>Oligo dT</td> <td style="text-align: center;">100 pmol</td> <td style="text-align: center;">200 pmol</td> </tr> <tr> <td>Random primer</td> <td style="text-align: center;">100 pmol</td> <td style="text-align: center;">200 pmol</td> </tr> <tr> <td>Gene specific primer</td> <td style="text-align: center;">10-30 pmol</td> <td style="text-align: center;">20-50 pmol</td> </tr> </tbody> </table> <p>2. Incubate the mixture at 70°C for 5 min and place it on ice directly.</p>				Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	0.5-1 µg	1-2 µg	Poly(A) RNA	0.05-0.1 µg	0.1-0.2 µg	Primers	Oligo dT	100 pmol	200 pmol	Random primer	100 pmol	200 pmol	Gene specific primer	10-30 pmol	20-50 pmol			
Components		20 µl reaction	50 µl reaction																										
Template RNA	Total RNA	0.5-1 µg	1-2 µg																										
	Poly(A) RNA	0.05-0.1 µg	0.1-0.2 µg																										
Primers	Oligo dT	100 pmol	200 pmol																										
	Random primer	100 pmol	200 pmol																										
	Gene specific primer	10-30 pmol	20-50 pmol																										
2	 cDNA synthesis	<p>3. Transfer the incubated primer annealing mixture to <i>AccuPower[®] RT PreMix</i> tube, and then fill up the nuclease-free water to make a total volume of 20 µl or 50 µl.</p> <p>4. Dissolve the vacuum-dried pellet by pipetting or vortexing, and briefly spin down.</p> <p>5. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td style="text-align: center;">42°C</td> <td style="text-align: center;">60 min</td> </tr> <tr> <td>Heat inactivation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> </tr> </tbody> </table> <p>6. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>				Step	Temperature	Time	cDNA synthesis	42°C	60 min	Heat inactivation	95°C	5 min															
Step	Temperature	Time																											
cDNA synthesis	42°C	60 min																											
Heat inactivation	95°C	5 min																											
	 Option	<ul style="list-style-type: none"> If PCR is needed, transfer 2-5 µl of reaction mixture (synthesized cDNA) into <i>AccuPower[®] PCR PreMix</i> tubes (K-2012, not provided), and perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">20 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td style="text-align: center;">45-65°C</td> <td style="text-align: center;">20 sec</td> <td style="text-align: center;">25-35 cycles</td> </tr> <tr> <td>Extension</td> <td style="text-align: center;">72°C</td> <td style="text-align: center;">0.5-1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td style="text-align: center;">72°C</td> <td style="text-align: center;">3-5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>* Note: For maximum yield and specificity, temperatures and cycling times should be optimized for each new template DNA or primers.</p>				Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing	45-65°C	20 sec	25-35 cycles	Extension	72°C	0.5-1 min/kb		Final extension	72°C	3-5 min	1 cycle
Step	Temperature	Time	Cycles																										
Pre-denaturation	95°C	5 min	1 cycle																										
Denaturation	95°C	20 sec																											
Annealing	45-65°C	20 sec	25-35 cycles																										
Extension	72°C	0.5-1 min/kb																											
Final extension	72°C	3-5 min	1 cycle																										