

[Cat. No.] **K-2217**

Introduction

AccuPower® RocketPlex™ RT-PCR Master Mix utilizes one-step RT-PCR performing cDNA synthesis and PCR of two or more target genes (up to 10 targets) in a single step, in a single tube. It can effectively synthesize cDNA with complex secondary RNA structures by using RocketScript™ Reverse Transcriptase (RocketScript™ RTase) developed by BIONEER. RocketScript™ RTase originated and engineered from *M-MLV* reverse transcriptase can provide increased thermal stability resulting in effective full length first strand cDNA synthesis. It can be used at various temperatures (42-70°C) allowing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptase. This product is a ready-to-use mixture containing RocketScript™ RTase, reaction buffer, DTT, dNTPs, RNase inhibitor, and HotStart Taq DNA Polymerase. It simplifies preparation of reverse transcription reaction mixture and PCR mixture by adding template RNA and primers without any extra process. After the reaction, samples can be applied directly on agarose gel for analysis due to included tracking dye.




Figure 1. Schematic diagram of 5' UTR of a gene, with complex secondary structure, at three different temperatures.

RocketScript™ Reverse Transcriptase is fully active at 70°C, so it can synthesize complete gene sequences that *M-MLV* and other reverse transcriptase can not synthesize.

Applications

- Multiplex RT-PCR
- Low copy detection
- Gene-expression analysis

Features & Benefits

- Multiplex RT-PCR: Simultaneous cDNA synthesis and amplification up to 10 target genes.
- Specificity: Minimized non-specific amplification and maximized PCR efficiency by using HotStart Taq DNA Polymerase.
- Ease-of-use: Reactants are included in a tube, it allows any user simply perform cDNA synthesis and PCR in one tube by adding template RNA and primers. In addition, tracking dye and sedimentation agents are included in products for electrophoresis.

Composition

Composition	Concentration
RocketScript™ Reverse Transcriptase	200 U
5X Reaction buffer	1X
DTT	0.25 mM
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 μM
RNase inhibitor	1 U
HotStart Taq DNA Polymerase	1 U
Stabilizer and tracking dye	1X

Specifications

HotStart Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 1 kb

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols

Ordering Information

Description	Cat. No.
1 ml of 2X Master Mix solution	1 ml x 1 ea K-2217





Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols

Batch Code	Biological Risks	Catalog Number	Caution
Consult Instructions For Use	Contains Sufficient for <n> tests	Do not Re-use	Manufacturer
Research Use Only	Temperature Limitation	Use-by Date	

Experimental Procedures

Steps		Procedure Details																												
1	 Thaw reagents	<p>1. Thaw <i>AccuPower® RocketPlex™</i> RT-PCR Master Mix on ice and mix thoroughly before use. Then, briefly spin down components.</p> <p>2. Dispense appropriate volumes of <i>AccuPower® RocketPlex™</i> RT-PCR Master Mix into PCR tubes (not provided). Use 10 µl and 25 µl of 2X Master Mix for 20 µl reaction and 50 µl reaction, respectively.</p>																												
2	 Preparation of reaction mixture	<p>3. Add template RNA, primers, and nuclease-free water into PCR tubes containing <i>AccuPower® RocketPlex™</i> RT-PCR Master Mix.</p> <ul style="list-style-type: none"> Amount of template RNA and primers <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Primers</td> <td>Gene specific primer</td> <td>5-10 pmol</td> <td>5-10 pmol</td> </tr> </tbody> </table> <p>4. Mix the reaction mixture by pipetting or vortexing, and briefly spin down.</p>	Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg	Primers	Gene specific primer	5-10 pmol	5-10 pmol													
Components		20 µl reaction	50 µl reaction																											
Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg																											
	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg																											
Primers	Gene specific primer	5-10 pmol	5-10 pmol																											
3	 RT-PCR	<p>5. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>42-70°C</td> <td>10-60 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>10-30 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>50-65°C</td> <td>10-30 sec</td> <td>30 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>1 kb/min</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: Recommended temperature of reverse transcription reaction is 50°C.</p>	Step	Temperature	Time	Cycles	cDNA synthesis	42-70°C	10-60 min	1 cycle	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	10-30 sec		Annealing	50-65°C	10-30 sec	30 cycles	Extension	72°C	1 kb/min		Final extension	72°C	5 min	1 cycle
Step	Temperature	Time	Cycles																											
cDNA synthesis	42-70°C	10-60 min	1 cycle																											
Pre-denaturation	95°C	5 min	1 cycle																											
Denaturation	95°C	10-30 sec																												
Annealing	50-65°C	10-30 sec	30 cycles																											
Extension	72°C	1 kb/min																												
Final extension	72°C	5 min	1 cycle																											
4	 Analyze with gel electrophoresis	<p>6. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p> <p>7. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																												