

[Cat. No.] Please refer to the **Ordering Information**

Introduction

AccuPower® RocketPlex™ RT-PCR PreMix utilizes one-step RT-PCR performing cDNA synthesis and PCR of two or more target genes (up to 10 targets) in a single step, in a single tube. It can effectively synthesize cDNA with complex secondary RNA structures by using RocketScript™ Reverse Transcriptase (RocketScript™ RTase) developed by BIONEER. RocketScript™ RTase originated and engineered from M-MLV reverse transcriptase can provide increased thermal stability resulting in effective full length first strand cDNA synthesis. It can be used at various temperatures (42-70°C) allowing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptase. This product contains vacuum-dried components including RocketScript™ RTase, reaction buffer, DTT, dNTPs, RNase inhibitor, and HotStart Taq DNA Polymerase. It simplifies preparation of reverse transcription reaction mixture and PCR mixture by adding template RNA and primers without any extra process. After the reaction, samples can be applied directly on agarose gel for analysis due to included tracking dye.




Figure 1. Schematic diagram of 5' UTR of a gene, with complex secondary structure, at three different temperatures.

RocketScript™ Reverse Transcriptase is fully active at 70°C, so it can synthesize complete gene sequences that M-MLV and other reverse transcriptase can not synthesize.

Applications

- Multiplex RT-PCR
- Low copy detection
- Gene-expression analysis

Features & Benefits

- Multiplex RT-PCR: Simultaneous cDNA synthesis and amplification up to 10 target genes.
- Specificity: Minimized non-specific amplification and maximized PCR efficiency by using HotStart Taq DNA Polymerase.
- Ease-of-use: Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform cDNA synthesis and PCR by adding template RNA and primers. In addition, tracking dye and sedimentation agents are included in products for electrophoresis.

Composition

Composition	Concentration
RocketScript™ Reverse Transcriptase	200 U
5X Reaction buffer	1X
DTT	0.25 mM
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM

RNase inhibitor	1 U
HotStart Taq DNA Polymerase	1 U
Stabilizer and tracking dye	1X

Specifications

HotStart Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 1 kb

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols


Ordering Information

Description	Cat. No.
0.2 ml thin-wall 96 tubes	20 µl/rxn K-2211
8-tube strips	50 µl/rxn K-2213
with attached cap	20 µl/rxn K-2212
480 tubes	50 µl/rxn K-2214




Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols



Experimental Procedures

Steps		Procedure Details																												
1	 Preparation of reaction mixture	<p>1. Add template RNA, primers and nuclease-free water into <i>AccuPower® RocketPlex™</i> RT-PCR PreMix tubes to make a total volume of 20 µl or 50 µl. Do not include the dried pellet.</p> <ul style="list-style-type: none"> Amount of template RNA and primers <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Primers</td> <td>Gene specific primer</td> <td>5-10 pmol</td> <td>5-10 pmol</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried pellet by pipetting or vortexing, and briefly spin down.</p>	Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg	Primers	Gene specific primer	5-10 pmol	5-10 pmol													
Components		20 µl reaction	50 µl reaction																											
Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg																											
	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg																											
Primers	Gene specific primer	5-10 pmol	5-10 pmol																											
2	 RT-PCR	<p>3. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>42-70°C</td> <td>10-60 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>10-30 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>50-65°C</td> <td>10-30 sec</td> <td>30 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>1 kb/min</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: Recommended temperature of reverse transcription reaction is 50°C.</p>	Step	Temperature	Time	Cycles	cDNA synthesis	42-70°C	10-60 min	1 cycle	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	10-30 sec		Annealing	50-65°C	10-30 sec	30 cycles	Extension	72°C	1 kb/min		Final extension	72°C	5 min	1 cycle
Step	Temperature	Time	Cycles																											
cDNA synthesis	42-70°C	10-60 min	1 cycle																											
Pre-denaturation	95°C	5 min	1 cycle																											
Denaturation	95°C	10-30 sec																												
Annealing	50-65°C	10-30 sec	30 cycles																											
Extension	72°C	1 kb/min																												
Final extension	72°C	5 min	1 cycle																											
3	 Analyze with gel electrophoresis	<p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p> <p>5. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																												