

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

AccuPower® Ligation PreMix is a New, Powerful, Ready-to-use Ligation reagent optimized for more accurate Ligation reactions.

Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme combines blunt-end and cohesive-end termini as well as repairs single-stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids.
Source	T4 DNA Ligase was isolated from a recombinant <i>E. coli</i> strain.
Storage Condition	4°C to -20°C (-20°C is recommended for longer storage.)
Applications	<ul style="list-style-type: none"> • Gene cloning: Joining double-stranded DNAs with cohesive or blunt ends. • Cloning into vectors • Library construction • TA cloning • Linker ligation • Recirculization of linear DNA
Advantages	
Speed	Only 5 minutes for ligation of cohesive end DNA and 10 minutes for ligation of blunt end DNA at room temperature
Stability	Stabilized for 4 months at room temperature and 3 years at -20°C
Simplicity	Ready to use, AccuPower® Ligation PreMix is a single-tube lyophilized premix including the necessary components for ligation reaction such as T4 DNA Ligase, ATP, reaction buffer and stabilizer. You can finish 96 samples ligation within 30 minutes by only adding your vector and insert DNA.
Reproducibility and Yield	Strict functional QC assays provide highly consistent and reproducible ligation performance. In most applications, an increase in yield is observed as compared to standard reactions.
Unit Definition	2 Weiss units / Reaction. One Weiss ⁽³⁾ unit is defined as the amount of enzyme required to catalyze the conversion of 1 nmole of ³² P pyrophosphate into Norit-adsorbable material in 20 minutes at 37°C. 0.01 Weiss unit of T4 DNA Ligase is the amount of enzyme required to catalyze ligation reaction of greater than 95 % of 1µg of Lambda DNA / Hind III fragment in 20 minutes at 16°C.
Quality Controls	<ul style="list-style-type: none"> • Cohesive-end Ligation: > 95% of 1µg Lambda / Hind III fragments are ligated within 5 minutes in 20µl reaction condition. • Blunt-end Ligation: > 95% of 1µg Lambda / EcoR V fragments are ligated within 10 minutes in 20µl reaction condition.
Experimental Protocol	<ol style="list-style-type: none"> 1. Prepare mixture of insert DNA and vector using the molar ratio specified by your standard ligation protocol. (We recommend the molar ratio of vector : insert DNA = 1 : 3) Molar ratio: $\frac{\square \text{ ng of vector} \times \square \text{ kb size of insert DNA}}{\square \text{ kb size of vector}} \times \frac{1}{3}^a = \square \text{ ng of insert DNA}$ ^a : molar ratio of insert DNA / vector 2. Add 20µl of the mixture to AccuPower® Ligation PreMix tube. If the volume of mixture is less than 20 µl, add distilled to water to adjust the final volume to 20µl. 3. Mix completely by pipetting and vortexing vigorously and then spin down briefly. 4. Incubate the reaction mixture at room temperature for 5~10 minutes. 5. Collect 10µl from tube and perform transformation with 100µl of competent cells. In case an electroporation of ligated products, salt precipitate the DNA before electroporation in order to remove the salts contained in ligation mixture.
Heat Inactivation	70°C for 10 minutes.
References	<ol style="list-style-type: none"> 1. Engler, M. J. and Richardson, C. C. (1982) In: <i>The Enzymes</i>, Boyer, P. D. ed., Academic Press, New York, NY. 2. Zimmerman, S. B. and Pfeiffer, B. H. (1983) <i>Proc. Natl. Acad. Sci. USA</i> 80, 5852-5856 3. Weiss, B., et al., (1968) <i>J. Biol. Chem.</i>, 243, 4543-4555

Ordering Information

Tube type	Reaction	Cat.No.	Description
0.2ml Tube	20ul	K-7103	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes

Related Products

Cat. No.	Description
E-3061	T4 DNA Ligase (20,000 unit)
E-3062	T4 DNA Ligase (100,000 unit)
E-3111	Thermostabe <i>Thermus filiformis</i> (<i>Tfi</i>) DNA Ligase (2,000 unit)
E-3112	Thermostabe <i>Thermus filiformis</i> (<i>Tfi</i>) DNA Ligase (10,000 unit)
K-7101	AccuRapid™ DNA Ligation Kit (50 ligations)
K-7102	AccuRapid™ DNA Ligation Kit (150 ligations)

A complete product list appears on our web site at www.bioneer.com.