

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

AccuPower® RT Master Mix is a new, powerful, ready-to-use RT kit for the synthesis of cDNA with superiority to other RT products.

Advantages

Speed

Substantial reduction in reaction setup time.

Stability

As each tube of *AccuPower®* RT Master Mix contains a stabilizer (patented in US and Korea), which can maintain the stability of the RTase up to 2 years at -20°C.

Reproducibility and Yield

The strict functional QC assays demonstrated highly consistent and reproducible RT performance. In most applications an increase in yield is observed as compared to the standard reactions.

Experimental Protocol

- Mix the template RNA and primer in a sterile tube as indicated below;

Concentration of template RNA and primer			
Reaction volume		20 µl reaction	50 µl reaction
Template RNA	Total RNA	0.5~1.0 µg	1.0~2.0 µg
	Poly(A) RNA	0.05~0.1 µg	0.1~0.2 µg
Primer	Oligo dT ₁₈	0.5 µg(100 pmole)	1.0 µg(200 pmole)
	Sequence specific	10~30 pmole	20~50 pmole

- Incubate the mixture at 70°C for 5 min and place it on ice.
- Transfer the incubated mixture and *AccuPower®* RT Master Mix(2X Master mix) 10 ul to an 0.2 ml Tube, and then fill up the reaction volume with DEPC – DW.
- Mix the reaction mixture by vortexing or pipetting, and briefly spin down.
- Add mineral oil to each tube (This step is unnecessary when using a thermal cycler with top heating).
- Perform cDNA synthesis reaction as follows :

42°C, 60 min. (cDNA synthesis)

94°C, 5 min. (RTase inactivation)

- If PCR is needed following RT reaction, perform the PCR using *AccuPower®* PCR PreMix as follows :

- i) Transfer 2~5 µl of the RT product (synthesized cDNA) to *AccuPower®* PCR PreMix tube.
- ii) Perform PCR cycles according to the PCR condition.

(Annealing temperature and time need to be optimized for each primer/template combination.)

Ordering Information

Tube type	Reaction	Cat.No	Description	Tube type	Reaction	Cat.No	Description
0.2 ml Tube	20 µl	K-2041	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes	96-well	20 µl	K-2261-4	thin-wall 96-well flat plate
		K-2041-B	0.2 ml thin-wall 8-strip tubes with attached cap/480 tubes			K-2261-5	thin-wall 96-well full-skirted plate
	50 µl	K-2043	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes			K-2261-6	thin-wall 96-well semi-skirted plate
		K-2043-B	0.2 ml thin-wall 8-strip tubes with attached cap/480 tubes	384-well	5 µl	K-2082-1	thin-wall 384-well full-skirted plate
96-well	10 µl	K-2261-1	thin-wall 96-well flat plate		10 µl	K-2082-2	thin-wall 384-well full-skirted plate
		K-2261-2	thin-wall 96-well full-skirted plate		20 µl	K-2082-3	thin-wall 384-well full-skirted plate
		K-2261-3	thin-wall 96-well semi-skirted plate		0.5ml Tube	50 µl	K-2042
1.5ml Tube	20 µl	K-2263	2X Master Mix				

Notice

Bioneer corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice. All information provided here is subject to change without notice.

A complete product list appears on our web site at www.bioneer.com