

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

I. Introduction

AccuPower® RocketScript™ RT Master Mix contains all components for first-strand cDNA synthesis from purified Poly(A) or total RNA template.

The AccuPower® RocketScript™ RT Master Mix contains RocketScript Reverse Transcriptase, a new M-MLV originated Reverse Transcriptase that has been engineered to provide increased thermal stability in order to synthesize full length first-strand cDNA more efficiently. The amount of starting material can vary from 1 pg to > 1 µg of total RNA and RNA targets from 100 bp to > 10 kb can be detected with the AccuPower® RocketScript™ RT Master Mix.

The AccuPower® RocketScript™ RT Master Mix can be used to synthesize cDNA at a temperature range of 42 – 70°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases.

II. Application

- Standard RT and RT-PCR
- Real-Time PCR
- Synthesis of double-stranded cDNA for cloning
- Gene expression level analysis

III. Contents

| Component | Amount |
|------------------------------------|-------------|
| RocketScript Reverse Transcriptase | 200 U |
| 5 x Reaction Buffer | 1 x |
| DTT | 0.25 mM |
| dNTP | 250 µM each |
| RNase Inhibitor | 1 U |

IV. Principle

RocketScript Reverse Transcriptase in the AccuPower® RocketScript™ RT Master Mix is genetically engineered thermal stable M-MLV Reverse Transcriptase with enhanced thermal stability and outstanding processivity. The enzyme also features increased specificity and improved efficiency allowing efficient reverse transcription of RNA molecules with complex secondary structures.

RocketScript Reverse Transcriptase at 70°C




Figure 1. Schematic representation of the 5'UTR of a gene, with complex secondary structure, at three different temperatures.

Note that RocketScript shows full activity at 70°C allowing it to synthesize the complete gene sequence where M-MLV and other Reverse Transcriptase's fail.

V. Storage

AccuPower® RocketScript™ RT Master Mix should be stored at -20°C upon receipt and is stable until the expiry date stated on the label.

VI. Notice to Purchaser

- AccuPower® RocketScript™ RT Master Mix exhibits RNase H activity.
- Bioneer corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice. All information provided here is subject to change without notice.

VII. Additional Required Materials & Devices

- Thermal cycler for PCR
- Target-specific primers, oligo dT, random hexamer, or nonamers
- Calibrated micropipette
- Sterilized micropipette tips with filters

VIII. General Precautions

- Wear gloves during experiments to prevent contamination.
- Store positive materials, such as samples and control templates, in separated freezer from freezers for the kit.
- Add templates to the reaction mixture in clean bench or a spatially separated facility.

IX. Protocol

1. Thaw total RNA, DEPC-water, and primer before use.
2. Mix the total RNA and primer (oligo dT, random primer, or specific primer) in a sterile tube as indicated below.

Recommended amount of template and primer

| Components | Amount |
|-------------------------|----------------|
| Total RNA | 10 pg – 5 µg |
| Oligo dT | 50 pmoles |
| Random primer(dN6,dN12) | 100 pmoles |
| Specific primer | 10 – 50 pmoles |

3. Fill up the reaction volume with DEPC-water.
4. Mix the reaction mixture by vortexing or pipetting and briefly spinning down.
5. Perform the reaction under the following conditions.

| Step | Temperature | | | | Time |
|-------------------|-----------------|------------------|------------------|--------------------------|--------|
| | dN ₆ | dN ₁₂ | dT ₂₀ | Specific primer | |
| Primer annealing | 15°C | 30°C | 37°C | T _m of primer | 10 min |
| cDNA synthesis | 42°C | | | | 30 min |
| Heat inactivation | 95°C | | | | 5 min |

Alternative protocol

| Step | Temperature | Time |
|-------------------|-------------|-------|
| cDNA synthesis | 42~70°C | 1 hr |
| Heat inactivation | 95°C | 5 min |

Note: for difficult or high GC-content templates, use a 55°C cDNA synthesis temperature.

6. Maintain the reaction at 4°C after amplification. The sample can be stored at -20°C until use.

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

X. Reaction Examples

1. Reaction mixture

| Component | Volume | Amount |
|--------------|--------|-----------|
| Template RNA | 1 µl | 100 pg |
| Oligo dT 20 | 1 ul | 50 pmoles |
| RT Mastermix | 10 ul | |
| DEPC-D.W | 8 µl | |
| Total | 20 µl | |

2. Reaction condition

| Step | Temperature | Time |
|-------------------------------|-------------|--------|
| Primer annealing (oligo dT20) | 37°C | 10 min |
| cDNA synthesis | 60°C | 1 hr |
| Heat inactivation | 95°C | 5 min |

XI. Experimental Data




Figure 1. Amplification results of AccuPower® RocketScript™ RT Master Mix using myc compared with Company I reverse transcriptase.

Reverse transcription condition: Incubation at each temperature 45, 50, 55, 60, 65, 70°C for 1 hr, inactivation at 95°C for 5 min / Primer set: human myc 495 bp set

Lane M: 1 kb DNA Ladder
Lane 1: 100 ng Human total RNA from HeLa cell
Lane 2: 10 ng Human total RNA from HeLa cell
Lane 3: 1 ng Human total RNA from HeLa cell
Lane 4: 100 pg Human total RNA from HeLa cell




Figure 2. Comparison of long kb amplification between AccuPower® RocketScript™ RT Master Mix and Company I reverse transcriptase.

Reverse transcription reactions were performed according to each manufacturer's recommendation. All cDNAs were amplified with AccuPower® HotStart PCR Master Mix (K-5050) from Bioneer. Note supplier I shows inhibition with high input concentration of total RNA.

Lane 1: 1 µg Human total RNA from HeLa cell
Lane 2: 100 ng Human total RNA from HeLa cell
Lane 3: 10 ng Human total RNA from HeLa cell




Figure 3. Performance on Comparison oligo (Oligo dT 20, dN 6, dN 12)

Template RNA: Human total RNA in Hela cell
Concentration 0.01 ng ~ 100 ng, 10 fold serial dilution.

XII. Notice

Bioneer corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice. All information provided here is subject to change without notice.

XIII. Ordering Information

| Cat. No. | Description |
|----------|---|
| K-2105 | AccuPower® RocketScript™ RT Master Mix (2X, 1 ml) |

A complete product list appears on our web site at www.bioneer.com