

Full Manual(EN)



For detailed information,
visit the website

Full Manual(KR)



Innovation • Value • Discovery

BIONEER

EcoQprep™ Genomic DNA Kit

• 50 rxn

K-3701

Name	Volume	Quantity	Name	Volume	Quantity
TL Buffer	15 mL	1 ea	Proteinase K powder, lyophilized	25 mg	1 ea
GB Buffer	15 mL	1 ea	RNaseA powder, lyophilized	24 mg	1 ea
WM1 Buffer	30 mL	1 ea	Poly(A), lyophilized	1 mg	1 ea
WB2 Buffer	10 mL	1 ea	Magnetic Nanobead	6 mL	1 ea
WE Buffer	40 mL	1 ea	1.5 mL Tube	50 ea	1 pack
EA Buffer	25 mL	1 ea	One Page Protocol	-	1 ea

EcoQprep™ Genomic DNA Kit

One Page Protocol

Before Start

Components	Before Use	Storage	Symbol
TL Buffer	If precipitation occurs, incubate on a heating block at 60°C before use.	Room temperature	Run centrifuge
WM1 and WB2 Buffer	As indicated on the bottle label, add the correct amount of absolute ethanol (not provided).	Room temperature	
Proteinase K	Dissolve in 1,250 µL of nuclease-free water.	Short-term : 4°C / Long-term : -20°C	Heating
RNase A	Dissolve in 600 µL of nuclease-free water.	Short-term : 4°C / Long-term : -20°C	Incubating (RT ~ 37°C)
Poly (A)	Dissolve in 1,000 µL of nuclease-free water.	Short-term : 4°C / Long-term : -20°C	Caution
Lysozyme (Not provided)	Prepare at a concentration of 100 mg/ml for gram-positive bacterial lysis.	-	Add the following solution
Lysis Buffer (Not provided)	For gram-positive bacteria, prepare as follows: 20 mM Tris-HCl (pH8.0), 2 mM Sodium-EDTA, 1.2% Triton X-100.	-	

Protocol

NOTE It is recommended to vortex after adding the solution.

Sample		Pre-treatment ①	Lysis ②	Precipitation ③	Binding ④	1 st Wash ⑤	2 nd Wash ⑥	3 rd Wash ⑦	Elution ⑧
Whole blood Buffy coat		Sample : 200 µL	⊕ 20 µL Proteinase K ⊕ 200 µL GB 🔥 60°C, 10 min	⊕ 400 µL EtOH	⊕ 100 µL Bead 1) Attach magnet	⊕ 500 µL WM1 Follow 1) ~ 4) as below. 🔴 Repeat whole steps again	⊕ 700 µL WB2 1) Attach magnet	⊕ 700 µL WE 🔴 Add to the Opposite side of the beads	⊕ 100 µL EA 🔥 60°C, 1 min
Cultured cells		Cell pellet : ~ 1 × 10 ⁶ cells 📺 300 xg, 10 min ⊕ 200 µL PBS(1X)	⊕ 20 µL Proteinase K ⊕ 10 µL RNase A 🌿 RT, 2 min ⊕ 200 µL GB 🔥 60°C, 10 min	⊕ 400 µL EtOH	 2) Invert 3~4X	⊕ 700 µL WM1 1) Attach magnet 2) Invert 3~4X	 2) Invert 3~4X	1) Gently invert 2X 	1) Attach magnet
Animal tissues		Homogenize sample : ~25 mg ⊕ 180 µL TL	⊕ 20 µL Proteinase K ⊕ 10 µL RNase A 🌿 RT, 2 min 🔥 60°C, until lysed	⊕ 200 µL GB ⊕ 400 µL EtOH	 3) Discard sup	2) Invert 3~4X 3) Discard sup	 3) Discard sup	2) Discard sup 	2) Transfer the eluent
Bacteria	Gram (-)	Cell pellet : ~ 1 × 10 ⁹ cells 📺 6000 xg, 10 min ⊕ 180 µL TL	⊕ 20 µL Proteinase K ⊕ 10 µL RNase A 🌿 RT, 2 min 🔥 60°C, until lysed	⊕ 200 µL GB ⊕ 400 µL EtOH	 4) Detach magnet	3) Discard sup 4) Detach magnet	 3) Discard sup	3) Detach magnet 	
	Gram (+)	Cell pellet : ~ 1 × 10 ⁹ cells 📺 6000 xg, 10 min ⊕ 180 µL Lysis Buffer	⊕ 20 µL Lysozyme ⊕ 10 µL RNase A 🌿 37°C, 30 min ⊕ 20 µL Proteinase K ⊕ 200 µL GB 🔥 60°C, until lysed	⊕ 400 µL EtOH		4) Detach magnet 	🔴 Do Not Detach magnet		