

Innovation • Value • Discovery

BIONEER

- 50 rxn

Name	Volume	Quantity	Name	Volume	Quantity
------	--------	----------	------	--------	----------

RD Buffer	30 mL	1 ea	1.5 mL Tube	50 ea	1 pack
RWMB1 Buffer	25 mL	1 ea	One Page Protocol	-	1 ea
RWB2 Buffer	15 mL	1 ea			
WE Buffer	40 mL	1 ea			
ER Buffer	25 mL	1 ea			
Magnetic Nanobead	6 mL	1 ea			

EcoQprep™ Universal RNA Kit

One Page Protocol

Before Start

Components	Before Use	Storage	Symbol
RD Buffer	Before each experiment, mix RD buffer with 10 µL of β-mercaptoethanol per 1 mL of RD buffer.	-	Run centrifuge Heating Add the following solution Caution
RWMB1 and RWB2 Buffer	As indicated on the bottle label, add the correct amount of absolute ethanol (not provided).	Room temperature	
80% EtOH (Not provided)	80% EtOH is needed only for the drying bead method.	Room temperature	

Protocol

NOTE It is recommended to vortex after adding the solution.

Sample		Pre-treatment ①	Lysis ②	Precipitation ③	Binding ④	1 st Wash ⑤	2 nd Wash ⑥	3 rd Wash ⑦	Elution ⑧
Cultured Cells	Suspension cell culture	Cell pellet ~1×10 ⁶ cells ☑ 300 xg, 5 min	⊕ 500 µL RD	⊕ 300 µL EtOH	⊕ 100 µL Bead 1) Attach magnet  2) Invert 3~4X 	⊕ 700 µL RWMB1 1) Attach magnet  2) Invert 3~4X 	⊕ 700 µL RWB2 1) Attach magnet  2) Invert 3~4X 	Perform either A or B A) Washing Bead ⊕ 700 µL WE ❶ Add to the Opposite side of the beads 1) Gently invert 2X 2) Discard sup 3) Detach magnet For 1) ~ 3), refer to 1 st and 2 nd Wash B) Drying Bead ⊕ 700 µL 80% EtOH 1) Detach magnet 2) Attach magnet 3) Invert 3~4X 4) Discard sup 5) Detach magnet 6) ⬇ 60°C, 5 min For 2) ~ 5), refer to 1 st and 2 nd Wash	⊕ 50 - 100 µL ER ⬇ 60°C, 1 min 1) Attach magnet  2) Transfer the eluent 
	Monolayer cell culture	Direct harvest 1) Discard cell culture medium Trypsin harvest 1) Detach cells 2) ☑ 300 xg, 5 min							
	Plant tissue	Tissue : ~100 mg 1) Grind in LN ₂	⊕ 500 µL RD ⬇ 60°C, 1~3 min ☑ full speed, 2 min Transfer the supernatant		3) Discard sup 	3) Discard sup 	3) Discard sup 		
Animal tissue		Homogenize sample : 20 ~ 30 mg	⊕ 500 µL RD ☑ full speed, 3 min Transfer the supernatant		4) Detach magnet 	4) Detach magnet 	❶ Do Not Detach magnet		