

[Cat. No.] **E-3131, E-3132**

## Introduction

*CycleScript™* Reverse Transcriptase is applied with BIONEER's patent technology called Cyclic Reverse Transcription (CRT), which not only increases the efficiency, but also is effective for full-length cDNA synthesis. The CRT reaction can perform homogeneous primer-annealing at low temperature as well as reverse transcription with high sensitivity at high temperature, which even complex RNA secondary structure is resolved. It is composed of 2 or 3 steps. The step 1 is performed at 15-25°C, at which short primers are fully annealed. Then, the step 2 is conducted at 42-48°C for cDNA synthesis and the step 3 is followed at high temperature 50-55°C at which secondary structure of template RNA obstructing reverse transcription is released.

## Applications

- First-strand cDNA synthesis from RNA molecules
- RT-PCR
- Random priming reaction
- Library construction
- Probe labeling
- mRNA 5' end mapping by primer extension analysis

## Components

Components	E-3131	E-3132
<i>CycleScript™</i> Reverse Transcriptase	10,000 U (50 µl)	50,000 U (50 µl x 5)
5X Reaction buffer	0.4 ml	0.4 ml x 5
10 mM dNTPs	0.2 ml	0.2 ml x 5
100 mM DTT	0.2 ml	0.2 ml x 5

\* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

## Specifications

<i>CycleScript™</i> Reverse Transcriptase	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	No
Strand displacement	Yes
Fragment size	Up to 9 kb

## Buffer Composition

5X Reaction buffer	150 mM Tris, 250 mM KCl, 10 mM MgCl <sub>2</sub> , etc, pH 8.1
--------------------	--

## Storage Buffer

*CycleScript™* Reverse Transcriptase is supplied in 50% (v/v) glycerol containing 20 mM Tris-HCl, 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, and stabilizer, pH 7.6.

## Unit Definition

One unit is defined as the amount of enzyme that incorporates 1 nmole of dTTP into acid-insoluble products in 10 min at 37°C using poly(A)-oligo(dT) as template primer.

## Quality Control

- Nuclease Contamination Assay: Nuclease activity is not detected after incubation of 1 µg of DNA and RNA with 200 U of *CycleScript™* Reverse Transcriptase at 37-42°C for 3 hrs.

## Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

## Online Resources



Korean



English

Visit our **product page** for additional information and protocols.

## Ordering Information

Description	Cat. No
<i>CycleScript™</i> Reverse Transcriptase	10,000 U (50 rxn) E-3131
	50,000 U (250 rxn) E-3132

## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols



Batch Code



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Research Use Only





Temperature Limitation



Use-by Date

**Experimental Procedures**

Steps		Procedure Details																																																	
1	 <p><b>Preparation of reaction mixture</b></p>	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> <li>Amount of template RNA and primers</li> </ul>	<table border="1"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>1 µg</td> </tr> <tr> <td>RNA</td> <td>5-100 ng</td> </tr> <tr> <td rowspan="2">Primers</td> <td>Oligo dT or Random primer</td> <td>10-100 pmol</td> </tr> <tr> <td>Sequence specific</td> <td>10-30 pmol</td> </tr> </tbody> </table>	Components		20 µl reaction	Template RNA	Total RNA	1 µg	RNA	5-100 ng	Primers	Oligo dT or Random primer	10-100 pmol	Sequence specific	10-30 pmol																																			
		Components		20 µl reaction																																															
Template RNA	Total RNA	1 µg																																																	
	RNA	5-100 ng																																																	
Primers	Oligo dT or Random primer	10-100 pmol																																																	
	Sequence specific	10-30 pmol																																																	
2	 <p><b>cDNA synthesis</b></p>	<p>2. Add all components for cDNA synthesis into PCR tubes (not provided) to a total volume of 20 µl.</p> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul>	<table border="1"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td colspan="2">Step 1 mixture</td> <td>Variable</td> </tr> <tr> <td colspan="2">5X Reaction buffer</td> <td>4 µl</td> </tr> <tr> <td colspan="2">10 mM dNTPs</td> <td>2 µl or Variable</td> </tr> <tr> <td colspan="2">100 mM DTT</td> <td>2 µl</td> </tr> <tr> <td colspan="2">RNase inhibitor (Not provided)</td> <td>Variable</td> </tr> <tr> <td colspan="2">CycleScript™ Reverse Transcriptase (200 U/ µl)</td> <td>200 U</td> </tr> <tr> <td colspan="2">Nuclease-free water</td> <td>Variable</td> </tr> <tr> <td colspan="2">Total volume</td> <td>20 µl</td> </tr> </tbody> </table>	Components		20 µl reaction	Step 1 mixture		Variable	5X Reaction buffer		4 µl	10 mM dNTPs		2 µl or Variable	100 mM DTT		2 µl	RNase inhibitor (Not provided)		Variable	CycleScript™ Reverse Transcriptase (200 U/ µl)		200 U	Nuclease-free water		Variable	Total volume		20 µl																					
		Components		20 µl reaction																																															
Step 1 mixture		Variable																																																	
5X Reaction buffer		4 µl																																																	
10 mM dNTPs		2 µl or Variable																																																	
100 mM DTT		2 µl																																																	
RNase inhibitor (Not provided)		Variable																																																	
CycleScript™ Reverse Transcriptase (200 U/ µl)		200 U																																																	
Nuclease-free water		Variable																																																	
Total volume		20 µl																																																	
		<p>3. Perform the reaction under the following conditions.</p> <p>3-1. Cyclic reverse transcription (example 1)</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td>15-25°C</td> <td>30 sec</td> <td></td> </tr> <tr> <td>cDNA synthesis</td> <td>42-45°C</td> <td>4 min</td> <td>12 cycles or less</td> </tr> <tr> <td>Melting secondary structure &amp; cDNA synthesis</td> <td>55°C</td> <td>30 sec</td> <td></td> </tr> <tr> <td>Heat inactivation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>3-2. Cyclic reverse transcription (example 2)</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td>15-25°C</td> <td>1 min</td> <td></td> </tr> <tr> <td>Melting secondary structure &amp; cDNA synthesis</td> <td>42-50°C</td> <td>4 min</td> <td>12 cycles or less</td> </tr> <tr> <td>Heat inactivation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: For difficult or high GC-content templates, perform cDNA synthesis at 55°C.</p> <p>3-3. Single temperature reaction (example 3)</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>37-50°C</td> <td>30-60 min</td> <td>1 cycle</td> </tr> <tr> <td>Heat inactivation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: Recommended temperature is range of 42-48°C.</p>	Step	Temperature	Time	Cycles	Primer annealing	15-25°C	30 sec		cDNA synthesis	42-45°C	4 min	12 cycles or less	Melting secondary structure & cDNA synthesis	55°C	30 sec		Heat inactivation	95°C	5 min	1 cycle	Step	Temperature	Time	Cycles	Primer annealing	15-25°C	1 min		Melting secondary structure & cDNA synthesis	42-50°C	4 min	12 cycles or less	Heat inactivation	95°C	5 min	1 cycle	Step	Temperature	Time	Cycles	cDNA synthesis	37-50°C	30-60 min	1 cycle	Heat inactivation	95°C	5 min	1 cycle	
Step	Temperature	Time	Cycles																																																
Primer annealing	15-25°C	30 sec																																																	
cDNA synthesis	42-45°C	4 min	12 cycles or less																																																
Melting secondary structure & cDNA synthesis	55°C	30 sec																																																	
Heat inactivation	95°C	5 min	1 cycle																																																
Step	Temperature	Time	Cycles																																																
Primer annealing	15-25°C	1 min																																																	
Melting secondary structure & cDNA synthesis	42-50°C	4 min	12 cycles or less																																																
Heat inactivation	95°C	5 min	1 cycle																																																
Step	Temperature	Time	Cycles																																																
cDNA synthesis	37-50°C	30-60 min	1 cycle																																																
Heat inactivation	95°C	5 min	1 cycle																																																
		<p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>																																																	