

[Cat. No.] **E-3161, E-3162**

## Introduction

RocketScript™ Reverse Transcriptase, RNase H Minus is an RNA and DNA-dependent polymerase removed with RNase H activity by introducing mutation in RNase H domain. This engineered enzyme provides enhanced extensibility, resulting in higher yields of cDNA and longer product (up to 13 kb) than RocketScript™ Reverse Transcriptase. It can be used to synthesize cDNA at a temperature range of 42-70°C, providing increased specificity and it can successfully synthesize cDNA even with 1 pg of human total RNA due to its sensitivity.

## Applications

- Gene synthesis
- First-strand cDNA synthesis from RNA molecules
- RT-PCR
- Random priming reactions
- Library construction
- Probe labeling
- mRNA 5'end mapping by primer extension analysis
- Real-time PCR

## Components

Components	E-3161	E-3162
RocketScript™ Reverse Transcriptase, RNase H Minus	10,000 U (50 µl)	50,000 U (50 µl x 5)
5X Reaction buffer	0.5 ml	0.5 ml x 5
10 mM dNTPs	0.2 ml	0.2 ml x 5
100 mM DTT	0.2 ml	0.2 ml x 5
RNase inhibitor	50 µl	50 µl x 5

\* Note: For research use only. Not for use in diagnostic or therapeutic procedures.

## Specifications

RocketScript™ Reverse Transcriptase, RNase H Minus	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	No
Fragment size	Up to 12.5 kb

## Buffer Composition

5X Reaction buffer	250 mM Tris-HCl, 375 mM KCl, 15 mM MgCl <sub>2</sub> , and stabilizer pH 8.3
--------------------	--

## Storage Buffer

RocketScript™ Reverse Transcriptase, RNase H Minus is supplied in 50% (v/v) glycerol containing 20 mM Tris-HCl, 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, and stabilizer, pH 7.6.

## Unit Definition

One unit is defined as the amount of enzyme that incorporates 1 nmole of dTTP into acid-insoluble products in 10 min at 37°C using poly(A)-oligo(dT) as template primer.

## Quality Control

- Nuclease Contamination Assay: Nuclease activity is not detected after incubation of 1 µg of DNA and RNA with 200 U of RocketScript™ Reverse Transcriptase, RNase H minus at 37-42°C for 3 hrs.

## Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

## Online Resources



Korean



English

Visit our [product page](#) for additional information and protocols.

## Ordering Information

Description	Cat. No
RocketScript™ Reverse Transcriptase, RNase H Minus	10,000 U (50 rxn) E-3161
	50,000 U (250 rxn) E-3162

## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols



Batch Code



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Research Use Only





Temperature Limitation



Use-by Date

**Experimental Procedures**

Steps		Procedure Details																																																																					
1	 <b>Preparation of reaction mixture</b>	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> <li>Amount of template RNA and primers</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2" style="text-align: center;">Template RNA</td> <td style="text-align: center;">Total RNA</td> <td style="text-align: center;">1-5 µg</td> </tr> <tr> <td style="text-align: center;">mRNA</td> <td style="text-align: center;">5-100 ng</td> </tr> <tr> <td rowspan="2" style="text-align: center;">Primers</td> <td style="text-align: center;">Oligo dT or Random primer</td> <td style="text-align: center;">10-100 pmol</td> </tr> <tr> <td style="text-align: center;">Sequence specific</td> <td style="text-align: center;">10-30 pmol</td> </tr> </tbody> </table> <p>2. Add all components for cDNA synthesis into PCR tubes (not provided) to a total volume of 20 µl.</p> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td colspan="2">Mixture of template RNA and primers</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td colspan="2">5X Reaction buffer</td> <td style="text-align: center;">4 µl</td> </tr> <tr> <td colspan="2">10 mM dNTPs</td> <td style="text-align: center;">2 µl or Variable</td> </tr> <tr> <td colspan="2">100 mM DTT</td> <td style="text-align: center;">2 µl</td> </tr> <tr> <td colspan="2">RNase inhibitor (100 ng/µl)</td> <td style="text-align: center;">0.5-1 µl</td> </tr> <tr> <td colspan="2"><i>RocketScript™</i> Reverse Transcriptase, RNase H Minus (200 U/µl)</td> <td style="text-align: center;">200 U</td> </tr> <tr> <td colspan="2">Nuclease-free water</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td colspan="2">Total volume</td> <td style="text-align: center;">20 µl</td> </tr> </tbody> </table>	Components		20 µl reaction	Template RNA	Total RNA	1-5 µg	mRNA	5-100 ng	Primers	Oligo dT or Random primer	10-100 pmol	Sequence specific	10-30 pmol	Components		20 µl reaction	Mixture of template RNA and primers		Variable	5X Reaction buffer		4 µl	10 mM dNTPs		2 µl or Variable	100 mM DTT		2 µl	RNase inhibitor (100 ng/µl)		0.5-1 µl	<i>RocketScript™</i> Reverse Transcriptase, RNase H Minus (200 U/µl)		200 U	Nuclease-free water		Variable	Total volume		20 µl																													
Components		20 µl reaction																																																																					
Template RNA	Total RNA	1-5 µg																																																																					
	mRNA	5-100 ng																																																																					
Primers	Oligo dT or Random primer	10-100 pmol																																																																					
	Sequence specific	10-30 pmol																																																																					
Components		20 µl reaction																																																																					
Mixture of template RNA and primers		Variable																																																																					
5X Reaction buffer		4 µl																																																																					
10 mM dNTPs		2 µl or Variable																																																																					
100 mM DTT		2 µl																																																																					
RNase inhibitor (100 ng/µl)		0.5-1 µl																																																																					
<i>RocketScript™</i> Reverse Transcriptase, RNase H Minus (200 U/µl)		200 U																																																																					
Nuclease-free water		Variable																																																																					
Total volume		20 µl																																																																					
2	 <b>cDNA synthesis</b>	<p>3. Perform the reaction under the following conditions.</p> <p>3-1. Cyclic reverse transcription (example 1)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Step</th> <th colspan="4">Temperature</th> <th rowspan="2">Time</th> <th rowspan="2">Cycles</th> </tr> <tr> <th>dN<sub>6</sub></th> <th>dN<sub>12</sub></th> <th>dT<sub>20</sub></th> <th>Sequence specific</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td style="text-align: center;">15°C</td> <td style="text-align: center;">30°C</td> <td style="text-align: center;">37°C</td> <td style="text-align: center;">T<sub>m</sub> of primers</td> <td style="text-align: center;">10-30 sec</td> <td rowspan="3" style="text-align: center;">12 cycles or less</td> </tr> <tr> <td>cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">50°C</td> <td style="text-align: center;">4 min</td> </tr> <tr> <td>Melting secondary structure &amp; cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">55°C</td> <td style="text-align: center;">30 sec</td> </tr> <tr> <td>Heat inactivation</td> <td colspan="3"></td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>3-2. Single temperature reaction (example 2)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Step</th> <th colspan="4">Temperature</th> <th rowspan="2">Time</th> <th rowspan="2">Cycles</th> </tr> <tr> <th>dN<sub>6</sub></th> <th>dN<sub>12</sub></th> <th>dT<sub>20</sub></th> <th>Sequence specific</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td style="text-align: center;">15°C</td> <td style="text-align: center;">30°C</td> <td style="text-align: center;">37°C</td> <td style="text-align: center;">T<sub>m</sub> of primers</td> <td style="text-align: center;">1 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">42-70°C</td> <td style="text-align: center;">10-60 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Heat inactivation</td> <td colspan="3"></td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>	Step	Temperature				Time	Cycles	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific	Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	10-30 sec	12 cycles or less	cDNA synthesis				50°C	4 min	Melting secondary structure & cDNA synthesis				55°C	30 sec	Heat inactivation				95°C	5 min	1 cycle	Step	Temperature				Time	Cycles	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific	Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	1 min	1 cycle	cDNA synthesis				42-70°C	10-60 min	1 cycle	Heat inactivation				95°C	5 min	1 cycle
Step	Temperature				Time	Cycles																																																																	
	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific																																																																			
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	10-30 sec	12 cycles or less																																																																	
cDNA synthesis				50°C	4 min																																																																		
Melting secondary structure & cDNA synthesis				55°C	30 sec																																																																		
Heat inactivation				95°C	5 min	1 cycle																																																																	
Step	Temperature				Time	Cycles																																																																	
	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific																																																																			
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	1 min	1 cycle																																																																	
cDNA synthesis				42-70°C	10-60 min	1 cycle																																																																	
Heat inactivation				95°C	5 min	1 cycle																																																																	