

# *ExiPrep*<sup>™</sup> Plus Viral DNA/RNA Kit

Cat. No. K-4271  
K-4272  
K-4273



***ExiPrep*<sup>™</sup> Plus Viral DNA/RNA Kit (K-4271)**  
***ExiPrep*<sup>™</sup> Plus Viral DNA Kit (K-4272)**  
***ExiPrep*<sup>™</sup> Plus Viral RNA Kit (K-4273)**

Kit for the extraction of DNA/RNA from virus

## User Guide

**K-4271**  
 **96**

**K-4272**  
 **96**

**K-4273**  
 **96**

Version No.: 5 (2022-06-10)

Please read all the information in booklet before using the unit



**BIONEER Corporation**

**Bioneer Global Center, 71, Techno-2-ro,  
Yuseong-gu, Daejeon, 34013, Republic of Korea**

Tel: 1588-9788

Email: [sales@bioneer.co.kr](mailto:sales@bioneer.co.kr)

[www.bioneer.com](http://www.bioneer.com)

## **Intended Use**

*ExiPrep*<sup>™</sup> Plus Viral DNA/RNA Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

## **Safety Warning and Precaution**

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information, please consult the appropriate Material Safety Data Sheet (MSDS).

## **Warranty and Liability**

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to BIONEER within 30 days.

## **Quality Management System ISO 9001 Certified**

Every aspect of our quality management system from product development, production to quality assurance and supplier qualification meets the world-class standards.

## **Patent**

*ExiPrep*<sup>™</sup> and its kits are protected by the patents KR10-2015-0089172.

## **Trademark**

*ExiPrep*<sup>™</sup> is a trademark of BIONEER Corporation.

## **Copyright**

Copyright 2022 BIONEER Corporation. All Rights Reserved.

## **Notice**

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Contents

|   |           |
|---|-----------|
| <b>Product Information</b> .....                  | <b>1</b>  |
| Components .....                                  | 1         |
| Storage .....                                     | 1         |
| Specifications.....                               | 2         |
| Precautions .....                                 | 2         |
| <b>Introduction</b> .....                         | <b>3</b>  |
| Product Description .....                         | 3         |
| Principle .....                                   | 4         |
| Features & Benefits .....                         | 5         |
| Components of Buffer Cartridges .....             | 6         |
| <b>Experimental Procedures</b> .....              | <b>7</b>  |
| Loading the Kit to the Instrument .....           | 7         |
| <b>Troubleshooting</b> .....                      | <b>11</b> |
| <b>Appendix A: List of Protocol Numbers</b> ..... | <b>12</b> |
| <b>Ordering Information</b> .....                 | <b>15</b> |
| <b>Related Products</b> .....                     | <b>15</b> |
| <b>Explanation of Symbols</b> .....               | <b>16</b> |

## Product Information

### Components

| Components            | K-4271<br>(Viral DNA/RNA) | K-4272<br>(Viral DNA)   | K-4273<br>(Viral RNA)   |
|-----------------------|---------------------------|-------------------------|-------------------------|
| Buffer Cartridge ①    | 6 ea                      | 6 ea                    | 6 ea                    |
| Buffer Cartridge ②    | 6 ea                      | 6 ea                    | 6 ea                    |
| Disposable Filter Tip | 3 packs<br>(32 ea/pack)   | 3 packs<br>(32 ea/pack) | 3 packs<br>(32 ea/pack) |
| Elution Tube          | 8-tube strips x 12 ea     | 8-tube strips x 12 ea   | 8-tube strips x 12 ea   |
| User Guide            | 1 ea                      | 1 ea                    | 1 ea                    |

\* **Note:** All provided consumables including disposable tips, reaction tubes, and elution tubes are DNase- and RNase-free.

### Storage

The kit will maintain performance for at least two years under standard storage conditions.

Each Buffer Cartridge is hermetically sealed with a three-ply sealing foil and then wrapped in film to protect against leakage, evaporation, and cross-contamination. The Buffer Cartridges can be stored dry at room temperature (15-25°C) for up to 2 years from the date of delivery, provided they remain sealed.

The kit also contains lyophilized enzymes (Proteinase K and carrier RNA), which are pre-loaded into Buffer Cartridge. They can be stored at room temperature up to 2 years without any reduction in activity, provided they remain unopened. Once dissolved, enzymes should be stored at -20°C for up to 6 months.

## Specifications

| Sample Type          | Amount of Starting Sample |
|----------------------|---------------------------|
| Serum                | 200 µl                    |
| Plasma               |                           |
| CSF                  |                           |
| Urine                |                           |
| BAL                  |                           |
| Cell free body fluid |                           |
| Saliva               |                           |
| Swab                 | 1 ea                      |

\* **Note:** There may be differences in measured values depending on the type of samples.

## Precautions

- RNA is highly susceptible to degradation by exogenous RNase that may be introduced during the handling steps, all the steps must be conducted under sterile, RNase-free condition.
- RNase-free reagents, pipette tips, and tubes must be used with gloved hands while handling them.

## Introduction

### Product Description

*ExiPrep*<sup>™</sup> Plus Viral DNA/RNA Kits are designed for extraction of highly purified DNA or RNA from serum, plasma, Cerebrospinal fluid (CSF), Bronchoalveolar lavage (BAL), urine, swab, or other cell-free body fluids. *ExiPrep*<sup>™</sup> Plus Viral DNA/RNA Kits provide total solution for accurate and rapid total DNA/RNA extraction. These kits employ our unique Buffer Cartridge system. The Buffer Cartridges contain all components for nucleic acid extraction, including: binding buffer, washing buffer, elution buffer, and magnetic nanobead solution. The Buffer Cartridges are key to extract total DNA/RNA with the aid of *ExiPrep*<sup>™</sup> 16 Plus (Cat. No. A-5030). *ExiPrep*<sup>™</sup> 16 Plus is designed for rapid extraction of nucleic acids delivering up to 16 extracted samples automatically. RNA extracted through this kit can be used for a variety of applications, including: reverse transcription PCR (RT-PCR), reverse transcription quantitative PCR (RT-qPCR), Northern blot analysis, and cDNA synthesis.

## **Principle**

*ExiPrep™* Plus Viral DNA/RNA Kits use Magnetic Nano Beads to extract nucleic acid. Buffers within the kit assist nucleic acid to bind to silica-coated magnetic nanobeads. As a result, high yield and highly purified nucleic acid is extracted from samples.

The Buffer Cartridges consist of binding buffer, washing buffer, elution buffer, and magnetic nanobeads. Samples are lysed in the presence of a guanidine-thiocyanate-containing buffer, which is a highly denaturing agent and inactivates RNase to isolate RNA. Extracted DNA/RNA is bound to silica-coated magnetic nanobeads. Proteins and other contaminants are eliminated by subsequent washing, and highly purified DNA or RNA is eluted in an elution buffer.

## Features & Benefits

- Convenient & Rapid: Uses a pre-filled buffer cartridge system in which enzymes and reagents for nucleic acid extraction are dispensed.
- Reproducible: Uses fully automatic nucleic acid extraction equipment, and reproducible results can be obtained.
- Efficient: Contains all required consumables such as Disposable Filter Tips and Elution Tubes.
- Ready-to-use: Extracted DNA/RNA is ready-to-use for various application.

### Components of Buffer Cartridges










Figure 1. Position of reagents, tubes, and starting material in Cartridges/Elution Tube Rack of *ExiPrep*™ Plus Viral DNA/RNA Kit (K-4271), *ExiPrep*™ Plus Viral DNA Kit (K-4272), *ExiPrep*™ Plus Viral RNA Kit (K-4273).

## Experimental Procedures


### Loading the Kit to the Instrument

| Steps  | Procedure Details  |
|--|--|
| <p>Example 1) For 1 sample</p>  <p>Example 2) For 8 samples</p>  <p>Example 3) For 16 samples</p>  | <p>1. Punch holes in the sealing films of Buffer Cartridge ① and ② using 6-Hole Punch (<i>ExiPrep</i><sup>™</sup> 16 Plus's accessory) according to the number of samples.</p> <p><b>* Note:</b> Before punching the hole, agitate the Buffer Cartridge gently to settle the beads and buffer.</p> |
|   | <p>2. Load 200 µl of sample into 'Sample loading well' of Cartridge ①.</p> <p><b>* Note:</b> Be careful not to contaminate any other wells.</p>  |
|   | <p>3. Open the door of <i>ExiPrep</i><sup>™</sup> 16 Plus and pull out the baseplate completely.</p>   |




4. Place the Buffer Cartridge ② onto the proper position of the baseplate.

\* **Note:** Please check the punched holes of the Buffer Cartridge ②.




5. Place the Buffer Cartridge ① onto the proper position of the baseplate.




6. Place the Elution Tube Rack including Elution Tubes onto the proper position of the baseplate. The Elution Tube Rack is slotted so it can only be placed in the correct orientation.

\* **Note:** Make sure the direction of the Elution Tube caps laid out as on the left when inserting into the Elution Tube Rack.



7. Place the Disposable Filter Tip Rack onto the proper position of the base plate.

\* **Note:** Tips should be placed in the corresponding positions with the punched holes of the Cartridges.



8. Place the Waste tray between Buffer Cartridge ① and ②.

9. Finally, confirm holes in the cartridges and position of samples and tips. Push the baseplate completely until you hear the click sound, then close the door.



10. Turn on the *ExiPrep*™ 16 Plus.

11. In the MENU screen, press 'Start' button to select a proper protocol.



12. The PREP SETUP screen appears as shown in the left, press protocol number according to the protocol number list (pp. 12-14). Confirm following information displayed on the screen, and then press the 'Enter' button.



13. Select the desired elution volume from the touch screen.

14. Press the 'ok' button to move to the next step.




15. Ensure that everything is correctly installed following the CHECK LIST, then choose "ok".



16. In the Running Mode screen, ensure that the protocol name appears on the screen.

17. Press the "RUN" button to initiate DNA extraction.



18. “Work Completion” screen appears when the protocol is completed. Open the door and collect final DNA from the Elution Tubes.

19. Remove all components used in the experiment, and choose 1, 2, or ok.

**\* Note:** If you want to quit and press the ‘ok’ button, the UV lamp will be turned on automatically.

## Troubleshooting

| Problem                                    | Comments   |
|--|--|
| <p><b>Low DNA/RNA yield or purity</b></p>  | <ul style="list-style-type: none"> <li>• <b>You may have used too much (or too little) starting material.</b><br/>The yield is dependent on the sample type and amount of starting sample. Appropriate amount of starting sample should be used for efficient extraction of nucleic acid.</li> <li>• <b>The lysis may have been incomplete.</b><br/>Centrifuge completely to obtain clear lysate.</li> <li>• <b>Incomplete suspension of the magnetic nanobeads may decrease the DNA/RNA yield or purity.</b><br/>You should agitate the Buffer Cartridge ① before use.</li> </ul> |
| <p><b>Co-eluted magnetic nanobeads</b></p> | <p>Sometimes magnetic nanobeads are eluted with your viral DNA/RNA. Magnetic nanobeads in the eluate will not affect the performance of the DNA/RNA in downstream applications. Furthermore, magnetic nanobeads cannot bind DNA/RNA in elution buffer, though it may affect readings on a spectrophotometer. Magnetic nanobeads that are carried over can be easily separated by centrifugation at 13,000 rpm for 1 minute.</p>  |

## Appendix A

## List of Protocol Numbers

| Protocol Number | Target    | Sample Source        |
|-----------------|-----------|----------------------|
| No. 401         | Viral DNA | Whole blood          |
| No. 402         | Viral DNA | Animal tissue        |
| No. 412         | Viral DNA | Plasma               |
| No. 413         | Viral DNA | Serum                |
| No. 414         | Viral DNA | Buffy coat           |
| No. 415         | Viral DNA | Sputum               |
| No. 416         | Viral DNA | BAL                  |
| No. 417         | Viral DNA | Saliva               |
| No. 418         | Viral DNA | Swab                 |
| No. 419         | Viral DNA | Urine                |
| No. 420         | Viral DNA | Stool                |
| No. 421         | Viral DNA | Cell-free body fluid |
| No. 422         | Viral DNA | Pleural fluid        |
| No. 423         | Viral DNA | CSF                  |
| No. 424         | Viral DNA | EPS                  |
| No. 425         | Viral DNA | Respiratory sample   |
| No. 426         | Viral DNA | Amniotic fluid       |
| No. 427         | Viral DNA | Forensic sample      |
| No. 428         | Viral DNA | Bone marrow          |
| No. 429         | Viral DNA | Bone                 |
| No. 430         | Viral DNA | Dried blood spot     |
| No. 431         | Viral DNA | Soil                 |
| No. 432         | Viral DNA | Hair                 |
| No. 433         | Viral DNA | Cell supernatant     |

|         |               |                      |
|---------|---------------|----------------------|
| No. 501 | Viral RNA     | Whole blood          |
| No. 502 | Viral RNA     | Animal tissue        |
| No. 512 | Viral RNA     | Plasma               |
| No. 513 | Viral RNA     | Serum                |
| No. 514 | Viral RNA     | Buffy coat           |
| No. 515 | Viral RNA     | Sputum               |
| No. 516 | Viral RNA     | BAL                  |
| No. 517 | Viral RNA     | Saliva               |
| No. 518 | Viral RNA     | Swab                 |
| No. 519 | Viral RNA     | Urine                |
| No. 520 | Viral RNA     | Stool                |
| No. 521 | Viral RNA     | Cell-free body fluid |
| No. 522 | Viral RNA     | Pleural fluid        |
| No. 523 | Viral RNA     | CSF                  |
| No. 524 | Viral RNA     | EPS                  |
| No. 525 | Viral RNA     | Respiratory sample   |
| No. 526 | Viral RNA     | Amniotic fluid       |
| No. 527 | Viral RNA     | Forensic sample      |
| No. 528 | Viral RNA     | Bone marrow          |
| No. 529 | Viral RNA     | Bone                 |
| No. 530 | Viral RNA     | Dried blood spot     |
| No. 531 | Viral RNA     | Soil                 |
| No. 532 | Viral RNA     | Hair                 |
| No. 533 | Viral RNA     | Cell supernatant     |
| No. 601 | Viral DNA/RNA | Whole blood          |
| No. 602 | Viral DNA/RNA | Animal tissue        |
| No. 612 | Viral DNA/RNA | Plasma               |
| No. 613 | Viral DNA/RNA | Serum                |

|         |               |                      |
|---------|---------------|----------------------|
| No. 614 | Viral DNA/RNA | Buffy coat           |
| No. 615 | Viral DNA/RNA | Sputum               |
| No. 616 | Viral DNA/RNA | BAL                  |
| No. 617 | Viral DNA/RNA | Saliva               |
| No. 618 | Viral DNA/RNA | Swab                 |
| No. 619 | Viral DNA/RNA | Urine                |
| No. 620 | Viral DNA/RNA | Stool                |
| No. 621 | Viral DNA/RNA | Cell-free body fluid |
| No. 622 | Viral DNA/RNA | Pleural fluid        |
| No. 623 | Viral DNA/RNA | CSF                  |
| No. 624 | Viral DNA/RNA | EPS                  |
| No. 625 | Viral DNA/RNA | Respiratory sample   |
| No. 626 | Viral DNA/RNA | Amniotic fluid       |
| No. 627 | Viral DNA/RNA | Forensic sample      |
| No. 628 | Viral DNA/RNA | Bone marrow          |
| No. 629 | Viral DNA/RNA | Bone                 |
| No. 630 | Viral DNA/RNA | Dried blood spot     |
| No. 631 | Viral DNA/RNA | Soil                 |
| No. 632 | Viral DNA/RNA | Hair                 |
| No. 633 | Viral DNA/RNA | Cell supernatant     |

## Ordering Information

| Description  | Cat. No |
|--|---------|
| <i>ExiPrep</i> <sup>™</sup> Plus Viral DNA/RNA Kit | K-4271  |
| <i>ExiPrep</i> <sup>™</sup> Plus Viral DNA Kit     | K-4272  |
| <i>ExiPrep</i> <sup>™</sup> Plus Viral RNA Kit     | K-4273  |

## Related Products

| Description                         | Cat. No |
|-------------------------------------|---------|
| <i>ExiPrep</i> <sup>™</sup> 16 Plus | A-5030  |
| <i>ExiProgen</i> <sup>™</sup>       | A-5041  |

### Explanation of Symbols

|   |  |   |   |
|---|--|---|---|
|  <p>Batch Code</p>                   |  <p>Biological Risks</p>                        |  <p>Catalog Number</p> |  <p>Caution</p>      |
|  <p>Consult Instructions For Use</p> |  <p>Contains Sufficient for &lt;n&gt; tests</p> |  <p>Do not Re-use</p>  |  <p>Manufacturer</p> |
|  <p>Research Use Only</p>            |  <p>Temperature Limitation</p>                  |  <p>Use-by Date</p>    |   |

#### BIONEER Corporation - HQ

**Address** 8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302, Republic of Korea  
**E-mail** [sales@bioneer.co.kr](mailto:sales@bioneer.co.kr)  
**Web** [www.bioneer.com](http://www.bioneer.com)

#### BIONEER Global Center

**Address** 71, Techno 2-ro, Yuseong-gu, Daejeon, 34013, Republic of Korea  
**E-mail** [sales@bioneer.co.kr](mailto:sales@bioneer.co.kr)  
**Web** [www.bioneer.com](http://www.bioneer.com)

#### BIONEER R&D Center

**Address** Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro, Bundang-gu, Seongnam-si  
Gyeonggi-do, 13488, Republic of Korea  
**E-mail** [sales@bioneer.co.kr](mailto:sales@bioneer.co.kr)  
**Web** [www.bioneer.com](http://www.bioneer.com)

#### BIONEER Inc. - USA Branch

**Address** 155 Filbert St. Suite 216 Oakland, CA 94607, USA  
**E-mail** [order.usa@bioneer.com](mailto:order.usa@bioneer.com)  
**Web** [us.bioneer.com](http://us.bioneer.com)

#### BIONEER Corp. - European Branch

**Address** Ludwig-Erhard-Strasse 30-34, 65760 Eschborn, Germany  
**E-mail** [euinfo@bioneer.com](mailto:euinfo@bioneer.com)  
**Web** [www.bioneer.com](http://www.bioneer.com)