

ExiProgen™ EC-Bulk Protein Synthesis Kit (Cat. No. K-7340)

Step I Protein Expression

A. Preparation of Experiment

- 1) Prepare the 'Template DNA' for protein expression.

(Note) The 'Template DNA' must contain the following : T7 Promoter-Ribosome binding site (RBS)-Target gene (with His-tag)-T7 terminator.

Use 5-20 µg of plasmid DNA per kilobase or 30-50 µg of PCR product.

(Purity: $A_{260}/_{280} > 1.7$, $A_{260}/_{230} > 1.5$).

- 2) Take out the 'Expression Cartridge ①, ②' from Box ② and thaw them at room temperature.

- 3) Take out 'E. coli extract', 'Master Mix', and 'DEPC DW' from a Box ② and thaw them on ice.

(Note) Please make sure that all solutions are completely thawed. It takes about 3-4 hours.

- 4) Punch all holes in the sealing films of 'Expression Cartridge ①, ②'.

- 5) Install 'Bulk SECF device' from a Box ① on a 'Bulk Reactor (Cat. No. KA-7340-1)' and

fill the 'Bulk SECF device' with the 'Reaction Mixture' as below.

(Note) 'Bulk SECF device' is filled with solution containing 0.05% Na₂S.

Remove the solution in the device before use.

Reaction Mixture	
Template DNA	X ml
E. coli extract	5.34 ml
Master Mix	9.34 ml
DEPC DW	(5.32-X) ml
Total volume	20 ml



B. Protein Expression with ExiProgen™

- 1) Turn on the ExiProgen™ and tap 'Press to Start' button.


Wait until the [MENU] screen appears.

- 2) Open the door of ExiProgen™ and pull out the baseplate.

- 3) Load 'Cartridge ①, ②', 'Waste Tray', and 'Bulk Reactor'.

Follow in order: 'Cartridge ②' → 'Cartridge ①' → 'Waste Tray' → 'Bulk Reactor'

, Place 8 'Disposable Filter Tips' on row B of 'Disposable Tip Rack'.



(Note I) Make sure that the first and second lines of 'Cartridge ②' are firmly fixed on the heating block.

(Note II) There are 'silicon rings' embedded on both sides of the 'Cartridge ①' installation position. Place the left side first, and then the right side.

- 4) Push the baseplate in until you hear the click sound, then close the door.

- 5) Press [MENU] 'Start' → [PREP SETUP], '906 (Bulk Expression)', 'Enter' → [PREP SETUP, Elution Volume] 'OK' → [PREP SETUP, Reaction Temperature] '30°C', 'OK' → [CHECK LIST] 'OK' → [Running Mode] 'RUN'.



- 6) The [Work Completion] screen will appear once the protocol is completed.

It takes approximately 25 hours.

- 7) Expressed protein samples will be collected in the 'Bulk SECF device'.

ExiProgen™ EC-Bulk Protein Synthesis Kit (Cat. No. K-7340)

Step II Protein Purification

A. Preparation of Experiment

- 1) Take out 'Purification Cartridge ①, ②' from a **Box ①** and punch all holes in the sealing films.
- 2) Load the **1 ml/well (Total 16 ml)** expressed protein samples (Final sample of Step I) on the '**Sample Loading Well (16 wells)**' of 'Purification Cartridge ①'.
- 3) Install the new 'Bulk SECF device' from a **Box ①** on a 'Bulk Reactor (Cat. No. KA-7340-1)' and fill **20 ml of sterile distilled water in the device**.
(Note) 'Bulk SECF device' is filled with solution containing 0.05% NaN₃.


Remove the solution in the device before use.

B. Protein Purification with ExiProgen™

- 1-4) 1-4 steps are the same as 'Protein Expression' protocol. Refer to front page.
- 5) Press [MENU]'Start' → [PREP SETUP], **'907 (Bulk_Purification)'**, 'Enter' → [PREP SETUP, Elution Volume]'OK' → [PREP SETUP, Reaction Temperature]'OK' → [CHECK LIST]'OK' → [Running Mode]'RUN'.
- 6) The [Work Completion] screen will appear once the protocol is completed.
It takes approximately **15 hours**.
(Note) The '907' protocol includes [Ni-NTA affinity purification] and [Dialysis to storage buffer] processes. If you want to separate two processes, use [908 (Bulk_Ni-NTA purify)] and [909 (Bulk_Storage)] protocol.
- 7) **Purified protein samples in storage buffer will be collected in the 'Bulk SECF device'**.
(Note I) Purified protein samples may contain trace amounts of magnetic beads. The beads do not influence to proteins and can be easily removed through spin-down.
(Note II) Composition of Storage Buffer: 50 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.05% NaN₃, 50% Glycerol, pH 7.6.

Step III Analysis of Sample

A. Position of each sample (in 'Purification Cartridge ①, ②')



B. Sampling for SDS-PAGE analysis

- 1) Prepare loading mixture as shown in the table.

	Expression/Unbound/1 st washing sample	Purified protein/Bead sample
Sample	5 µl	15 µl
4X Loading dye	5 µl	5 µl
Sterile distilled water	10 µl	-
Total volume	20 µl	20 µl

(Note) Add 360-400 µl/well of sterile distilled water to the well containing the beads for suspension. Proceed with the sampling afterwards.

- 2) Incubate the samples at 95°C for 5-10 min.
- 3) Load 5 µl of each sample to a well of SDS-PAGE gel [10 x 8 (cm), 0.75 mm thick, 10-well].

※ For more information, visit our **website (www.bioneer.com)**.