


## ExiProgen™ His-tagged Protein Purification Kit (Cat.No. K-7220, K-7221)

### Step I Sample Preparation (e.g. *E. coli* cell lysate)

- 1) Harvest the cells, including the target proteins by centrifuging at 3,000 rpm for 10 min.
- 2) Resuspend the cells with lysis buffer (not provided).  
(Note I) It is recommended to use 10 ml of lysis buffer per 0.1 g of wet cells.  
(Note II) Example of lysis buffer: 20 mM Tris-HCl, 1 mM 2-Mercaptoethanol, pH 7.6
- 3) Lyse the cells on ice with a sonicator equipped with a microtip.
- 4) Centrifuge the cell lysate at 13,000 rpm for 1 min to separate supernatant and pellet.
- 5) Check the presence of target proteins in supernatant with SDS-PAGE.

### Step II Preparation of Experiment

- 1) Take out 'Cartridge ①, ②', 'Disposable Filter Tip', 'Elution Tube', and 'Protection Cover' from the kit box.
- 2) Punch holes in the sealing films of 'Cartridge ①, ②'. Add 700-750 µl of samples (e.g. *E. coli* cell lysate) to 'Sample Loading Well' of 'Cartridge ②'.  
(Note) The number of columns punched should match with of the samples.



- 3) Load 'Elution Tube' on the 'Elution Tube Rack' (*ExiProgen™* accessory).  
(Rows in alphabet: Elution Tube)  
Place the 'Protection cover' on the 'Elution Tube Rack'.



### Step III Protein Purification with *ExiProgen™*

- 1) Turn on the *ExiProgen™* and tap 'Press to Start' button.  
Wait until the [MENU] screen appears.
- 2) Open the door of *ExiProgen™* and pull out the baseplate.
- 3) Load 'Cartridge ①, ②', 'Waste Tray', and 'Elution Tube Rack'.  
Follow in order: 'Cartridge ②' → 'Cartridge ①' → 'Waste Tray' → 'Elution Tube Rack'.  
(Note I) Make sure that the first and second lines of 'Cartridge ②' are firmly fixed on the heating block.  
(Note II) There are 'silicon rings' embedded on both sides of the 'Cartridge ①' installation position. Place the left side first, and then the right side.
- 4) Place the 'Disposable Filter Tips' on row A or B of 'Disposable Tip Rack'.  
(Note) Tips should be placed in the same columns with the punched holes of the Cartridges.



- 5) Push the baseplate in until you hear the click sound, then close the door.

## ExiProgen™ His-tagged Protein Purification Kit (Cat.No. K-7220, K-7221)

- 6) Tap the following: 'Start' → [PREP SETUP] **901 (Protein\_Purification)**, 'Enter' → [PREP SETUP, Elution Volume] 'OK' → [PREP SETUP, Reaction Temperature] 'OK' → [CHECK LIST] 'OK' → [Running Mode] 'RUN'.



- 7) The [Work Completion] screen will appear once the protocol is completed. It takes approximately **2 hours**.
- 8) **Purified protein samples** in the elution buffer can be collected from the elution tubes in the rack, located at the **rows labeled with alphabets**.

### Step IV Analysis of Sample

#### A. Position of each sample (in 'Cartridge ②')



- Bead sample
- Unbound sample
- 1<sup>st</sup> washing sample

#### Unbound sample

: Samples not bound to Ni-NTA magnetic bead

#### 1<sup>st</sup> washing sample

: Samples after 1<sup>st</sup> washing step in purification process

#### Bead sample

: Used bead samples for purification

#### B. Sampling for SDS-PAGE analysis

- 1) Prepare the loading mixture as shown in the table. Incubate the samples at 95°C for 5-10 min.

	Unbound/1 <sup>st</sup> washing sample	Purified protein/Bead sample
Sample	5 µl	15 µl
4X Loading dye	5 µl	5 µl
Sterile distilled water	10 µl	-
Total volume	20 µl	20 µl

(Note) Add 200 µl of sterile distilled water to the well containing the beads for suspension. Proceed with the sampling afterwards.

- 2) Load each sample to the wells of SDS-PAGE gel [10 x 8 (cm), 0.75 mm thick, 10-well].

(Note) Unbound and 1<sup>st</sup> washing sample: 5 µl/well,

Purified proteins and Bead samples: 10 µl/well