

ExiProgen™ EC-Bulk Protein Synthesis Kit (Cat. No. K-7340)

Step I Protein Expression

A. Preparation of Experiment

- 1) Prepare the 'Template DNA' for protein expression.

Note) The template DNA should have the structure of "T7 Promoter-RBS-Target gene (with His-tag)-T7 terminator".

Use 5~20 µg of plasmid DNA(per kilobase) or 30~50µg of PCR product.

(Purity-A_{260/280} ; >1.7, A_{260/230} ; >1.5).

- 2) Take out 'Expression Cartridge ①, ②' from a **Box ②** and thaw them at the room temperature.

Note) Please make sure that all solutions are completely thawed. It takes about 3~4 hours.


- 3) Punch all holes in the sealing films of 'Expression Cartridge ① and ②'.

- 4) Take out 'E. coli extract', 'Master mix' and 'DEPC DW' from a **Box ②** and thaw them on ice.

- 5) Install 'Bulk SECF device from a **Box ①**' on a 'Bulk Reactor (Cat. No. KA-7340-1)' and **Fill the Reaction mixture into the device.**

Note) 'Bulk SECF device', is filled with distilled water containing 0.05% NaN₃. Remove the water in the device before use.

Reaction mixture	
Template DNA	X mL
E. coli extract	5.34 mL
Master mix	9.34 mL
DEPC DW	(5.32-X) mL
Total	20 mL



B. Protein Expression with ExiProgen™

- 1) Turn on the ExiProgen™ and press the 'Press to start' button.


A loading bar will appear on the screen. Afterwards, the [MENU] screen will appear.

- 2) Open the door of the ExiProgen™ and pull out its base plate.

- 3) Install 'Cartridge ①, ②', 'Waste tray', and 'Bulk Reactor'.

Follow the exact order: 'Cartridge ②' → 'Cartridge ①' → 'Waste tray' → 'Bulk Reactor'.

And place **8 tips** of 'Disposable filter tips' on **B row** of 'Disposable Tip Rack'.



Note I) Make sure that the first and second lines of 'Cartridge ②' are firmly fixed on the heating block.

Note II) There are 'silicon rings' embedded in both sides of the installation position of 'Cartridge ①'. Install the left side first and then the right side.

- 4) Push the Base plate completely, until you hear 'click' sound or feel clicking, and close the door.

- 5) Press [MENU] 'Start' → [PREP SETUP], **'906 (Bulk_Expression)'**, 'Enter' → [PREP SETUP, Elution Volume] 'OK' → [PREP SETUP, Reaction Temperature] '30 °C', 'OK' → [CHECK LIST] 'OK' → [Running Mode] 'Run'.



- 6) About **25 hours** later, [Work Completion] screen will appear when the protocol is Completed.

- 7) Expressed protein samples will be collected in the 'Bulk SECF device'.

ExiProgen™ EC-Bulk Protein Synthesis Kit (Cat. No. K-7340)

Step II Protein Purification

A. Preparation of Experiment


- 1) Take out 'Purification Cartridge ①, ②' from a **Box ①** and Punch all holes in the sealing films of them.
- 2) Load the **1 mL/well (Total 16 mL)** expressed protein samples (Final sample of Step I) on the '**Sample Loading Well (16 wells)**' of 'Purification Cartridge ①'.
- 3) Install 'Bulk SECF device' from a **Box ①** on a 'Bulk Reactor (Cat. No. KA-7340-1)' and **Fill 20 mL of sterile distilled water in the device.**
Note) 'Bulk SECF device', is filled with distilled water containing 0.05% NaN₃.
Remove the water in the device before use.

B. Protein Purification with ExiProgen™

- 1~4) 1~4 steps are the same as 'Protein Expression' protocol. Refer to front page.
- 5) Press [MENU] 'Start' → [PREP SETUP, '907 (Bulk_Purification)', 'Enter'] → [PREP SETUP, Elution Volume] 'OK' → [PREP SETUP, Reaction Temperature] 'OK' → [CHECK LIST] 'OK' → [Running Mode] 'Run'.
- 6) About **15 hours** later, [Work Completion] screen will appear when the protocol is completed.
Note) The '907' protocol includes [Ni-NTA affinity purification] and [Dialysis to storage buffer] processes. If you want to separate two processes, use [908 (Bulk_Ni-NTA purify)] and [909 (Bulk_Storage)] protocol.
- 7) Purified protein samples in storage buffer will be collected in the 'Bulk SECF device'.
Note I) Purified protein samples may be contain trace amount magnetic beads. The beads have no influence to proteins and can easily removed through spin-down.
Note II) Composition of Storage buffer : [50mM Tris-Cl(pH7.6), 100mM NaCl, 1mM DTT, 0.05% NaN₃, 50% Glycerol].

Step III Analysis of sample

A. Position of each sample (in Purification Cartridges)



B. Sampling for SDS-PAGE analysis

- 1) Prepare loading mixture as shown in the table.

	Expression/Unbound/1 st washing sample	Purified protein/Bead sample
Sample	5 μ L	15 μ L
4x loading dye	5 μ L	5 μ L
Sterile distilled water	10 μ L	-
Total	20 μ L	20 μ L

Note) Bead sample was resuspended in 500 μ L/well sterile distilled water.

- 2) Treat the samples at 95 °C, for 5~10 min.
- 3) Load 5 μ L of each sample to a well of SDS-PAGE gel (10x8(cm), 0.75mm thick, 10 wells).

※ For more information on how to use this kit, please see our [website](#).