
User's Guide

ExiProgen™ **EC-Disulfide Protein Synthesis Kit**

REF **K-7330**

***ExiProgen*[™] EC-Disulfide Protein Synthesis Kit**

User's Guide

Version No.: 2.1 (2016-02)

Please read all the information in booklet before using the kit

**Bioneer Corporation
8-11, Munpyeongseo-ro
Daedeuk-gu, Daejeon 34302
Korea**

Tel: +82-42-930-8777

Fax: +82-42-930-8688

Email: sales@bioneer.com

www.bioneer.com

Safety warning and Precaution

*ExiProgen*TM EC-Disulfide Protein Synthesis Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to BIONEER within 30 days.

Patent

*ExiProgen*TM and its kits are protected by the patents KR10-2011-0085824, PCT/KR2012/006715, and KR 10-2012-0138335.

Trademark

*ExiProgen*TM is trademark of Bioneer Corporation.
Copyright © 2012 by Bioneer Corporation, All rights reserved.

Contents

1. Introduction	2
2. The <i>ExiProgen™</i> protein synthesis system	3
3. <i>ExiProgen™</i> EC-Disulfide Protein Synthesis Kit	5
3.1. Overview	5
3.2. Preparation of template DNA	6
3.3. Kit contents and storage	7
3.4. Information of components	8
3.4.1. Cartridges	8
3.4.2. Other components	9
3.5. Experimental procedure	10
3.5.1. Preparation of experiment	10
3.5.2. Preparation of protein expression mixture	11
3.5.3. Protein synthesis with <i>ExiProgen™</i>	12
3.5.4. Analysis of sample.....	16
4. Maintenance	19
5. Troubleshooting guide	20
6. Appendix: pBIVT vectors	23
7. References	24
8. Related products	28

1. Introduction

Proteins are one of the essential components in the vital reactions in vivo and in vitro, with various functions as enzymes, hormones, and structural proteins. The research on the roles and structures of proteins has been done actively in the post-genomic era, starting with the production of proteins of interest by using recombinant DNA technology.

In order to produce recombinant proteins people generally transformed various host cells, such as *E. coli*, yeast, and mammalian cells with a vector that contains a piece of DNA necessary for the expression of a target protein. The transformed cells that demonstrate a stable expression of the recombinant protein can be selected by the selection marker co-inserted into the host cells with a target DNA. The recombinant proteins are produced during the culture process of the selected transformed cells and after the culture collected and purified from either the cell lysates or the culture solution. The whole process of cell-based protein synthesis comprising from recombination and transformation up to purification is usually a tedious and labor intensive job. Additionally for the toxic proteins to the host cells it is very hard to express those proteins in a cell-based system. Therefore those toxic proteins can be obtained through repeated trials and errors under various conditions, which takes several days or even several months.

To overcome these obstacles, a cell-free protein expression method and its related products have been developed to accomplish protein synthesis inside a reaction vessel in a short period of time. This method only requires to add cell lysate and template DNA, either of an expression vector or a PCR product, to protein expression solution that contains amino acids, nucleic acids and energy source in the reaction vessel. Then the target recombinant protein is expressed during a reaction at appropriate temperature. This method has the advantage of quite shorter processing time, over the cell based method, as well as of expression of toxic proteins to the cells.

Cell-free protein expression method has been used in functional research on new genes and comparison studies on functional activities of various mutant proteins. It also has broad applications such as research on the mechanism of protein-protein interaction and on the active site of proteins.


2. The *ExiProgen™* protein synthesis system

Bioneer has developed the *ExiProgen™* that combines cell-free protein expression system and highly efficient purification system. The *ExiProgen™* expresses and purifies various proteins within a day in a fully automated fashion and provides the target protein with high purity and amount of from micrograms up to milligrams. Furthermore, the *ExiProgen™* can extract DNA and RNA from various biological samples within an hour.

The *ExiProgen™* expressed a target protein from a template DNA and reagents provided by Bioneer and purifies the expressed protein using affinity reactions between the Ni-NTA magnetic nanoparticles and the 6x histidine tag attached at the end of the protein. The entire process of expression and purification is automated. The only requirement is to construct a user's own template DNA and to install a protein synthesis kit in the *ExiProgen™*.

Bioneer currently provides various protein synthesis kits that have been optimized for many different proteins. Therefore, researchers can simply execute a protocol already programmed in the *ExiProgen™* along with installing the kit of interest.

Bioneer's cell-free protein synthesis products are divided into three groups (a figure on page 4): 1) products for the generation of template DNA; 2) *AccuRapid™* products for manual expression of proteins; 3) *ExiProgen™* products for automated protein expression and purification using *ExiProgen™* instrument.


3. ExiProgen™ EC-Disulfide Protein Synthesis Kit

3.1. Overview

Cell-free protein expression method provides several advantages over cell-based method such as an expression of toxic proteins and a production of small quantities of proteins with ease and speed. However, in a cell-free system using *E. coli* extract, it is hard to produce proteins containing disulfide bonds.

ExiProgen™ EC-Disulfide Protein Synthesis Kit allows disulfide-bonded proteins to be produced in a cell-free system using a buffer solution that has high redox potential. This kit adopts Bioneer's stepwise exchange cell-free protein expression method (Figure below).

ExiProgen™ EC-Disulfide Protein Synthesis Kit can synthesize proteins containing up to nine disulfide bonds. With this kit, the entire processes from protein expression and purification to dialysis are conducted in a fully automated manner and the final products are harvested in a storage buffer.

3.2. Prepare of template DNA

A protein expression vector can be used as template DNA for the *ExiProgen™* EC-Disulfide Protein Synthesis Kit. Figure 1 shows the typical structure of template DNA for the *ExiProgen™*. The template DNA should have the structure of “T7 promoter – Ribosome binding site (RBS) – Target gene – T7 terminator”. The target gene should also have start codon and stop codon for protein translation, as well as 6x histidine tag at either 5’ – or 3’-terminal for the ensuing purification.

Figure 1. Structure of Template DNA

1) NH constructs (5’ end 6x histidine tagged template)

2) CH constructs (3’end 6x histidine tagged template)

For cloning of a target gene, one can use pBIVT vector from Bioneer, which is a protein expression vector for *in vitro* translation (Cat. No. K-7350). A detailed explanation on pBIVT vector is provided in Appendix of the manual. pK7, pIVEX, and pET-series vectors could be used for this purpose.

In order to obtain the maximum amount of synthesized protein, it is recommended that the DNA sequence of target gene be optimized to the codon usage for *E. coli*. Bioneer provides Gene Synthesis Service to a customer if the target gene is not optimized to *E. coli*. Throughout the service a customer can get codon-optimized target gene cloned into pBIVT vector (or other *in vitro* translation vector a customer selects). Please refer to Bioneer’s homepage (www.bioneer.com) for more information.

3.3. Kit contents and storage

ExiProgen™ EC-Disulfide Protein Synthesis Kit ①		
Cartridge ①	96 well x 1 ea	Store at 4 ~ 8℃
Dialysis tube	1 pack (16 ea/pack)	
Disposable filter tip	1 pack (8 ea/pack)	Store at Room temperature
Protection cover	1 ea	

ExiProgen™ EC-Disulfide Protein Synthesis Kit ②		
Cartridge ②	96 well x 1 ea	Store at -20~ -70℃
<i>E. coli</i> extract	8-tube strip (Yellow) x 1 ea	
Master mix	8-tube strip (Violet) x 1 ea	
DEPC DW	8-tube strip (White) x 1 ea	
Storage buffer	35 mL x 2 bottle	
Positive Control DNA	1.5 mL tube x 1 ea	

ExiProgen™ EC-Disulfide Protein Synthesis Kit consists of Kit ① and Kit ②.


The Kit ① includes Cartridge ①, Disposable filter tip, Dialysis tube, and a Protection cover. Cartridge ① contains buffers to be used for protein purification and should be stored in a refrigerator (4–8℃) along with Dialysis tube.

The Kit ② includes Cartridge ②, *E. coli* extract, Master mix, DEPC DW, Storage buffer and a tube of Positive Control DNA. Cartridge ② contains Feeding solution for protein expression. The Kit ② should be stored at -20℃ or below, and *E. coli* extract is strongly recommended to locate at a -70℃ deep freezer.

The Cartridges of this kit are covered with 3-ply sealing films in order to prevent cross-contamination, evaporation, or leakage of solutions inside. All of the plastic products and buffers in this kit are provided under Nuclease- and Protease-free condition. Please be careful not to contaminate any part of the kit with nucleases or proteases.

3.4. Information of components

3.4.1. Cartridges

3.4.2. Other components

A. *E. coli* extract

The *E. coli* extract in the Kit ② is a cell lysate which includes T7 RNA polymerase, ribosome, tRNA and disulfide bond isomerase required for protein expression. It is provided as an 8-tube strip, and each tube contains **130 μ L**.

B. Master mix

The Master mix, which is also provided in the Kit ②, includes NTPs, amino acids, energy sources and salts required for protein expression. It is provided as an 8-tube strip, and each tube contains **220 μ L**.

C. DEPC DW

It is used to adjust final volume of protein expression solution. It is provided as an 8-tube strip, and each tube contains **130 μ L**.

D. Dialysis tube

Dialysis tube, which is specifically manufactured by Bioneer, is used to supply the protein source materials from the Feeding buffer to the protein expression solution, and to exchange the purification buffer to Storage buffer after purification of a target protein. It is provided in 20% ethanol and individually wrapped can be used as many as needed. It is required to remove 20% ethanol and rinse out with distilled water once before use.

E. Storage buffer

Composition of Storage buffer in the Kit ② is as follows;

50mM Tris-Cl (pH7.6), 100mM NaCl, 1mM DTT, 0.05%(v/v) Na ₂ S ₂ O ₃ , 50%(v/v) glycerol
--

Please add 1 ml of the storage buffer provided in two separate bottles in the kit to each well of G1 to I2 in Cartridge ② right before use. When necessary, a user may construct another buffer with different composition and use it instead of the storage buffer provided by the manufacturer. Just keep in mind that a user's own buffer should contain 10% of glycerol not to overflow out of the Reaction block.


3.5. Experimental procedure

3.5.1. Preparation of experiment

Supplies

; latex gloves, forceps, sterile distilled water in squeeze bottle, scissors, micro pipettes, Reaction block (provided with the *ExiProgen™* instrument)

- ① Take out the Cartridge ② from a freezer and thaw it at room temperature. (Please ensure that all solutions are completely thawed before the experiment.)
- ② Get the Cartridge ① in the Kit ① out of a refrigerator and leave it at room temperature.
- ③ Prepare the Reaction block.
- ④ Take out twice as many Dialysis tubes as the number of samples, remove the 20% ethanol solution and rinse out the inside and the outside of the tube with sterile distilled water.
- ⑤ Install the DW-rinsed Dialysis tube on a Reaction Block as the diagram below indicates. **Fill 500 uL of sterile distilled water in the Dialysis tube located in the row B.**

- ⑥ Take Disposable filter tip and a Protection cover out of the Kit box ①.

3.5.2. Preparation of protein expression mixture

The first step of protein synthesis using the *ExiProgen™* is to prepare the protein expression solution containing template DNA.

- ① Cut off quickly 8-tube strip of *E. coli* extract, Master mix, and DEPC DW as many as the number of samples and thaw them on ice.
- ② When using Positive Control DNA, take out the Positive Control DNA tube from the kit② and thaw it on ice.

Note) The scFv (recombinant single chain variable fragment) of 34 kDa size is provided as a Positive Control DNA.

- ③ Prepare protein expression solution using pre-thawed *E. coli* extract, Master mix and DEPC DW. It is necessary to spin-down each tube with a micro centrifuge and then mix with a pipette before use.

	Positive Control	Sample
Template DNA	2 uL	X uL
<i>E. coli</i> extract	120 uL	120 uL
Master mix	210 uL	210 uL
DEPC DW	118 uL	(120-X) uL
Total	450 uL	450 uL

Note) The amount of template DNA, in microgram quantity, can be determined as follows:

Use 1 µg of plasmid DNA per kb (kilo base) of the insert size of DNA. (For example, 1.5 µg of Plasmid DNA is appropriate for 1.5 kb-sized insert).


- ④ Remove water completely from the inside of Dialysis tube in the **row A of the reaction block**, and add protein expression solution.
- ⑤ Cover the Protection Cover on a Reaction Block.

3.5.3. Protein synthesis with *ExiProgen™*


- ① Punch holes in the sealing films of Buffer Cartridge ① and ② using a 6 hole punch before the installation of the Cartridges on the *ExiProgen™* instrument, and add the Storage buffer on G1~I2 rows of Cartridge ② (Refer to the page 8~9.)

Note) The number of columns punched should be matched with the sample number.

Example1) For 1 sample

Example 2) For 8 sample

- ② Install the Buffer Cartridges and related components on the *ExiProgen™* instrument as follows:

1. Open the door of *ExiProgen™* instrument and pull out the Base Plate completely.

2. Install the Cartridge ② in the position of ② on the Base Plate.

* Ensure that Cartridge ② is firmly installed without any movement and that the row L of the Cartridge ② is located on the top of the heating block of the Base Plate.

	<p>3. Install the Cartridge ① in the position of ① on the Base Plate.</p> <p>* There are 'silicon rings' embedded in both sides of the installation position of Cartridge ①. Install the left side of the Cartridge ① first and then the right side. Ensure that the Cartridge ① is firmly installed, without any</p>
	<p>4. Install the Waste Tray between the Cartridge ② and ①.</p> <p>* Follow the exact installation order of Cartridge ② → Cartridge ① → Waste Tray. Ensure that the Cartridges and the Waste Tray are firmly installed, without any movement.</p>
	<p>5. Install the Reaction Block on the magnetic part of the Base Plate (see the left figure).</p> <p>* Ensure that the row A of the Reaction block is located toward the two Cartridges (inside) and the row B toward a user (outside).</p>
	<p>6. Load as many Disposable tip(s) as the number of samples into the same location of column with punched hole(s) in the <u>row B of the Disposable Tip Rack</u></p> <p>* Ensure that 1) tips are located at the same columns with the punched holes of the Cartridges, and 2) there should be no tip at the columns which will not be used.</p>
	<p>7. Push the Base Plate completely, until you hear 'click' sound or feel clicking, and close the door of ExiProgen™.</p>


- ③ After completion of setup, run the *ExiProgen™* instrument as follows.

 <p style="text-align: center;">↓</p> 	<p>1. Turn on the <i>ExiProgen™</i> instrument and press the 'Press to start' button in the center of the screen.</p> <p>Soon the <i>ExiProgen™</i> screen with moving scroll bar appears, and after a short period of time it moves to the "MENU" screen.</p> <p>The short pause is for the initialization process of the instrument. If it does not move to the next screen, turn off the instrument and call the customer center for A/S.</p>
	<p>2. In the MENU screen, press 'Start' button to select a proper protocol.</p>
	<p>3. In the PREP SETUP screen shown as left, input '<u>905</u>' of protocol number. Ensure that the statement</p> <p>Prep type: Protein Sample SRC: Synthesis_DS</p> <p>appears on the screen. Then select 'Enter'.</p>

	<p>4. Choose 'ok' to go to the next step.</p> <p>* This step is for nucleic acid extraction only.</p>
	<p>5. Choose '<u>26°C</u>' and "ok" to go to the next step.</p>
	<p>6. Ensure that everything is correctly installed following the CHECK LIST, then choose "ok". If not, make an appropriate change.</p>
 	<p>7. In the Running Mode screen shown as left, ensure that the statement Prep type: Protein Sample SRC: Synthesis_DS appears on the screen. Then press "RUN" to initiate protein synthesis.</p> <p>"Work Completion" screen appears when the Disulfide protocol is completed. Open the door and collect purified proteins. Remove all components used in the experiment, and choose 1, 2, or OK;</p> <p>* UV lamp will be on automatically if 'ok' button is selected to finish the use of the ExiProgen™.</p>


3.5.4. Analysis of sample

After the synthesis of protein using the *ExiProgen™* instrument, the final target protein in about 120~130 μ L of Storage buffer is collected from the Dialysis tube in the **row B of Reaction Block**. (It may contain trace amount of Ni-NTA magnetic beads, which has no influence to the protein and can be easily removed through centrifugation at 10,000 rpm for 1 min).

Some of reaction samples from protein expression and/or purification processes could also be collected from the wells of row L1, K1, or K2 at the Cartridge ② as shown in Figure 2.

Figure 2. Samples at the selected rows of Cartridge ②


Through SDS-PAGE or other methods a user is able to check whether the desired protein is synthesized accordingly. Here is the protocol for SDS-PAGE:

- a. Prepare loading mixture as shown in the table.

	Final target protein solution	Expression/Unbound/1 st washing samples
sample	10 μ l	5 μ l
4x loading dye	5 μ l	5 μ l
Sterile DW	5 μ l	10 μ l
Total	20 μ l	20 μ l

- b. Treat the samples at 95 °C for 5–10 min.
- c. Prepare SDS-PAGE mini-gel either 10% or 12% SDS-PAGE gel (10x8 (cm), 10 wells) is typically used.
Load 10 μ l of “Final target protein solution” and 5 μ l each of “Expression/Unbound samples/1st washing samples” to a well and run the SDS-PAGE.
- d. After staining of the gel with Coomassie blue solution, a band of target protein is detected as in Figure 3.


Figure 3. Final target protein on SDS-PAGE

M: *AccuLadder™* Protein Size Marker (Low), 1: scFv (34kDa) 2: rtPA (42kDa)

Note) A protein band shown around 24 kDa on the gel is a disulfide bond isomerase that does not interfere in the activity of the target protein.

Figure 4. Activity of the final protein produced by the *ExiProgen™* EC-Disulfide Protein Synthesis Kit

Note) The graph shown above is the results of activity measurements of rPA (recombinant plasminogen activator). rtPA was incubated with S-2288™, a substrate from *Chromogenix*, and the absorbance at 405nm was detected in every minute for 10 minutes.

4. Maintenance

After protein synthesis, wash and store the Reaction Block and other accessories as follows.

A. Reaction Block

Rinse all the used wells of the Reaction Block with sterile distilled water. Remove any residual water completely and store it in 70% ethanol solution.

B. Waste tray and Disposable tip rack

Discard all the solution in the waste tray, wash it in the running water and rinse it with 20% ethanol. If there remains any dirty sludge on the Disposable tip rack, clean it with 20% ethanol solution. And store at room temperature.

C. Cartridge ①, ②

The Cartridges with unused wells should be covered with their lid and be stored in a refrigerator (Cartridge ①) and in a freezer (Cartridge ②).

5. Troubleshooting guide

Please refer to the troubleshooting guide below if there is a problem in protein synthesis/purification using Bioneer's *ExiProgen™* protein synthesis kits and the *ExiProgen™* instrument. Keep in mind those solutions in the guide are for most of proteins, but not for all proteins.

1. No protein synthesis, including positive control protein

Causes	Solutions
Contamination of nucleases (DNase or RNase)	Wear gloves at all times. Use DNase- & RNase-free pipette tips.
Pipetting error or no reagent added	Check accuracy of pipettes. Add exact amount of DNA and/or reagents. Mix every reagents exactly according to the protocol.
Inappropriate storage of reagents	Store every reagents/components in the kit at the recommended temperature. No repeated freeze-and-thaw of <i>E. coli</i> extract.

2. No protein synthesis, except positive control protein

Causes	Solutions
Error(s) in gene sequence of Template DNA	Check that the target gene contains start codon (ATG) and stop codon (TAA, TAG, or TAG) at the right position. Check that the target gene codon is in-frame. Check if there is a mutation in the ORF (Open Reading Frame) of template DNA (Translation may be stopped.).

Structure of expression vector	Check that the gene of a target protein is cloned into the expression vector with T7 expression system containing T7 promoter, T7 terminator, and RBS. Add appropriate amount of IPTG (Isopropyl β -D-1-thiogalactopyranoside) to the reaction mixture if the expression vector can produce endogenous lac repressor.
Contamination of template DNA	Prepare and use fresh, uncontaminated template DNA.

3. Low yield of target protein

Causes	Solutions
Low purity of template DNA	Use template DNA of which $A_{260/280}$ is between 1.7 and 2.0, and $A_{260/230}$ is over 1.5.
Concentration of template DNA	Determine and use the optimal concentration of template DNA (see Appendix 2.).
No codon optimization of template DNA to <i>E. coli</i>	Do the codon optimization of target gene to <i>E. coli</i> .
Location of histidine tag	<p>Try to locate the histidine tag on the different location of target protein.</p> <p>* Histidine tag may inhibit expression of target protein due to its location. The expressed target protein may not be purified through Ni-NTA beads, if histidine tag is not exposed outside effectively.</p>


4. Low activity and/or solubility of target protein

Causes	Solutions
Requiring posttranslational modification for activity and/or solubility	Not available yet, due to the lack of posttranslational modification in cell-free protein synthesis using <i>E. coli</i>
Requiring cofactor for protein activity	Add the cofactor required to the final target protein.
Aggregation of protein due to low solubility	Synthesize protein at lower temperature. Add chaperone protein to help protein folding.

6. Appendix: pBIVT vectors

pBIVT vector set (Cat. No. K-7350) is designed for Bioneer's cell-free protein synthesis system. It fits well into all kits of *ExiProgen™* protein synthesis system.

Using this vector set it is possible to add 6x histidine tag to either N-terminal or C-terminal of target protein. The structure of vector is as follows.

① pBIVT-1's MCS (Multi Cloning Site) sequences

CATATGCACCAACCACCAACCACCGGATCCGAGCTCAAGCTTGCGGCCGC

Nde I 6x histidine tag BamH I Sac I Not I

ATAGGTCGAC

Sal I

② pBIVT-2's MCS (Multi Cloning Site) sequences

CATATGGGATCCGAGCTCAAGCTTGCGGCCGCACACCAACCACCAACC

Nde I BamH I Sac I Not I 6x histidine tag

ACTAGGTCGAC

Sal I

7. References

Kim HC, Kim TW, Kim DM (2011) Prolonged production of proteins in a cell free protein synthesis system using polymeric carbohydrates as an energy source. *Process Biochemistry* 46: 1366–1369.

Chursov A, Walter MC, Schmidt T, Mironov A, Shneider A, et al. (2011) Sequence–structure relationships in yeast mRNAs. *Nucleic Acids Res* In press.

Trotta E (2011) The 3–base periodicity and codon usage of coding sequences are correlated with gene expression at the level of transcription elongation. *PLoS ONE* 6: e21590.

Oh IS, Lee JC, Lee MS, Chung JH, Kim DM (2010) Cell-free production of functional antibody fragments. *Bioproc. Biosyst. Eng.*, 33(1): 127–132

Park CG, Kwon MA, Song JK, Kim DM (2010) Cell-free synthesis and multi-fold screening of *Candida antarctica* lipase B (CalB) variants after combinatorial mutagenesis of hot spots. *Biotechnol. Prog.*, 27(1): 47–53.

Park S, Hamad–Schifferli K (2010) Enhancement of in vitro translation by gold nanoparticle–DNA conjugates. *ACS Nano* 4: 2555–2560.

Park CG, Ahn MA, Song JK, Kim DM (2009) Expression of functional *Candida antarctica* lipase B in a cell-free protein synthesis system derived from *Escherichia coli*. *Biotechnol. Prog.*, 25(2): 589–593.

Keum JW, Ahn JH, Kang TJ, Kim DM (2009) Combinatorial, selective and reversible control of gene expression using oligodeoxynucleotides in a cell-free protein synthesis system. *Biotechnol Bioeng* 102: 577–582.

Ahn JH, Kang TJ, Kim DM (2008) Tuning the expression level of recombinant proteins by modulating mRNA stability in a cell-free protein synthesis system. *Biotechnol Bioeng* 101: 422–427.

Ahn JH, Keum JW, Kim DM (2008) High-throughput, combinatorial engineering of initial codons for tunable expression of recombinant proteins. *J Proteome Res* 7: 2107–2113.

Hino M, Kataoka M, Kajimoto K, Yamamoto T, Kido J, Shinohara Y, Baba Y (2008)

Efficiency of cell-free protein synthesis based on a crude cell extract from *Escherichia coli*, wheat germ, and rabbit reticulocytes. *J Biotechnol*, 133(2):183–189.

Oh IS, Kim TW, Ahn JH, Keum JW, Choi CY, Kim DM (2007) Use of L-buthionine sulfoximine for the efficient expression of disulfide-containing proteins in cell-free extracts of *Escherichia coli*. *Biotechnol Bioproc Eng*, 12: 574–578

Kim TW, Oh IS, Keum JW, Kwon YC, Byun JY, et al. (2007) Prolonged cell free protein synthesis using dual energy sources: Combined use of creatine phosphate and glucose for the efficient supply of ATP and retarded accumulation of phosphate. *Biotechnol Bioeng* 97: 1510–1515.

Ohashi H, Shimizu Y, Ying BW, Ueda T (2007) Efficient protein selection based on ribosome display system with purified components. *Biochem Biophys Res Commun*, 352(1):270–276.

Oh IS, Kim DM, Kim TW, Park CG, Choi CY (2006) Providing an oxidizing environment for the cell-free expression of disulfide-containing proteins by exhausting the reducing activity of *Escherichia coli* S30 extract. *Biotechnol Prog*, 22: 1225–1228

Keum JW, Ahn JH, Choi CY, Lee KH, Kwon YC, et al. (2006) The presence of a common downstream box enables the simultaneous expression of multiple proteins in an *E. coli* extract. *Biochem Biophys Res Commun* 350: 562–567.

Kim TW, Keum JW, Oh IS, Choi CY, Park CG, et al. (2006) Simple procedures for the construction of a robust and cost-effective cell-free protein synthesis system. *J Biotechnol* 126: 554–561.

Swartz J (2006) Developing cell-free biology for industrial applications. *J Ind Microbiol Biotechnol*, 33(7):476–485.

Villemagne D, Jackson R, Douthwaite JA (2006) Highly efficient ribosome display selection by use of purified components for in vitro translation. *J Immunol Methods*, 313(1–2):140–148.

Katzen F, Chang G, Kudlicki W (2005) The past, present and future of cell-free protein synthesis. *Trends Biotechnol*, 23(3):150–156.

Josephson K, Hartman MC, Szostak JW (2005) Ribosomal synthesis of unnatural

peptides. *J Am Chem Soc*, 127(33):11727–11735.

Shimizu Y, Kanamori T, Ueda T (2005) Protein synthesis by pure translation systems. *Methods*, 36(3):299–304.

Underwood KA, Swartz JR, Puglisi JD (2005) Quantitative polysome analysis identifies limitations in bacterial cell-free protein synthesis. *Biotechnol Bioeng*, 91(4):425–435.

Forster AC, Cornish VW, Blacklow SC (2004) Pure translation display. *Anal Biochem* 2004, 333(2):358–364.

Torizawa T, Shimizu M, Taoka M, Miyano H, Kainosho M (2004) Efficient production of isotopically labeled proteins by cell-free synthesis: a practical protocol. *J Biomol NMR* 30: 311–325.

Ying BW, Taguchi H, Ueda H, Ueda T (2004) Chaperone-assisted folding of a single-chain antibody in a reconstituted translation system. *Biochem Biophys Res Commun* 2004, 320(4):1359–1364.

Gualerzi CO, Giuliodori AM, Pon CL (2003) Transcriptional and posttranscriptional control of cold-shock genes. *J Mol Biol*, 331(3):527–539.

Rungpragayphan S, Nakano H, Yamane T (2003) PCR-linked in vitro expression: a novel system for high-throughput construction and screening of protein libraries. *FEBS Lett* 540: 147–150.

Kigawa T, Yamaguchi-Nunokawa E, Kodama K, Matsuda T, Yabuki T, et al. (2002) Selenomethionine incorporation into a protein by cell-free synthesis. *J Struct Funct Genomics* 2: 29–35.

Frydman J (2001) Folding of newly translated proteins *in vivo*: the role of molecular chaperones. *Annu Rev Biochem*, 70:603–647.

Shimizu Y, Inoue A, Tomari Y, Suzuki T, Yokogawa T, Nishikawa K, Ueda T (2001) Cell-free translation reconstituted with purified components. *Nat Biotechnol*, 19(8):751–755.

Ramachandiran V, Kramer G, Hardesty B (2000) Expression of different coding sequences in cell-free bacterial and eukaryotic systems indicates translational pausing on *Escherichia coli* ribosomes. *FEBS Lett*, 482(3):185–188.

Jermutus L, Ryabova LA, Pluckthun A (1998) Recent advances in producing and selecting functional proteins by using cell-free translation. *Curr Opin Biotechnol* , 9(5):534–548.

Netzer WJ, Hartl FU (1997) Recombination of protein domains facilitated by co-translational folding in eukaryotes. *Nature* , 388(6640):343–349.

Roberts RW, Szostak JW (1997) RNA-peptide fusions for the *in vitro* selection of peptides and proteins. *Proc Natl Acad Sci U S A*, 94(23):12297–12302.

Frydman J, Hartl FU (1996) Principles of chaperone-assisted protein folding: differences between *in vitro* and *in vivo* mechanisms. *Science* , 272(5267):1497–1502.

Kuriki Y (1986) Stimulation *in vitro* of expression of the amp gene of pBR322 by soluble protein fractions isolated from *E. coli*. *Biochem Int* , 12(4):593–602.

Zubay G (1973) *In vitro* synthesis of protein in microbial systems. *Annu Rev Genet*, 7:267–287.

8. Related products

Cat. No.	Product	Size
K-7250	<i>AccuRapid™</i> Cell-Free Protein Expression Kit	45 μ L x 24 reactions
K-7260	<i>AccuRapid™</i> Midi Protein Expression Kit	1 mL x 5 reactions
K-7270	<i>AccuRapid™</i> Maxi Protein Expression Kit	10 mL x 1 reaction
K-7280	<i>AccuRapid™</i> Protein Synthesis Kit	5 reactions
K-7300	<i>ExiProgen™</i> EC Protein Synthesis Kit	16 reactions
K-7301		32 reactions
K-7302		96 reactions
K-7310	<i>ExiProgen™</i> EC-Maxi Protein Synthesis Kit	8 reactions
K-7320	<i>ExiProgen™</i> EC-Tagfree Protein Synthesis Kit	8 reactions
K-7350	pBIVT vector set	Each 5 μ g
K-7400	<i>ExiProgen™</i> ProXpress PCR Template Kit	16 reactions
K-7401		32 reactions
D-2010	<i>AccuLadder™</i> Protein Size Marker (Broad)	500 μ L
D-2020	<i>AccuLadder™</i> Protein Size Marker (Low)	500 μ L
S-2041	Gene Synthesis Service	–
S-2500	Protein Synthesis Service	–
A-5041	<i>ExiProgen™</i>	–

ExiProgen™ EC-Disulfide Protein Synthesis Kit

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Version No.: 2.1 (2016-02)

Please read all the information in booklet before using the kit

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Safety warning and Precaution

ExiProgen[™] EC–Disulfide Protein Synthesis Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to BIONEER within 30 days.

Patent

ExiProgen[™] and its kits are protected by the patents KR10–2011–0085824, PCT/KR2012/006715, and KR 10–2012–0138335.

Trademark

ExiProgen[™] is trademark of Bioneer Corporation.
Copyright © 2012 by Bioneer Corporation, All rights reserved.

Contents

1. Introduction	2
2. The <i>ExiProgen™</i> protein synthesis system	3
3. <i>ExiProgen™</i> EC-Disulfide Protein Synthesis Kit	5
3.1. Overview	5
3.2. Preparation of template DNA	6
3.3. Kit contents and storage	7
3.4. Information of components	8
3.4.1. Cartridges	8
3.4.2. Other components	9
3.5. Experimental procedure	10
3.5.1. Preparation of experiment	10
3.5.2. Preparation of protein expression mixture	11
3.5.3. Protein synthesis with <i>ExiProgen™</i>	12
3.5.4. Analysis of sample.....	16
4. Maintenance	19
5. Troubleshooting guide	20
6. Appendix: pBIVT vectors	23
7. References	24
8. Related products	28

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2. The *ExiProgen*™ protein synthesis system

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Park CG, Kwon MA, Song JK, Kim DM (2010) Cell-free synthesis and multi-fold screening of *Candida antarctica* lipase B (CalB) variants after combinatorial mutagenesis of hot spots. *Biotechnol Prog.*, 27(1): 47–53.

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Hino M, Kataoka M, Kajimoto K, Yamamoto T, Kido J, Shinohara Y, Baba Y (2008) Efficiency of cell-free protein synthesis based on a crude cell extract from

Escherichia coli, wheat germ, and rabbit reticulocytes. *J Biotechnol* , 133(2):183–189.

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Ohashi H, Shimizu Y, Ying BW, Ueda T (2007) Efficient protein selection based on ribosome display system with purified components. *Biochem Biophys Res Commun* , 352(1):270–276.

Oh IS, Kim DM, Kim TW, Park CG, Choi CY (2006) Providing an oxidizing environment for the cell-free expression of disulfide-containing proteins by exhausting the reducing activity of *Escherichia coli* S30 extract. *Biotechnol Prog* , 22: 1225–1228

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Katzen F, Chang G, Kudlicki W (2005) The past, present and future of cell-free protein synthesis. *Trends Biotechnol* , 23(3):150–156.

Josephson K, Hartman MC, Szostak JW (2005) Ribosomal synthesis of unnatural peptides. *J Am Chem Soc*, 127(33):11727–11735.

Shimizu Y, Kanamori T, Ueda T (2005) Protein synthesis by pure translation systems. *Methods* , 36(3):299–304.

Underwood KA, Swartz JR, Puglisi JD (2005) Quantitative polysome analysis identifies limitations in bacterial cell-free protein synthesis. *Biotechnol Bioeng* , 91(4):425–435.

Forster AC, Cornish VW, Blacklow SC (2004) Pure translation display. *Anal Biochem* 2004, 333(2):358–364.

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Frydman J (2001) Folding of newly translated proteins *in vivo*: the role of molecular chaperones. *Annu Rev Biochem* , 70:603–647.

Shimizu Y, Inoue A, Tomari Y, Suzuki T, Yokogawa T, Nishikawa K, Ueda T (2001) Cell-free translation reconstituted with purified components. *Nat Biotechnol* , 19(8):751–755.

Ramachandiran V, Kramer G, Hardesty B (2000) Expression of different coding sequences in cell-free bacterial and eukaryotic systems indicates translational pausing on *Escherichia coli* ribosomes. *FEBS Lett* , 482(3):185–188.

Jermutus L, Ryabova LA, Pluckthun A (1998) Recent advances in producing and selecting functional proteins by using cell-free translation. *Curr Opin Biotechnol* ,

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8. Related products

Cat. No.	Product	Size
K-7250	<i>AccuRapid™</i> Cell-Free Protein Expression Kit	45 uL x 24 reactions
K-7260	<i>AccuRapid™</i> Midi Protein Expression Kit	1 mL x 5 reactions
K-7270	<i>AccuRapid™</i> Maxi Protein Expression Kit	10 mL x 1 reaction
K-7280	<i>AccuRapid™</i> Protein Synthesis Kit	5 reactions
K-7300	<i>ExiProgen™</i> EC Protein Synthesis Kit	16 reactions
K-7301		32 reactions
K-7302		96 reactions
K-7310	<i>ExiProgen™</i> EC-Maxi Protein Synthesis Kit	8 reactions
K-7320	<i>ExiProgen™</i> EC-Tagfree Protein Synthesis Kit	8 reactions
K-7350	pBIVT vector set	Each 5ug
K-7400	<i>ExiProgen™</i> ProXpress PCR Template Kit	16 reactions
K-7401		32 reactions
D-2010	<i>AccuLadder™</i> Protein Size Marker (Broad)	500 uL
D-2020	<i>AccuLadder™</i> Protein Size Marker (Low)	500 uL
S-2041	Gene Synthesis Service	–
S-2500	Protein Synthesis Service	–
A-5041	<i>ExiProgen™</i>	–

● Bioneer Worldwide

Bioneer Corporation

Address 8-11 Munpyeongseoro, Daedeok-gu, Daejeon, 34302, Republic of Korea
Tel +82-42-930-8777 (Korea: 1588-9788)
Fax +82-42-930-8688
E-mail sales@bioneer.com
Web www.bioneer.com

Bioneer Inc.

Address 1301 Marina Village PKWY, Suite 110, Alameda, CA 94501, USA
Tel +1-877-264-4300 (Toll free)
Fax +1-510-865-0350
E-mail order.usa@bioneer.com
Web us.bioneer.com

Bioneer R&D Center

Address Korea Bio Park BLDG #B702, 700 Daewangpangyoro, Bundang-gu, Seongnam-si
Gyeonggi-do, 13488, Republic of Korea
Tel +82-31-628-0500
Fax +82-31-628-0555
E-mail sales@bioneer.co.kr
Web www.bioneer.co.kr