

**AccuPrep® Universal RNA Extraction Kit (K-3140, K-3141)****I Before You Begin**

- 1) Add 10  $\mu$ l  $\beta$ -mercaptoethanol per 1 ml RB Buffer.
- 2) Prepare additional ethanol(80% and 100%) that is not included.

**II Cultured Cell Collection**

- 1) Cells grown in suspension:  
Count the cell number and centrifuge the cultured cells ( $10^4 \sim 10^6$ ) at 300  $\times g$  for 5 min. Discard supernatant and go to step III. RNA Extraction from Cultured Cell.
- 2) Cells grown in a monolayer: There are 2 different ways to collect cells grown in a monolayer.
  - a. Direct cell lysis from culture dish: Completely remove cell culture medium and go to step III. RNA Extraction from Cultured Cell.
  - b. Harvesting cells with trypsin: Remove cell culture medium and wash the monolayer with DPBS. Add 0.1%~0.25% trypsin to the washed cell monolayer. When the cells are detached, add cell culture medium to inactivate trypsin. Transfer the cells into a RNase-free tube and centrifuge at 300  $\times g$  for 5 min. Discard supernatant carefully and go to step III. RNA Extraction from Cultured Cell.

**III RNA Extraction from Cultured Cell**

- 1) Add 400  $\mu$ l of RB Buffer to the cell pellet and mix by vortex mixer.
- 2) Add 300  $\mu$ l of ethanol (80%) and mix immediately by using pipette.
- 3) Transfer the sample to a Binding column in a 2 ml collection tube.
- 4) Close the lid and centrifuge at  $\geq 14,000$  rpm for 20 sec.
- 5) Discard the flow-through from the collection tube and reuse the collection tube.
- 6) Add 700  $\mu$ l of RWA1 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 20 sec.
- 7) Discard the solution from the collection tube and reuse the collection tube.
- 8) Add 500  $\mu$ l of RWA2 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 20 sec.
- 9) Discard the solution from the collection tube and reuse the collection tube.
- 10) Add 500  $\mu$ l of RWA2 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 2 min.
- 11) Discard the solution from the collection tube and reuse the collection tube.
- 12) Centrifuge once more at 14,000 rpm for 1 min to completely remove ethanol, and check that there is no droplet clinging to the bottom of

**AccuPrep® Universal RNA Extraction Kit (K-3140, K-3141)**

binding column tube.

- 13) Transfer the Binding column tube to a new 1.5 ml tube for elution, add **50~200 µl** of **ER Buffer** onto Binding column tube, and wait for at least 1 min at RT (15~25°C).
- 14) Centrifuge at **10,000 rpm** for **1 min** to elute.

**IV RNA Extraction from Plant Tissue**

- 1) Grind Sample (up to 100 mg) under liquid nitrogen. Transfer the tissue powder to an appropriately sized tube and add **500 µl** of **RB Buffer** to a maximum of 100 mg tissue powder and vortex vigorously.
- 2) Incubate at **60°C** for **1~3 min**.
- 3) Centrifuge at full speed for **2 min**.
- 4) Transfer the supernatant, aqueous phase to a new microcentrifuge tube.
- 5) Add **0.5 sample volume** of **ethanol (96~100%)** and mix immediately by using pipette.
- 6) Go to **step 3** of "RNA Extraction from Cultured Cell" in page 1 and follow the instructions accordingly.

**V RNA Extraction from Animal Tissue**

- 1) Homogenize the sample (20~30 mg) with a homogenizer, place them in a new 1.5 ml tube, and add **500 µl** of **RB Buffer**.
- 2) Centrifuge the lysate for 3 min at full speed, and transfer the supernatant to a new 1.5 ml tube.
- 3) Add **200 µl** of **ethanol (96~100%)** and mix immediately by using pipette.
- 4) Go to **step 3** of "RNA Extraction from Cultured Cell" in page 1 and follow the instructions accordingly.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User's Guide of this kit.