

## MagListo™ 5M Viral DNA/RNA Extraction Kit (K-3623, K-3624)

### I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of **ER Buffer**. For short term storage, dissolved Proteinase K can be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) Dissolve **Poly(A)** with **500 µl** of **ER Buffer**. Gently mix with vortex mixer.  
Mix dissolved Poly(A) solution into **VB Buffer**. Shake it thoroughly.
- 3) Add correct amount of **absolute ethanol** to **VWM1 Buffer**.

### II Viral DNA/RNA Extraction

- 1) Add **10 µl** of **Proteinase K** solution to a 1.5 ml or 2 ml tube.
- 2) Add **200 µl** of **Serum, Plasma** or **CSF** to the tube.  
(Note) Serum and Plasma can be used below **200 µl**.
- 3) Add **200 µl** of **VB Buffer** in the tube and mix by vortexing for 10 sec.  
To ensure efficient lysis, the sample should be mixed thoroughly with VB Buffer.
- 4) Incubate at **56~60°C** for **10 min**.
- 5) Add **400 µl** of **absolute ethanol** and mix thoroughly by vortexing.
- 6) Add **200 µl** of **Magnetic NanoBeads Solution** and mix thoroughly by vortexing.  
(Note) Please vortex Magnetic NanoBeads solution well before use.
- 7) Place the tube in **MagListo™ -2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.
- 8) Without removing the tubes from rack, carefully pour the supernatant out and completely remove the remaining supernatant using paper towel by blotting.
- 9) Detach the magnet plate from **MagListo™** stand. Add **700 µl** of **VWM1 Buffer** to the tube.  
Close the cap and mix by vortexing or shaking until the beads are fully resuspended.
- 10) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 11) Without removing the tubes from **MagListo™** rack, remove the supernatant.
- 12) Repeat the above step 9~11.
- 13) Detach the magnet plate from **MagListo™** stand. Add **700 µl** of **RWA2 Buffer**.  
Mix thoroughly by vortexing, and repeat the above step 10~11.
- 14) Choose from two ways to remove residual ethanol.  
14-1. (Washing Bead) Without removing the tube from **MagListo™** rack with the magnet plate attached, add **700 µl** of **WE Buffer** to "the opposite side of bead". Close the cap and gently invert the tube twice. Discard the supernatant and completely remove the remaining supernatant by blotting.  
(Note) Do not vortex or tap the mixture after adding WE Buffer, because vigorous mixing can reduce the yield of nucleic acid.

## **MagListo™ 5M Viral DNA/RNA Extraction Kit (K-3623, K-3624)**


- 14-2. **(Drying Bead)** Add **700 µl** of **80% Ethanol**, mix thoroughly by vortexing, and repeat the above step 10~11. Completely dry the beads with the tube open at 60°C for at least 5 min. Remove the remaining supernatant by using pipette.
- 15) Detach the magnet plate from *MagListo™* stand. Add **50~100 µl** of **ER Buffer**. Mix thoroughly by vortexing, and incubate the tube at 56~60°C at least 1 min.
- 16) Attach the magnet plate to *MagListo™* stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 17) Without removing the tube from *MagListo™* rack, transfer the supernatant containing nucleic acid carefully to a sterile microcentrifuge tube.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User's Guide of this kit.

# MagListo™ 5M Viral DNA/RNA Extraction Kit (K-3623, K-3624)

## Before You Begin

- 1) Completely dissolve **Proteinase K** in **1,250 µl** of **ER Buffer**.
- 2) Dissolve **Poly(A)** with **500 µl** of **ER Buffer**. Gently mix with vortex mixer.  
Mix dissolved Poly(A) solution into **VB buffer** and shake it thoroughly.
- 3) Add correct amount of **absolute ethanol** to **VWM1 Buffer**.

| Step                 | Image                                                                                                 | Description                                                                                                                                                           |
|----------------------|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Lysis                |                      | Add <b>VB Buffer</b> to sample<br>(Sample type: Serum, Plasma, CSF)<br>- Sample ~200 µl<br>- Proteinase K 10 µl<br>- VB Buffer 200 µl                                 |
| Incubation           |                                                                                                       | 60°C Heating block for 10 min                                                                                                                                         |
| Precipitation        |                      | Add <b>absolute ethanol 400 µl</b> and mix                                                                                                                            |
| Binding              |                      | Add <b>Magnetic NanoBead 200 µl</b> and mix until the beads are fully resuspended                                                                                     |
|                      |                      | Follow these 3 steps Attach/Detach Step<br>1. Attach magnet<br>2. Discard the supernatant<br>3. Detach magnet                                                         |
| 1 <sup>st</sup> Wash |                      | Add <b>VWM1 Buffer 700 µl</b> and mix until the beads are fully resuspended                                                                                           |
|                      |                                                                                                       | Repeat the above magnet attach/detach step (step 1, 2 and 3)                                                                                                          |
| 2 <sup>nd</sup> Wash |                      | Add <b>VWM1 Buffer 700 µl</b> and mix until the beads are fully resuspended                                                                                           |
|                      |                                                                                                       | Repeat the above magnet attach/detach step (step 1, 2 and 3)                                                                                                          |
| 3 <sup>rd</sup> Wash |                     | Add <b>RWA2 Buffer 700 µl</b> and mix until the beads are fully resuspended                                                                                           |
|                      |                                                                                                       | Repeat the above magnet attach step (step 1 and 2)                                                                                                                    |
| 4 <sup>th</sup> Wash | (Washing Bead)<br> | Add <b>WE Buffer 700 µl</b> to "the opposite side of the bead pellet"<br><b>Close the cap and gently invert the rack twice</b>                                        |
|                      |                                                                                                       | Repeat the above magnet detach step (step 2 and 3)                                                                                                                    |
|                      | (Drying Bead)<br>  | Add <b>80 % Ethanol 700 µl</b> and mix<br>Repeat the above magnet attach/detach step (step 1, 2 and 3)<br>Dry the beads with the tube open at 60°C for at least 5 min |
| Elution              |                    | Add <b>ER Buffer 50~100 µl</b> and mix                                                                                                                                |
|                      |                                                                                                       | 60°C Heating block at least 1 min                                                                                                                                     |