

MagListo™ 5M PCR/Gel Purification Kit (K-3627)

I Before You Begin

- 1) Prepare the absolute isopropanol and absolute ethanol (not provided in the kit).
- 2) Gels containing $\leq 1\%$ agarose for the purification are required.
- 3) Ensure to add **29 mL** of **absolute isopropanol** to **PB Buffer** before the purification.

II Fragment DNA purification from PCR product sample

- 1) Add **5 volumes** of **PB Buffer** to the PCR product and mix thoroughly by vortexing.
(e.g., For 20 μl of PCR product, add 100 μl of PB Buffer in 1.5 ml tube.)
- 2) Add **100 μl** of **Magnetic Nano Beads Solution** and mix thoroughly by vortexing.
(**Note**) Vortex mix the Magnetic Nano Beads Solution before use.
- 3) Place the tube in **MagListo™-2 Magnetic Separation Rack**. The **MagListo™** stand is attached with the magnet plate. Then invert the tube 3~4 times gently.
- 4) Without removing the tubes from the rack, carefully pour the supernatant out and remove the remaining supernatant using the paper towel by blotting.
- 5) Detach the magnet plate from the stand. Add **700 μl** of **W2 Buffer** to the tube and mix thoroughly by vortexing.
- 6) Attach the magnet plate to the stand and invert the tube 3~4 times gently until the beads are tightly bound to the magnet.
- 7) Remove the supernatants without removing the tubes from the rack.
- 8) Detach the magnet plate from the stand. Add **700 μl** of **Absolute ethanol**, and mix thoroughly by vortexing, then repeat the above steps 6~7.
- 9) Detach the magnet plate from the stand. Completely dry the beads with the tube opened at 60°C using the Heating Block for 10 min. Remove the remaining supernatant using the pipette.
- 9) Add **30~50 μl** of **EA Buffer** and mix thoroughly by vortexing. Incubate the tube at 60°C for more than 1 min.
- 10) Attach the magnet plate to the stand and carefully transfer the supernatant containing the nucleic acid to a sterile microcentrifuge tube.








MagListo™ 5M PCR/Gel Purification Kit (K-3627)

III Fragment DNA purification from Agarose gel









- 1) Visualize the band in agarose gel stained with any nucleic acid staining chemicals and cut the gel around the DNA band of interest using a scalpel blade.
- 2) The maximum amount of gel slice in each sample is **approximately 400 mg**.
Add **3 volumes** of **FB Buffer** to the gel slice (e.g., For 200 mg gel slice, add 600 μ l of FB Buffer).
(Note) Ensure the FB Buffer color remains yellow.
- 3) Incubate the sample at **60°C** for **10 min** and invert the tube every 2~3 min during the incubation to mix.
(Note) Ensure complete gel dissolution.
- 4) Add **100 μ l** of **Magnetic Nano Beads Solution** and mix thoroughly by vortexing.
(Note) Vortex mix the Magnetic Nano Beads Solution before use.
- 5) Place the tube in **the rack** in which the stand is attached with the magnet plate and invert the tube 3~4 times gently.
- 6) Without removing the tubes from the rack, carefully pour the supernatant out and remove the remaining supernatant using the paper towel by blotting.
- 7) Detach the magnet plate from the stand. Add **700 μ l** of **FWM1 Buffer** to the tube, and mix thoroughly by vortexing.
- 8) Attach the magnet plate to the stand and invert the tube 3~4 times gently until the beads are tightly bound to the magnet.
- 9) Remove the supernatants without removing the tubes from the rack.
- 10) Detach the magnet plate from the stand. Add **700 μ l** of **W2 Buffer** to the tube and mix thoroughly by vortexing. Repeat the above steps 8~9.
- 11) Detach the magnet plate from the stand. Add **700 μ l** of **Absolute ethanol** to the tube, mix thoroughly by vortexing, and repeat the above steps 8~9.
- 12) Detach the magnet plate from the stand. Completely dry the beads with the tube opened at 60°C using the Heating Block for 10 min. Remove the remaining supernatant using the pipette.
- 13) Add **30~50 μ l** of **EA Buffer** and mix thoroughly by vortexing. Incubate the tube at 60°C for more than 1 min.
- 14) Attach the magnet plate to the stand and carefully transfer the supernatant containing the nucleic acid to a sterile microcentrifuge tube.

※ For more information, please visit www.bioneer.com and refer to the User's Guide of this kit.

MagListo™ 5M PCR/Gel Purification Kit (K-3627)

Step	Image	Description
Binding		Add PB Buffer to the PCR product and mix. - Add 5 volumes of PB buffer to the PCR product. (ex. For 20 µl of PCR product, add 100 µl of PB buffer)
		Add 100 µl of Magnetic Nano Bead and mix
		Follow these 3 Magnet Attach/Detach Steps 1. Attach the magnet 2. Discard the supernatant 3. Detach the magnet
1 st Wash		Add 700 µl of W2 Buffer and mix
		Repeat the above "magnet attach/detach step" (step 1,2 and 3)
2 nd Wash		Add 700 µl of Absolute ethanol and mix
		Repeat the above "magnet attach/detach step" (step 1,2 and 3)
Drying	 Drying Bead	Dry the beads with the tube opened at 60°C for 10 min
Elution		Add 30~50 µl of EA buffer and mix
		Use the 60°C Heating block for more than 1 min

MagListo™ 5M PCR/Gel Purification Kit (K-3627)

Step	Image	Description
Gel solubilization		Add FB Buffer to the gel slice - Gel sample: approx. 400 mg - Add 3 volumes of FB buffer to the gel slice (e.g., for 200 mg of gel slice, add 600 µl of FB buffer)
		Use the 60°C Heating block for 10 min and invert the tube every 2-3 min during incubation
Binding	 	Add 100 µl of Magnetic Nano Bead and mix
		Follow these 3 Magnet Attach/Detach Steps 1. Attach the magnet 2. Discard the supernatant 3. Detach the magnet
1 st Wash		Add 700 µl of FWM1 Buffer and mix
		Repeat the above "magnet attach/detach step" (step 1,2 and 3)
2 nd Wash		Add 700 µl of W2 Buffer and mix
		Repeat the above "magnet attach/detach step" (step 1,2 and 3)
3 rd Wash		Add 700 µl of Absolute ethanol and mix
		Repeat the above "magnet attach/detach step" (step 1,2 and 3)
Drying	 Drying Bead	Dry the beads with the tube opened at 60°C for 10 min
Elution		Add 30~50 µl of EA buffer and mix
		Use the 60°C Heating block for more than 1 min