

## MagListo™ 5M Plasmid Extraction Kit (K-3600, K-3601)

### I Before You Begin

- 1) Completely dissolve RNase A powder in PM1 Buffer.  
After adding RNase A, PM1 Buffer should be stored at 4°C.
- 2) PB Buffer contains chaotropic salt. You should take the appropriate laboratory safety precautions and wear gloves when handling.
- 3) If there is any precipitate in P2 Buffer (Lysis), incubate on a heating block at 56°C.

### II DNA Extraction from Plasmid

- 1) Harvest the 1~5 ml (high-copy number plasmid) or 2~10 ml (low-copy number plasmid) of cultured *E. coli* cells by centrifugation at 6,000 x g for 15 min at 4°C and completely remove the media by pipetting.
- 2) Add **200 µl** of **PM1 Buffer** to the harvested cell pellet and resuspend by pipetting.  
(Note) PM1 Buffer contains Magnetic Nanobeads, Please shake well the bottle before use.
- 3) Add **200 µl** of **P2 Buffer** and invert the tube 3~4 times gently.  
(Note) Do not vortex but just invert gently. Vortexing can cause shearing of genomic DNA.
- 4) Add **200 µl** of **PC Buffer** and invert the tube 3~4 times gently.  
(Note) Again, be cautious not to shear genomic DNA. Genomic DNA and cell debris will form an insoluble complex.
- 5) Place the tube in **MagListo™ -2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.  
(Optional) For better yield, centrifuge the tube at >13,000 rpm for 1 min in a micro-centrifuge. Transfer cleared lysate to a new **2 ml (or 1.5 ml)** tube and go step 8.
- 6) Without removing the tube from **MagListo™** rack, take the supernatant only and transfer into a new 1.5 ml tube.
- 7) Repeat the above step 5 and 6 for clearer supernatant.
- 8) Add **600 µl** of **PB Buffer** and then **50 µl** of **Magnetic NanoBead solution**.  
Invert the tube 3~4 times gently.  
(Note) Please vortex Magnetic NanoBead Solution well before use.
- 9) Place the tube in **MagListo™ -2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.
- 10) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 11) Detach the magnet plate from the **MagListo™** stand. Add **1 ml** of **80% ethanol**.  
Invert the tube until the beads are fully resuspended.
- 12) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 13) Without removing the tube from **MagListo™** rack, remove the supernatant.

## **MagListo™ 5M Plasmid Extraction Kit (K-3600, K-3601)**











- 14) Detach the magnet plate from the *MagListo™* stand. Add **1 ml of 80% ethanol**.  
Invert the tube until the beads are fully resuspended. Repeat to the above step 12~13.
- 15) Without removing the tube from *MagListo™* rack, add **700 µl of WE Buffer** to “the opposite side of bead”. Close the cap and gently invert the tube twice.
- 16) Discard the supernatant and completely remove the remaining supernatant by blotting.
- 17) Detach the magnet plate from the *MagListo™* rack. Add **100 µl of EA Buffer**.  
Mix thoroughly by vortexing and incubate at 60°C for 1 min.
- 18) Attach the magnet plate to *MagListo™* stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 19) Without removing the tube from *MagListo™* rack, transfer the supernatant containing DNA carefully to a sterile microcentrifuge tube.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User's Guide of this kit.

# MagListo™ 5M Plasmid Extraction Kit (K-3600, K-3601)

## Before You Begin

- 1) Completely dissolve RNase A powder in PM1 Buffer.
- 2) If there is any precipitate in P2 Buffer, incubate on a heating block at 56°C.

Step	Image	Description
Resuspension		Add PM1 Buffer to sample and resuspend it completely - Sample : Cell pellet - PM1 Buffer : 200 µl
Lysis		Add P2 Buffer 200 µl and mix by inverting
Neutralization		Add PC Buffer 200 µl and mix by inverting
		Attach magnet and transfer the supernatant into a new 2 ml tube
Binding		Add PB Buffer 600 µl and Magnetic NanoBead 50 µl and mix until the beads are fully resuspended
		Follow these 3 Magnet Attach/Detach Step 1. Attach magnet 2. Discard the supernatant 3. Detach magnet
1 <sup>st</sup> Wash		Add 80% Ethanol 1 ml and mix by inverting until the beads are fully resuspended
		Repeat the above magnet attach/detach step (step 1, 2 and 3)
2 <sup>nd</sup> Wash		Add 80% Ethanol 1 ml and mix by inverting until the beads are fully resuspended
		Repeat the above magnet attach step (step 1, 2)
3 <sup>rd</sup> Wash		Add WE Buffer 700 µl to "the opposite side of the bead pellet". <b>Close the cap and gently invert the rack twice</b>
		Repeat the above magnet detach step (step 2 and 3)
Elution		Add EA Buffer 100 µl and mix
		60°C Heating block for 1 min