

## MagListo™ 5M Universal RNA Extraction Kit (K-3613)

### I Before You Begin

- 1) Add 10 µl β-mercaptoethanol (not provided) per 1 ml RD Buffer.
- 2) Prepare additional ethanol (80 % and 100 %) that is not included.

### II RNA Extraction from Cultured Cell

- 1) Centrifuge the cultured cells ( $10^4$ ~ $10^6$ ) for 5 min at 300 x g.  
Discard the supernatant carefully without disturbing the pellet.
- 2) Add **500 µl** of **RD Buffer** to the cell pellet and mix thoroughly by vortexing.
- 3) Add **300 µl** of **ethanol (96~100 %)** and mix immediately by using pipette.
- 4) Add **100 µl** of **Magnetic NanoBeads Solution** and mix thoroughly by vortexing.  
(Note) Please vortex Magnetic NanoBeads solution well before use.
- 5) Place the tube in **MagListo™-2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.
- 6) Without removing the tubes from rack, carefully pour the supernatant out and completely remove the remaining supernatant using paper towel by blotting.
- 7) Detach the magnet plate from **MagListo™** stand. Add **800 µl** of **RWM1 Buffer** to the tube.  
Close the cap and mix by vortexing or shaking until the beads are fully resuspended.
- 8) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 9) Without removing the tubes from **MagListo™** rack, remove the supernatant.
- 10) Add **800 µl** of **RWA2 Buffer**, mix thoroughly by vortexing, and repeat the above step 8~9.
- 11) Choose from two ways to remove residual ethanol.  
11-1. (**Washing Bead**) Without removing the tube from **MagListo™**, Add **700 µl** of **WE Buffer** to “the opposite side of bead”. Close the cap and gently invert the tube twice. Discard the supernatant and completely remove the remaining supernatant by blotting.  
11-2. (**Drying Bead**) Add **800 µl** of **80 % ethanol**, mix thoroughly by vortexing, and repeat the above step 8~9. Completely dry the beads with the tube open at 60 °C for at least 5 min. Remove the remaining supernatant by using pipette.
- 12) Detach the magnet plate from **MagListo™** stand. Add **50~100 µl** of **ER Buffer**.  
Mix thoroughly by vortexing, and incubate the tube at 60 °C for 1 min.
- 13) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- 14) Without removing the tube from **MagListo™** rack, transfer the supernatant containing RNA carefully to a sterile microcentrifuge tube.

## MagListo™ 5M Universal RNA Extraction Kit (K-3613)

### III RNA Extraction from Plant Tissue

- 1) Add **500 µl** of **RD Buffer** to a maximum of 100 mg tissue powder and vortex vigorously.
- 2) Incubate at 60°C for 1~3 min.
- 3) Centrifuge at full speed for 2 min.
- 4) Transfer the supernatant to a new microcentrifuge tube carefully without disturbing the pellet.
- 5) Add **300 µl** of **ethanol (96~100 %)** and mix immediately by using pipette.
- 6) Go to **step 4** of "RNA Extraction from Cultured Cell" in page 1 and continue the extraction process.

### IV RNA Extraction from Animal Tissue












- 1) Homogenize the sample (20~30 mg) with a homogenizer, place them in a clean 1.5 ml tube, and add **500 µl** of **RD Buffer**.
- 2) Centrifuge the lysate at full speed for 3 min.
- 3) Transfer the supernatant, aqueous phase to a new microcentrifuge tube.
- 4) Add **300 µl** of **ethanol (96~100 %)** and mix immediately by using pipette.
- 5) Go to **step 4** of "RNA Extraction from Cultured Cell" in page 1 and continue the extraction process.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User's Guide of this kit.

# MagListo™ 5M Universal RNA Extraction Kit (K-3613)

## Before You Begin

- 1) Add 10  $\mu$ l  $\beta$ -mercaptoethanol per 1 ml RD Buffer.
- 2) Prepare additional ethanol (80% and 100%) that is not included.

Step	Image	Description
Lysis		Add Reagents according to sample type <b>1. HeLa cell</b> - Cell pellet (~1x10 <sup>6</sup> ) - RD Buffer 500 $\mu$ l
		<b>2. Plant Tissue</b> - Tissue ~100 mg - RD Buffer 500 $\mu$ l - 60°C Heating block for 1~3 min
		<b>3. Animal Tissue</b> - Tissue ~20 mg - RD Buffer 500 $\mu$ l
Precipitation		Add absolute ethanol 300 $\mu$ l and mix
Binding		Add Magnetic NanoBead 100 $\mu$ l and mix until the beads are fully resuspended
		Follow these 3 Magnet Attach/Detach Step 1. Attach magnet 2. Discard the supernatant 3. Detach magnet
1 <sup>st</sup> Wash		Add RWM1 Buffer 800 $\mu$ l and mix until the beads are fully resuspended
		Repeat the above magnet attach/detach step (step 1, 2 and 3)
2 <sup>nd</sup> Wash		Add RWA2 Buffer 800 $\mu$ l and mix until the beads are fully resuspended
		Repeat the above magnet attach step (step 1 and 2)
3 <sup>rd</sup> Wash		Add WE Buffer 800 $\mu$ l to "the opposite side of the bead pellet". Close the cap and gently invert the rack twice
		Repeat the above magnet detach step (step 2 and 3)
		Add 80% Ethanol 800 $\mu$ l and mix Repeat the above magnet attach/detach step (step 1,2 and 3) Dry the beads with the tube open at 60°C for at least 5 min
Elution		Add ER Buffer 50~100 $\mu$ l and mix
		60°C Heating block for 1 min