

MagListo™ Universal RNA Extraction Kit for ExiPrep™96 Lite

I Before You Begin

- 1) Add 10 µl of β-mercaptoethanol per 1 ml of RD Buffer.
- 2) Prepare 80% and 100% of additional ethanol (not provided).

II Sample Preparation

A. RNA Extraction from Cultured Cell

- 1) Centrifuge the cultured cells (10^4 - 10^6 cells) for 5 min at 300 x g. Discard the supernatant carefully without disturbing the pellet.
- 2) Add **500 µl of RD Buffer** to the cell pellet and mix thoroughly by vortexing.
- 3) Transfer the lysate to a new 96-well dome plate carefully without disturbing the pellet.

B. RNA Extraction from Plant Tissue

- 1) Add **500 µl of RD Buffer** to a maximum of 100 mg tissue powder and vortex vigorously.
- 2) Incubate at 60°C for 1-3 min.
- 3) Centrifuge at full speed for 2 min.
- 4) Transfer the supernatant to a new 96-well dome plate carefully without disturbing the pellet.

C. RNA Extraction from Animal Tissue

- 1) Disrupt (or homogenize) the sample (20–30 mg), place them in a 1.5 ml tube, and add **500 µl of RD Buffer**.
- 2) Centrifuge the lysate at full speed for 3 min.
- 3) Transfer the supernatant to a new 96-well dome plate carefully without disturbing the pellet.

III Loading the Kit to the Instrument

- 1) Add **300 µl of 100% Ethanol** to a 96-well dome plate containing lysate using multichannel pipette.
- 2) Aliquot the solution from *MagListo™* Universal RNA Extraction Kit to each of the new 96-well dome plate using multichannel pipette.

Cartridge No.	Solution	Volume
①	Lysate + Ethanol	Up to 800 µl
②	Magnetic Nano Bead solution	100 µl
③	RWM1 Buffer	800 µl
④	RWA2 Buffer	800 µl
⑤	80% Ethanol	800 µl
⑥	100% Ethanol	1 ml
⑦	ER Buffer	100 µl

MagListo™ Universal RNA Extraction Kit for *ExiPrep™* 96 Lite

- 3) Press the 'Plate' Button on the instrument.
- 4) Place the Magnetic Rod Cover to the Magnetic Rod.
- 5) Place the plate onto the proper position of the base plate.
- 6) Press the 'Standard Protocol' Button and select '**K-3613_Universal RNA (V1.0)**'.
- 7) Press the 'Run' Button to start the selected protocol.