

MagListo™ Forensic Sample DNA Extraction Kit for ExiPrep™ 96 Lite

I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease-free water. For short term storage, dissolved Proteinase K should be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) Completely dissolve **Poly(A)** in **1,000 µl** of nuclease-free water. Divide the solution into aliquots of appropriate amount and store the aliquots at **-20°C**. Repeated freeze-thawing of the aliquots more than 3 times will result in the degradation of RNA.
- 3) Prepare 80% and 100% of additional ethanol (not provided).
- 4) **(Optional)** Prepare 1 M DTT (Dithiothreitol, not provided) for lysis of samples including hair, nail clipping or semen stains.

II Sample Preparation

A. DNA Extraction from Buccal swab, hair, nail clipping or semen stain

- 1) Apply the forensic sample to a 1.5 ml tube.
- 2) Add **250 µl** of **TL Buffer** and **10 µl** of **Proteinase K** and mix thoroughly by vortexing.
- 3) **(Optional)** Add **20 µl** of 1 M DTT and mix thoroughly by vortexing. **If the sample is hair, nail clipping or semen stains, this step is necessary to increase the yield.**
- 4) Incubate at **60°C** for **at least 1 hour**.

(Note) Incubation time for complete lysis varies depending on the type of sample used and the age of starting materials. If the sample is nail clipping, bone or is considered old, extend incubation time up to overnight.

- 5) Add **250 µl** of **GB Buffer** and mix thoroughly by vortexing.
- 6) **(Optional)** Add **2 µl** of **dissolved Poly(A)** to improve extraction efficiency of DNA.
- 7) Incubate at **60°C** for **20 min**.
- 8) Centrifuge at 13,000 rpm for 1 min to obtain clear lysate.
- 9) Transfer **500 µl** of **lysate** to a new 96-well dome plate.

B. DNA Extraction from Small amount of Tissue

- 1) Disrupt (or homogenize) **the sample (~10 mg)** and place them in a 1.5 ml tube.
- 2) Add **90 µl** of **TL Buffer** and **10 µl** of **Proteinase K** and mix thoroughly by vortexing.
- 3) Incubate at **60°C** **until the sample is completely lysed**.
- 4) Add **100 µl** of **GB Buffer** and mix thoroughly by vortexing.
- 5) Transfer the lysate to a new 96-well dome plate.

C. DNA Extraction from Urine

- 1) Centrifuge the urine and discard the supernatant.
 - a. (~ 2 ml) Centrifuge at 8,000 rpm (6,000 x g) for 2 min.
 - b. (~ 15 ml) Centrifuge at 3,500 rpm (2,000 x g) for 10 min.

MagListo™ Forensic Sample DNA Extraction Kit for ExiPrep™96 Lite

- 2) Add 500 µl of 1X PBS and mix thoroughly by vortexing. Transfer the sample to a 1.5 ml or 2 ml tube.
- 3) Centrifuge at 8,000 rpm for 2 min and discard the supernatant.
- 4) Add 250 µl of TL Buffer and 10 µl of Proteinase K and mix thoroughly by vortexing.
- 5) (Optional) Add 10 µl of 1 M DTT and mix thoroughly by vortexing. **If the urine is expected to contain sperm cells, this step will improve yield.**
- 6) Incubate at 60°C for at least 1 hour.
- 7) Add 250 µl of GB Buffer and mix thoroughly by vortexing.
- 8) Incubate at 60°C for 10 min.
- 9) Transfer the lysate to a new 96-well dome plate.

D. DNA Extraction from Small Volumes of Blood and Saliva

- 1) Apply 1-100 µl of sample to a 1.5 ml tube.
- 2) Add 10 µl of Proteinase K.
- 3) Add 100 µl of GB Buffer and mix thoroughly by vortexing.
- 4) Incubate at 60°C for 10 min.
- 5) Transfer the lysate to a new 96-well dome plate.

III Loading the Kit to the Instrument

- 1) Add a proper volume of 100% Ethanol to the 96-well dome plate containing lysate using multichannel pipette.

| Sample | Lysate | | | | | 100% Ethanol |
|--|-----------|------------|-----------|---------------------|--------------------|--------------|
| | TL Buffer | Protease K | GB Buffer | 1 M DTT | Poly(A) | |
| Buccal swab/hair/nail clipping/ semen stain | 250 µl | 10 µl | 250 µl | 20 µl (optional) | 2 µl (optional) | 500 µl |
| Small amount of tissue | 90 µl | | 100 µl | - | - | 200 µl |
| Urine | 250 µl | | 250 µl | 10 µl (optional) | - | 500 µl |
| Small volume of Blood and Saliva | - | | 100 µl | - | - | 200 µl |

- 2) Aliquot the solution from MagListo™ Forensic Sample DNA Extraction Kit to each of the new 96-well dome plate using multichannel pipette.

| Cartridge No. | Solution | Volume |
|---------------|-----------------------------|------------|
| ① | Lysate + Ethanol | Up to 1 ml |
| ② | Magnetic NanoBeads solution | 100 µl |
| ③ | WM1 Buffer | 700 µl |
| ④ | W2 Buffer | 700 µl |
| ⑤ | 80% Ethanol | 700 µl |
| ⑥ | 80% Ethanol | 700 µl |
| ⑦ | EA Buffer | 100 µl |

MagListo™ Forensic Sample DNA Extraction Kit for *ExiPrep™*96 Lite

- 3) Press the 'Plate' Button on the instrument.
- 4) Place the Magnetic Rod Cover to the Magnetic Rod.
- 5) Place the plate onto the proper position of the base plate.
- 6) Press the 'Standard Protocol' Button and select '**K-3603/K-3615 Genomic DNA (V1.0)**'.
- 7) Press the 'Run' Button to start the selected protocol.