

AccuPrep® Viral RNA Extraction Kit (K-3033)

I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease free water. For short term storage, dissolved Proteinase K should be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) Dissolve **Poly (A)** with **500 µl** of **ER Buffer**. Gently mix with vortex mixer. Mix dissolved Poly (A) solution into **VB buffer**. Shake it thoroughly.
- 3) Add correct amount of **absolute ethanol** to **VW1 Buffer**.
- 4) Before starting extraction process, heat the **ER Buffer** at **56~60°C**.
- 5) The protective seal in **BST Solution** should be completely removed. **BST Solution** may be discolored, but it does not affect nucleic acid extraction.

II Viral RNA Extraction

- 1) Add **10 µl** of **Proteinase K** to a 1.5 ml or 2 ml tube.
- 2) Add **200 µl** of **Serum, Plasma, Urine** or **CSF** to the tube. If the sample is swab, then add 1X PBS (phosphate buffered saline) and vortexing it to use only the supernatant.
- 3) Add **300 µl** of **VB Buffer** in the tube and mix by vortexing for 10 sec. To ensure efficient lysis, the sample should be mixed thoroughly with VB Buffer.
- 4) Incubate at **56~60°C** for **10 min**.
- 5) Add **300 µl** of **isopropanol**, lightly vortex for about 10 sec. Spin down for 5 sec to down the lysate clinging to the walls and lid of the tube.
- 6) Add **100 µl** of **BST Solution** to the Binding column tube (fit in a collection tube) and centrifuge for **30 sec** at **13,000 rpm**.
- 7) Discard the solution from the collection tube and reuse the collection tube.
- 8) Transfer the lysate into the Binding column tube not getting the lid wet.
- 9) Close the tube and centrifuge at **13,000 rpm** for **1 min**.
(Option) If the liquid has not completely passed the column following centrifugation, then centrifuge again until the liquid completely passes through.
- 10) Discard the solution from the collection tube and reuse the collection tube.
- 11) Add **500 µl** of **VW1 Buffer** to the column, close the lid, and centrifuge at **13,000 rpm** for **1 min**.
- 12) Discard the solution from the collection tube and reuse the collection tube.
- 13) Add **600 µl** of **RWA2 Buffer**, to the column, close the lid, and centrifuge at **13,000 rpm** for **1 min**.
- 14) Discard the solution from the collection tube and reuse the collection tube.

AccuPrep[®] Viral RNA Extraction Kit (K-3033)

- 15) Centrifuge once more at **13,000 rpm** for **1 min** to remove ethanol completely.
(**Caution**) Make sure that there is no droplet hanging from the bottom of the Binding column. Residual RWA2 Buffer left in the Binding column may cause problems in later applications.
- 16) Transfer the Binding column to a new 1.5 ml tube for elution, add **100 µl** of **ER Buffer**, and let stand for 1 min to allow the buffer to permeate the column.
(**Option**) We recommend letting stand for about 5 min to increase RNA yield. You can also increase yield by heating the ER Buffer at about 60°C before adding to the column.
- 17) Elute by centrifuge at **13,000 rpm** for **1 min**. The eluted RNA solution can directly be used, or stored at -70°C for longer storage.

III Viral DNA Extraction

- 1) Add **10 µl** of **Proteinase K** to a 1.5 ml or 2 ml tube.
- 2) Add **200 µl** of **Serum, Plasma, Urine** or **CSF** to the tube. If the sample is swab, then add 1X PBS (phosphate buffered saline) and vortexing it to use only the supernatant.
- 3) Add **200 µl** of **VB Buffer** in the tube and mix by vortexing for 10 sec. To ensure efficient lysis, the sample should be mixed thoroughly with VB Buffer.
- 4) Incubate at **56~60°C** for **10 min**.
- 5) Add **400 µl** of **isopropanol**, lightly vortex for about 10 sec. Spin down for 5 sec to down the lysate clinging to the walls and lid of the tube.
- 6) Add **100 µl** of **BST Solution** to the Binding column tube (fit in a collection tube) and centrifuge for **30 sec** at **13,000 rpm**.
- 7) Discard the solution from the collection tube and reuse the collection tube.
- 8) Go to **step 8** of “Viral RNA Extraction” in page 1 and continue the instructions accordingly.

※ For more information, please visit www.bioneer.com and refer to the User's Guide of this kit.