

[Cat. No.] **K-2036, K-2037**

Introduction

AccuPower® PCR PreMix (Negative Dye) is the powerful technology for convenient and easy performance of DNA amplification. This product contains vacuum-dried components including *Top* DNA Polymerase, dNTPs, reaction buffer, and stabilizer. It simplifies preparation of reaction mixture by adding template DNA and primers without any extra process. After the reaction, loading-dye mixture must be added to samples, and then loaded on agarose gel for analysis.

Applications

- Conventional PCR
- Primer extension
- TA cloning
- Gene sequencing

Features & Benefits

- **Stability:** Included stabilizer enables delivery at room temperature and provides increased stability compared to solution-type products.
- **User-friendly:** Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform PCR by adding template DNA and primers.
- **Sensitivity:** Excellent sensitivity and amplification efficiency even with small amounts of DNA.
- **Reproducibility:** Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

Composition

Composition	Concentration
<i>Top</i> DNA Polymerase	1 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM
Reaction buffer with 1.5 mM MgCl ₂	1X
Stabilizer	O

Specifications

<i>Top</i> DNA Polymerase	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 10 kb

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Precautions

- This enzyme is specifically optimized for increasing base incorporation rate by inactivating 5' to 3' exonuclease activity. Therefore, this product is not recommended to use for real-time PCR using hydrolysis probe method.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols


Ordering Information

Description	Cat. No.
0.2 ml thin-wall 8-tube strips with attached cap	20 µl/rxn (-dye) K-2036
	480 tubes 20 µl/rxn (-dye) K-2037


Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols



Experimental Procedures

Steps		Procedure Details																								
1	 Preparation of reaction mixture	<p>1. Add template DNA, primers, and nuclease-free water into <i>AccuPower®</i> PCR PreMix (Negative Dye) tubes to make a total volume of 20 µl. Do not include the dried pellet.</p> <ul style="list-style-type: none"> Amount of template <table border="1"> <thead> <tr> <th>Template DNA</th> <th>Amount of template</th> </tr> </thead> <tbody> <tr> <td>Bacteriophage λ, Plasmid DNA</td> <td>100 fg-200 ng</td> </tr> <tr> <td>Total genomic DNA</td> <td>1-500 ng</td> </tr> </tbody> </table> <ul style="list-style-type: none"> Preparation of reaction mixture <table border="1"> <thead> <tr> <th>Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA</td> <td>Variable (1-10 µl)</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td>0.5-2 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td>0.5-2 µl</td> </tr> <tr> <td>Nuclease-free water</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>20 µl</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried pellet by pipetting or vortexing, and briefly spin down.</p>	Template DNA	Amount of template	Bacteriophage λ, Plasmid DNA	100 fg-200 ng	Total genomic DNA	1-500 ng	Components	20 µl reaction	Template DNA	Variable (1-10 µl)	Forward primer (10 pmol/µl)	0.5-2 µl	Reverse primer (10 pmol/µl)	0.5-2 µl	Nuclease-free water	Variable	Total volume	20 µl						
Template DNA	Amount of template																									
Bacteriophage λ, Plasmid DNA	100 fg-200 ng																									
Total genomic DNA	1-500 ng																									
Components	20 µl reaction																									
Template DNA	Variable (1-10 µl)																									
Forward primer (10 pmol/µl)	0.5-2 µl																									
Reverse primer (10 pmol/µl)	0.5-2 µl																									
Nuclease-free water	Variable																									
Total volume	20 µl																									
2	 Incubate reactions in a thermal cycler	<p>3. Perform the reaction under the following conditions.</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>20 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>45-65°C</td> <td>20 sec</td> <td>25-35 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>0.5-1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> </tbody> </table>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing	45-65°C	20 sec	25-35 cycles	Extension	72°C	0.5-1 min/kb		Final extension	72°C	3-5 min	1 cycle
Step	Temperature	Time	Cycles																							
Pre-denaturation	95°C	5 min	1 cycle																							
Denaturation	95°C	20 sec																								
Annealing	45-65°C	20 sec	25-35 cycles																							
Extension	72°C	0.5-1 min/kb																								
Final extension	72°C	3-5 min	1 cycle																							
3	 Analyze with gel electrophoresis	<p>4. After the reaction, maintain the reaction mixture at 4-8°C. The samples can be stored at -20°C until use.</p> <p>5. Load samples on agarose gel with adding a loading-dye mixture (Cat. No. C-9029, not provided), and perform gel electrophoresis for analysis.</p>																								
	 Option	<ul style="list-style-type: none"> If primer's T_m value is more than 65°C or PCR product size is more than 5 kb, follow the conditions as below. <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>20 sec</td> <td></td> </tr> <tr> <td>Annealing/Extension</td> <td>68°C</td> <td>1 min/kb</td> <td>30-35 cycles</td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> </tbody> </table>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	20 sec		Annealing/Extension	68°C	1 min/kb	30-35 cycles	Final extension	72°C	3-5 min	1 cycle				
Step	Temperature	Time	Cycles																							
Pre-denaturation	95°C	5 min	1 cycle																							
Denaturation	95°C	20 sec																								
Annealing/Extension	68°C	1 min/kb	30-35 cycles																							
Final extension	72°C	3-5 min	1 cycle																							