

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

mRNA Synthesis Service

1. Precautions

RNA products are highly susceptible to degradation due to both the instability of RNA and the presence of exogenous RNase. All steps must be conducted under a sterile, RNase-free conditions. Disinfect the pipette, reagents, tips, and tubes prior to the procedures and always wear gloves during the handling steps.

- The RNA oligonucleotides can be safely stored up to 6 months at -70°C.

2. Products information

- Dissolve the delivered RNA at 1 µg/µl concentration in DEPC treated D.W. Be sure to store the solution at -70°C. The solution may be aliquoted before its usage.

- Standard mRNA Synthesis Service

Normalized amount (µg)	Volume (µl)	Concentration (µg/µl)
100	100	1

- Complete mRNA Synthesis Service

Normalized amount (µg)	Volume (µl)	Concentration (µg/µl)
50	50	1

- Cloning linked mRNA Synthesis Service

Normalized amount (µg)	Volume (µl)	Concentration (µg/µl)
50	50	1

3. Notice

Bioneer corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice. All information provided here is subject to change without notice.